



Sep 07, 2021

PBMC Thawing Protocol

Fang Zhang¹¹Fang Zhang

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.bm9ck92w

Fang.Zhang

ABSTRACT

This protocol details methods for thawing peripheral blood mononuclear cells (PBMC).

For a protocol detailing Culture and Stimulation, please view the following: [PBMCs Culture and Stimulation](#).

For a protocol detailing Cell Staining for Flow Cytometry Assay, please view the following: [Cell Staining for Flow Cytometry Assay](#).

ATTACHMENTS

[PBMC_Thawing_Protocol.pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.bm9ck92w

PROTOCOL CITATION

Fang Zhang 2021. PBMC Thawing Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bm9ck92w>

KEYWORDS

PBMC, thawing, cell culture, RPMI

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 10, 2020

LAST MODIFIED

Sep 07, 2021

OWNERSHIP HISTORY

Oct 10, 2020 Liz Brydon Protocols.io

Feb 23, 2021 Fang.Zhang

PROTOCOL INTEGER ID

43012

MATERIALS TEXT

Reagents:

- *PBMC washing medium*
 - A: RPMI-1640 with 5 to 10% FBS, 100 U/ml penicillin, 100 U/ml streptomycin, 1 mM sodium pyruvate, 2 mM L-glutamine
 - B: 50% X-vivo 15 medium (Lonza) + 25U/ benzonase
 - C: 1X CTL-Anti-Aggregate-Wash™ (CTL)
- *PBMC complete culture medium*
 - RPMI-1640 with 10% FBS, 100 U/ml penicillin, 100 U/ml streptomycin, 1 mM sodium pyruvate, 2 mM L-glutamine

Materials:

- 50 mL falcon tube (Fisher scientific #14-432-22)
- Trypan blue
- Water bath
- Dry ice
- 70% EtOH

SAFETY WARNINGS

Please refer to Safety Data Sheets (SDS) for health and environmental hazards.

- 1 Warm Washing buffer and medium to **37 °C** in a water bath.
- 2 Remove vials from liquid nitrogen and transport them to the lab on dry ice.
- 3 Thaw frozen vials, only 1 vial at a time, in a **37 °C** water bath. When cells are nearly completely thawed, carry the vials to the hood and swab them with 70% EtOH.
- 4 Gently remove PBMCs (avoid pipetting up and down, as the cells are very fragile at this stage) and transfer the cells into a 50 mL falcon tube (Fisher scientific #14-432-22) containing **25 mL warmed washing buffer**.



Use **1 mL washing medium** to rinse out the cryovial and gently mix the cells by inverting the 50 mL Falcon tube ~5x.






Wash 1:

Spin the cells: **400 x g, Room temperature , 00:08:00** . Pour off the supernatant.



Wash 2:

Suspend the cell pellet in  **1 mL prewarmed medium** (dropped slowly along the side of the tube) and resuspend the cell pellets, add  **9 mL complete medium** . Spin the cells:  **400 x g, Room temperature , 00:08:00** .

8 

If cells were thawed in the presence of benzonase, perform an additional wash with culture medium in the absence of benzonase.

9 

Count cells and determine viability by Trypan blue staining.