

Nov 12, 2020

low input ChIP-seq of immune cells

Wiebke Nahrendorf¹, Philip J Spence¹¹Institute of Immunology & Infection Research, University of Edinburgh, Edinburgh, United Kingdom

1 Works for me dx.doi.org/10.17504/protocols.io.bja3kign

Wiebke Nahrendorf

ABSTRACT

Cells can stably (and heritably) alter their gene expression profile through epigenetic modifications. Histones package DNA into chromatin and can be post-translationally modified - most prominently by methylation and acetylation. These histone modifications alter chromatin structure and DNA accessibility. We optimised a protocol for reliable high quality chromatin immunoprecipitation followed by DNA sequencing (ChIPseq) starting with just 60,000 monocytes isolated directly from mouse tissues by flow sorting. Our protocol can easily be adapted to other mouse or human cell types to interrogate the genome-wide distribution of histone modifications or transcription factor binding sites in immune cells directly *ex vivo*.

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Wiebke Nahrendorf, Alasdair Ivens and Philip J Spence Inducible mechanisms of disease tolerance provide an alternative strategy of acquired immunity to malaria. *bioRxiv* 2020 doi.org/10.1101/2020.10.01.322180

DOI

dx.doi.org/10.17504/protocols.io.bja3kign

PROTOCOL CITATION

Wiebke Nahrendorf, Philip J Spence 2020. low input ChIP-seq of immune cells. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bja3kign>

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Wiebke Nahrendorf, Alasdair Ivens and Philip J Spence Inducible mechanisms of disease tolerance provide an alternative strategy of acquired immunity to malaria. *bioRxiv* 2020 doi.org/10.1101/2020.10.01.322180

KEYWORDS

ChIPseq, histone modification, epigenetic reprogramming, low input

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

IMAGE ATTRIBUTION

icons credit: the noun project - thenounproject.com

CREATED

Aug 03, 2020

LAST MODIFIED

Nov 12, 2020

PROTOCOL INTEGER ID

39995

MATERIALS TEXT

MATERIALS

MATERIALS

- [Chloroform-Isoamyl Alcohol](#) **Sigma**
- Aldrich Catalog #25666**
- [20X](#)
- [EvaGreen](#) **Biotium Catalog #31000**
- [Agencourt AMPure XP magnetic beads](#) **Beckman**
- Coulter Catalog #A63880** **Step 117**
- [Agilent High Sensitivity DNA Kit](#) **Agilent**
- Technologies Catalog #5067-4626** **Step 113**
- [cOmplete ULTRA Tablets, Mini, EDTA-free,](#)
- EASYPack Roche Catalog #05 892 791 001**
- [Sodium Butyrate 500 mg](#) **Stemcell**
- Technologies Catalog #72242**
- [Glycine](#) **Sigma Catalog #50046**
- [Qubit® Assay Tubes](#) **Life**
- Technologies Catalog #Q32856**
- [DynaMag™-2 Magnet](#) **Life**
- Technologies Catalog #12321D** **Step 39**
- [4% Paraformaldehyde in PBS](#) **Alfa**
- Aesar Catalog #J61899-AK**
- [DNA LoBind Tubes, 1.5](#)
- mL Eppendorf Catalog #0030108051**
- [UltraPure®; DEPC-treated Water](#) **Thermo**
- Fisher Catalog #10813012**
- [IMDM](#) **Thermo**
- Fisher Catalog #12440053**
- [UltraPure®; Phenol:Chloroform:Isoamyl Alcohol \(25:24:1, v/v\)](#) **Thermo**
- Fisher Catalog #15593049**
- [PBS, pH 7.2](#) **Thermo**
- Fisher Catalog #20012019**
- [RNase Cocktail®; Enzyme Mix](#) **Thermo**
- Fisher Catalog #AM2286**
- [DNA Zap®; PCR DNA Degradation Solutions](#) **Thermo**
- Fisher Catalog #AM9890**
- [Qubit®; dsDNA HS Assay Kit](#) **Thermo**
- Fisher Catalog #Q32851** **Step 112**
- [Eppendorf Safe-Lock Tubes 1.5 mL PCR clean colorless 500](#)
- tubes Eppendorf Catalog #022363212**
- [True MicroChIP](#)
- kit Diagenode Catalog #C01010130** **Step 19**

[1.5 ml Bioruptor Pico Microtubes with Caps](#) **Diagenode Catalog #C30010016**
[Corning 15mL PP Centrifuge Tubes with CentriStar Cap Sterile](#) **Corning Catalog #430791**
[Falcon® 5 mL Round Bottom High Clarity PP Test Tube with Snap Cap Sterile](#) **Corning Catalog #352063**
[Drosophila spike-in chromatin](#) **Active Motif Catalog #61686**
[spike-in antibody](#) **Active Motif Catalog #61686**
[H3K27ac Antibody - ChIP-seq](#) **Grade Diagenode Catalog #C15410196**
[H3K4me1 Antibody - ChIP-seq](#) **Grade Diagenode Catalog #C15410037-50**
[H3K9me3 Antibody - ChIP-seq](#) **Grade Diagenode Catalog #C15410193**
[DiaMag protein A-coated magnetic beads \(ChIP-seq grade\)](#) **Diagenode Catalog #C03010020**
[MicroPlex Library Preparation Kit v2 \(12 indexes\)](#) **Diagenode Catalog #C05010012** Step 99
[MicroChIP DiaPure](#) **columns Diagenode Catalog #C03040001** Step 82
[Lightcycler 480 multiwell plate 96](#) **clear Roche Catalog #05102413001**
[MicroAmp Optical 8-Cap Strip lids](#) **Thermo Fisher Catalog #4323032**
[TE Buffer Tris-EDTA 1X Solution pH 8.0](#) **Fisher Scientific Catalog #10224683**
[Premium Fetal Bovine Serum \(FBS\)](#) **Thermo Fisher Catalog #16000044**
[HBSS no calcium no magnesium no phenol red](#) **Thermo Fisher Catalog #14175053**
[Certified Molecular Biology Agarose](#) **BIO-RAD Catalog #1613100**
FLOW BUFFER
 PBS
 2 % filtered (0.22 µm) heat-inactivated FBS
 5 mM EDTA

 > can be stored for 2 weeks at **4 °C**

FIXATION BUFFER – 1 ml per 1x10⁷ cells
 PBS
 10 % filtered (0.22 µm) heat-inactivated FBS

1 % Paraformaldehyde

> make fresh just before use, equilibrate to **Room temperature** before use

QUENCH BUFFER – use at final concentration of 125 mM Glycine

2.5 M Glycine in PBS (1.87 g in 10 ml PBS; gently heat to 37 °C in water bath, then rotate for **03:00:00** at

Room temperature to dissolve completely)

> can be stored for 2 weeks at **4 °C**, equilibrate to **Room temperature** before use

COLLECTION MEDIA

IMDM (specifically formulated for mouse cells - if working with human cells use RPMI instead)

5 % filtered (0.22 µm) heat-inactivated FBS

HBSS + PIC (protease inhibitor cocktail)

1 complete ULTRA Protease Inhibitor Cocktail tablet per 10 ml HBSS, rotate for 3 h at **Room temperature** to dissolve completely

> store at **4 °C** for no longer than 1 week

100 mM sodium butyrate – use at final concentration of 5 mM sodium butyrate

dissolve in DEPC water

> store aliquots (single use – do not freeze/thaw) at **-20 °C** for no longer than 3 months

***Drosophila* spike in chromatin**

stock: 10 ng/ µl

dilute in DEPC water to 40 pg/ µl

> store aliquots (single use – do not freeze/thaw) at **-80 °C** for no longer than 6 months

ABSTRACT

Cells can stably (and heritably) alter their gene expression profile through epigenetic modifications. Histones package DNA into chromatin and can be post-translationally modified - most prominently by methylation and acetylation. These histone modifications alter chromatin structure and DNA accessibility. We optimised a protocol for reliable high quality chromatin immunoprecipitation followed by DNA sequencing (ChIPseq) starting with just 60,000 monocytes isolated directly from mouse tissues by flow sorting. Our protocol can easily be adapted to other mouse or human cell types to interrogate the genome-wide distribution of histone modifications or transcription factor binding sites in immune cells directly *ex vivo*.

DNA-Protein crosslinking and cell sorting

1d

- 1 prepare single cell suspensions and lyse erythrocytes at **Room temperature**

RECOMMENDATION: work in 15 ml sterile, RNase/ DNase free, non pyrogenic polypropylene conical tubes until step 11.

NOTE: this protocol was optimised to isolate monocytes from mouse spleens and bone marrow - but it can be easily adapted to work with most other mouse and human tissues and cell types. adhere to best practice for your tissue when preparing single cells suspensions and lysing erythrocytes.

- 2 count cells, then Fc block and antibody stain in **FLOW BUFFER*** (scale appropriately: stain max 2×10^7 cells in **1 mL**) for **00:20:00** at **Room temperature**

* Buffers and Solutions in bold capitals are described in detail in Materials section

NOTE: design, titrate and test your antibody panel carefully beforehand.

3 wash cells twice in PBS 🌀 **350 x g, Room temperature , 00:05:00**

4 gently resuspend cells in **FIXATION BUFFER*** (📏 **1 mL** for every 1×10^7 cells)

CRITICAL: warm FIXATION BUFFER to 🌡 **Room temperature** before use.

5 incubate for exactly ⌚ **00:10:00** at 🌡 **Room temperature** , gently flick to mix occasionally

6 add 🌡 **Room temperature** **QUENCH BUFFER** (final Glycine concentration [M] **125 Milimolar (mM)**) : for every 📏 **1 mL** **FIXATION BUFFER** added in step 4 add 📏 **50 µl** **QUENCH BUFFER**)

7 incubate for ⌚ **00:05:00** at 🌡 **Room temperature** , gently flick to mix occasionally

8 🌀 **450 x g, 4°C, 00:10:00 , slow brake**

NOTE: faster, longer centrifugations going forward, since cell velocity changes after fixation.

9 aspirate supernatant carefully (leave approx. 📏 **100 µl**) and resuspend cells in 📏 **12 mL** cold PBS

10 🌀 **450 x g, 4°C, 00:10:00 , slow brake**

11 aspirate supernatant carefully (leave approx. 📏 **100 µl**) and resuspend cells in 📏 **3 mL** cold **FLOW BUFFER**

12 ⚠

SORT 60,000 desired cells on BD FACS Aria III or similar cell sorter (85 µm nozzle, sort precision mode: purity, sample and collection chamber 🌡 **4 °C**) into 5 ml polypropylene FACS tubes with 📏 **2.5 mL** **COLLECTION BUFFER**.

RECOMMENDATION: sort several technical replicates from one biological sample and chromatin-immunoprecipitate each replicate with an antibody against a different histone modifications. in this way you will get a more detailed picture of the epigenetic landscape within each biological sample.

NOTE: always perform a **test sort** beforehand, where you set up all parameters, gates and compensation ready for your big ChIPseq sort day. check cell recovery - some cell types are very fragile and may require alteration of sort or collection parameters for optimal viability. our recovery was 40 - 60 % of sorted cells i.e. we continue the protocol with

approx 30,000 cells. always check the purity of your sort before and after your last sample (and in between if you encountered any problems): > 95% of sorted cells should fall in the gates for your population of interest and debris should be minimal.

13

🌀 **450 x g, 4°C, 00:10:00 , slow acceleration, slow brake**

14

aspirate supernatant carefully (leave approx  **50 µl** behind)


15

resuspend in  **2 mL** cold **HBSS + protease inhibitors (PIC) + $[M]5$ Milimolar (mM) sodium butyrate**

16


🌀 **450 x g, 4°C, 00:10:00 , slow acceleration, slow brake**

17

aspirate supernatant carefully (leave approx  **50 µl** behind)

18



flash freeze cell pellet in methanol bath on dry ice and store at  **-80 °C** for up to 3 months

cell lysis and chromatin shearing

4h

19

*the following section uses reagents from the **True MicroChIP kit** (Diagenode, #C01010130) with a modified protocol*


NOTE: work in area designated for low input DNA work, use designated pipettes with sterile RNase and DNase free tips, clean area and pipettes as well as all other equipment with 1 % Distel, 70 % Ethanol and DNAZap before starting.

 [True MicroChIP](#)

[kit Diagenode Catalog #C01010130](#)



 [TrueMicroChIP-kit-manual.pdf](#)

20

equilibrate lysis buffer tL1 to  **Room temperature** (all crystals should be dissolved) and add protease inhibitor cocktail (PIC, from True MicroChIP Kit)






for 1 sample:  **25 µl** tL1 +  **0.125 µl** PIC (1:200)

21

thaw samples slowly  **On ice** and add  **1 mL** ice-cold **HBSS + PIC**

22

🌀 **450 x g, 4°C, 00:10:00 , slow acceleration, slow brake**


- 23 aspirate supernatant carefully (leave as little behind as possible), keep pellets  **On ice**
- 24 add  **25 µl** tL1 with PIC to the cell pellet – gently vortex to resuspend and flick until bubbles form
- 25 incubate for  **00:05:00**  **On ice**
- 26 add  **75 µl** HBSS with PIC, mix by pipetting and transfer to 1.5 ml Bioruptor Pico microtubes

27 

using the Bioruptor Pico sonicate for 5 cycles 30 sec ON 30 sec OFF to shear the chromatin.


Bioruptor Pico sonication device
Sonicator
Diagenode B01060010

NOTE: sonication time and intervals are unique for each cell type. for optimal ChIPseq chromatin should be sheared into 100 - 300 bp fragments. see [QC chromatin shearing] in step 30.1 for how to optimise shearing.

- 28 briefly vortex and place  **On ice**

- 29  **14000 x g, 4°C, 00:10:00**

30 

transfer supernatant (=  **100 µl** sheared chromatin) to 1.5 ml DNA LoBind tube

sheared chromatin can be stored at  **-80 °C** for up to 8 weeks or immediately immunoprecipitated (see next section)

30.1

8h

OPTIONAL: **optimise Chromatin shearing** for cell type of interest [uses one flow-sorted technical replicate]

*use reagents from the **True MicroChIP kit** (Diagenode, #C01010130) with a modified protocol. during optimisation it may pay off to use the designated **Chromatin shearing optimization kit** – **high SDS** (Diagenode, #C01020012)*

1. start with **100 µl** sheared chromatin in 1.5 ml DNA LoBind tube
2. dilute RNase cocktail (**1 µl** + **150 µl** DEPC-treated water) and add **2 µl** to the sheared chromatin
3. incubate **01:00:00** at **37 °C**
4. add **100 µl** elution buffer tE1 and **8 µl** elution buffer tE2, mix thoroughly by pipetting
5. decrosslink proteins from DNA for at least **04:00:00** or **Overnight** in a ThermoMixer (1300 rpm) at **65 °C**
6. spin tubes briefly
7. add **200 µl** **Room temperature** Phenol/Chloroform/Isoamyl alcohol 25:24:1
8. vortex for **00:00:15**, incubate for **00:10:00** on rotating wheel at **Room temperature**
9. **14000 x g, Room temperature, 00:02:00**
10. transfer aqueous phase to new 1.5 ml DNA LoBind tube
11. add **200 µl** **Room temperature** Chloroform/Isoamyl alcohol 24:1
12. vortex for **00:00:15**, incubate for **00:10:00** on rotating wheel at **Room temperature**
13. **14000 x g, Room temperature, 00:02:00**
14. transfer aqueous phase to new 1.5 ml DNA LoBind tube (approx **150 µl**)
15. to precipitate the DNA add:
 - 15 µl** tP1
 - 2 µl** tCP1
 - 2 µl** tCP2
 - 1 µl** ice cold 100 % Ethanol
16. incubate at **-80 °C** for **00:30:00**
17. **14000 x g, 4°C, 00:25:00**
18. carefully discard supernatant and add **500 µl** ice cold 70 % Ethanol
19. **14000 x g, 4°C, 00:10:00**
20. carefully remove all supernatant and allow pellet to air-dry for max. **00:05:00**
21. resuspend pellet in **12 µl** TE: to assess DNA fragment size distribution and integrity use **1 µl** for Bioanalyzer HS DNA Chip (see step 113 for details) and **10 µl** to run on a 1.5 % Agarose gel at 100 V for **01:00:00** (use a 100 kb ladder)
22. stain and assess gel image



sheared chromatin fragments should be between 100 and 300 bp

chromatin immunoprecipitation (ChIP)

1d


31 *the following section uses reagents from the True MicroChIP kit with a modified protocol*

NOTE: work in area designated for low input DNA work, use designated pipettes with sterile RNase and DNase free tips, clean area and pipettes as well as all other equipment with 1 % Distel, 70 % Ethanol and DNAzap before starting.

32 add protease inhibitor cocktail (True MicroChIP Kit) to Chip buffer tC1

for 1 sample:  100 µl tC1 +  0.5 µl PIC

33 add  100 µl tC1 + PIC to  100 µl sheared chromatin

34 



OPTIONAL: to normalise for technical variation between samples from this point onwards spike in a small amount of *Drosophila melanogaster* chromatin into each sample. an antibody against the *Drosophila*-specific histone variant H2Av then reliably pulls down a fraction of the *Drosophila* chromatin. this should happen consistently across all samples. after sequencing, the ratio of data mapping to the *Drosophila* genome vs your organisms genome creates a normalisation factor for each sample. you can then normalise your experimental tag counts by this factor. for more information: <https://www.activemotif.com/catalog/1091/chip-normalization>

add  140 pg *Drosophila* spike-in chromatin

(if you do this also add  0.3 µg spike-in antibody in step 36)




Egan B, Yuan CC, Craske ML, Labhart P, Guler GD, Arnott D, Maile TM, Busby J, Henry C, Kelly TK, Tindell CA, Jhunjhunwala S, Zhao F, Hatton C, Bryant BM, Classon M, Trojer P (2016). An Alternative Approach to ChIP-Seq Normalization Enables Detection of Genome-Wide Changes in Histone H3 Lysine 27 Trimethylation upon EZH2 Inhibition.. PLoS one.
<https://doi.org/10.1371/journal.pone.0166438>

35 

CRITICAL: remove  10 µl (5 %) of sheared chromatin as **input sample**, store at  4 °C in a 1.5 ml DNA LoBind tube.

NOTE: *input samples* (one for each biological replicate) are essential to analyse ChIPseq. it is the measurement of epigenetic landscape in your cells before immunoprecipitation - all enrichment is measured relative to it.

36 add your antibody of interest to the remaining 95% of sheared chromatin for immunoprecipitation:

we used antibodies against H3K27ac ( 2 µg), H3K4me1 ( 5 µg) and H3K9me3 ( 1 µg) to investigate activation and repression of transcription as well as the future potential to respond to stimuli.

H3K27ac marks transcription start sites to activate transcription



H3K4me1 marks enhancers and superenhancers to promote gene expression



H3K9me3 condenses DNA into heterochromatin to silence gene expression


NOTE: other ChIPseq-grade antibodies (for example against transcription factors) can be used. titrate all antibodies for optimal ChIPseq. some protocol recommend qPCR for validating titrations; however, we find that qPCR results do not predict ChIP sequencing outcome. we instead recommend a test sequencing run to validate antibodies and the concentrations they are used at.

OPTIONAL: add  0.3 µg spike-in antibody (see step 34)

37 

incubate  **Overnight** on a rotating wheel (40 rpm) in the cold room  **4 °C**

38 next morning: prepare DiaMag Protein A-coated magnetic beads
for 1 sample mix  **10 µl** beads (pipette up and down > 20 times to get an even suspension) with  **50 µl** beads
wash buffer tBW1 in a 1.5 ml tube

39 place in the DynaMag- 2 magnet and wait for  **00:01:00**

 [DynaMag™-2 Magnet Life](#)

[Technologies Catalog #12321D](#)

40 discard the supernatant (keep tube in magnet)

41 take tube out of magnet and gently resuspend the beads in  **50 µl** tBW1

42 place in the magnet and wait for  **00:01:00**

43 discard the supernatant (keep tube in magnet)

44 take tube out of magnet and gently resuspend the beads in  **10 µl** tBW1

45 remove samples from rotating wheel (keep  **On ice**) and spin briefly to collect all liquid in the bottom of the tube

46 add  **10 µl** of washed beads





47 incubate for  **06:00:00** on a rotating wheel (40 rpm) in the cold room  **4 °C**






washes 2h



48 *the following section uses reagents from the **True MicroChIP kit** with a modified protocol*

thorough, careful washing is key for high quality ChIPseq, since it removes non-antibody bound chromatin fragments and therefore reduces background.


NOTE: work in area designated for low input DNA work, use designated pipettes with sterile RNase and DNase free tips, clean area and pipettes as well as all other equipment with 1 % Distel, 70 % Ethanol and DNAzap before starting.

- 49 place magnet  **On ice** and keep samples and all buffers ice-cold throughout
- 50 remove samples from rotating wheel (keep  **On ice**) and spin briefly to collect all liquid in the bottom of the tube
- 51 place your samples in the magnet and wait for  **00:01:00** - the beads (and the immunoprecipitated chromatin bound to them) will bind to the side of the tube facing the magnet
- 52 discard the supernatant (keep tube in magnet)
- 53 take tube out of magnet and gently resuspend the beads in  **100 µl** ice cold wash buffer tW1

NOTE: do not create bubbles.
- 54 incubate for  **00:04:00** on a rotating wheel (40 rpm) in the cold room  **4 °C**
- 55 remove samples from rotating wheel (keep  **On ice**) and spin briefly to collect all liquid in the bottom of the tube
- 56 place your samples in the magnet and wait for  **00:01:00**
- 57 discard the supernatant (keep tube in magnet)
- 58 take tube out of magnet and gently resuspend the beads in  **100 µl** ice cold wash buffer tW2

NOTE: do not create bubbles.
- 59 incubate for  **00:04:00** on a rotating wheel (40 rpm) in the cold room  **4 °C**

60 remove samples from rotating wheel (keep  **On ice**) and spin briefly to collect all liquid in the bottom of the tube

61 place your samples in the magnet and wait for  **00:01:00**


62 discard the supernatant (keep tube in magnet)

63 take tube out of magnet and gently resuspend the beads in  **100 µl** ice cold wash buffer tW3

NOTE: do not create bubbles.

64 incubate for  **00:04:00** on a rotating wheel (40 rpm) in the cold room  **4 °C**

65 remove samples from rotating wheel (keep  **On ice**) and spin briefly to collect all liquid in the bottom of the tube

66 place your samples in the magnet and wait for  **00:01:00**


67 discard the supernatant (keep tube in magnet)

68 take tube out of magnet and gently resuspend the beads in  **100 µl** ice cold wash buffer tW4

NOTE: do not create bubbles.

69 incubate for  **00:04:00** on a rotating wheel (40 rpm) in the cold room  **4 °C**

70 remove samples from rotating wheel (keep  **On ice**) and spin briefly to collect all liquid in the bottom of the tube

71 place your samples in the magnet and wait for  **00:01:00**

1m



72 discard the supernatant (keep tube in magnet)

DNA decrosslinking

18h

73 *the following section uses reagents from the **True MicroChIP kit** with a modified protocol*

NOTE: work in area designated for low input DNA work, use designated pipettes with sterile RNase and DNase free tips, clean area and pipettes as well as all other equipment with 1 % Distel, 70 % Ethanol and DNAzap before starting.

74 after removing wash buffer tW4 take tube out of magnet and gently resuspend the beads in  **200 µl** elution buffer tE1 (equilibrate to  **Room temperature** before use - tE1 should be a clear solution)

75 take the *input samples* (10 µl) you saved in step 35 out of the fridge and add  **190 µl** elution buffer tE1

76 incubate *both ChIP and input samples* for  **00:30:00** on a rotating wheel (40 rpm) at  **Room temperature**

77 remove samples from rotating wheel and spin briefly to collect all liquid in the bottom of the tube



78 place *ChIP samples* in the  **Room temperature** magnet and wait for  **00:01:00**

79 

ChIP samples: keep the tube in the magnet and transfer the supernatant (= your immunoprecipitated chromatin) to a new 1.5 ml DNA LoBind tube

80 add  **8 µl** elution buffer tE2 to *both ChIP and input samples*

81 

decrosslink proteins from DNA (for both *ChIP and input samples*)  **Overnight** in a ThermoMixer (1300 rpm) at  **65 °C**

DNA purification using Micro ChIP DiaPure columns

1h









82 *the following section uses **Micro ChIP DiaPure columns** (Diagenode, #C03040001) according to manufacturers instructions*




NOTE: work in area designated for low input DNA work, use designated pipettes with sterile RNase and DNase free tips, clean area and pipettes as well as all other equipment with 1 % Distel, 70 % Ethanol and DNAzap before starting.

 [MicroChIP DiaPure](#)

[columns Diagenode Catalog #C03040001](#)

 [MicroChIP_DiaPure_manual.pdf](#)

- 83 spin decrosslinked samples (*ChIP and input samples*: 200 µl Volume) briefly to collect all liquid in the bottom of the tube
- 84 add  **1000 µl** (5 Vol)  **Room temperature** ChIP DNA binding buffer and mix gently by pipetting
- 85 transfer  **600 µl** to the spin column in its collection tube
- 86  **10000 x g, Room temperature , 00:00:30**
- 87 discard the flow-through
- 88 transfer the remaining  **600 µl** to the spin column in its collection tube
- 89  **10000 x g, Room temperature , 00:00:30**
- 90 discard the flow-through
- 91 add  **200 µl**  **Room temperature** DNA wash buffer

CRITICAL: make sure Ethanol was added to the buffer.
- 92  **10000 x g, Room temperature , 00:00:30**
- 93 add  **200 µl**  **Room temperature** DNA wash buffer

94  **10000 x g, Room temperature , 00:00:30**

95 transfer the column to a new 1.5 ml DNA LoBind tube

96 to elute the DNA add  **15.2 µl**  **Room temperature** DNA elution buffer directly onto the column matrix and incubate for  **00:03:00** at  **Room temperature**

97  **10000 x g, Room temperature , 00:00:30** then discard the column and transfer tube with DNA  **On ice**

98 

DNA can be stored at  **-20 °C** before library preparation for up to 2 weeks

Library preparation

1d

99 *the following section uses the***Diagenode MicroPlex Library preparation kit v2** (*Diagenode, #C05010012*) ^{1d}
according to manufacturers instructions




NOTE: work in area designated for library preparation (distinct from low input DNA area), use designated library preparation pipettes with sterile RNase and DNase free tips, clean area and pipettes as well as all other equipment with 1 % Distel, 70 % Ethanol and DNA-ZAP before starting.

 [MicroPlex Library Preparation Kit v2 \(12](#)

[indexes\)](#) **Diagenode Catalog #C05010012**

 [MicroPlex-Library-Prep-Kit-v2-manual.pdf](#)



NOTE: keep samples and reagents  **On ice** throughout.



100 in clear 96 well Lightcycler 480 plate mix  **10 µl** ChIP-ed, purified DNA with  **2 µl** template preparation buffer and  **1 µl** template preparation enzyme

101 gently mix by pipetting, cap using strip lids and spin briefly to collect all liquid in the bottom of the wells

102 run on a standard PCR machine (settings: plate, 13 µl Volume, heated lid)

temperature	time
22 °C	25 min
55 °C	20 min
4 °C	∞



103 transfer plate back  **On ice** as soon as PCR machine has cooled to  **4 °C** , spin briefly to collect all liquid in the bottom of the wells

104 carefully open lids and add  **1 µl** library synthesis buffer and  **1 µl** library synthesis enzyme

105 gently mix by pipetting, cap using strip-lids and spin briefly to collect all liquid in the bottom of the tube

106 incubate once more using the same standard PCR machine (settings: plate, 15 µl Volume, heated lid)


temperature	time
22 °C	40 min
4 °C	∞

107 transfer plate back  **On ice** as soon as samples PCR machine has cooled to  **4 °C** , spin briefly to collect all liquid in the bottom of the wells

108 carefully open lids and add  **30 µl** library amplification master mix:

reagent	volume/reaction
Library amplification buffer	25 µl
Library amplification enzyme	1 µl
EvaGreen	2.5 µl
Nuclease free water	1.5 µl


109 

add  **5 µl** indexing reagent (total volume 50 µl) - to avoid cross-contamination spray index lid with DNAzap and wipe dry, change gloves after each index

NOTE: carefully consider your sequencing requirements and plan which/how many libraries you are going to pool in each lane and index samples accordingly. see below for details of the standard Illumina indices supplied with the MicroPlex Library preparation kit v2 (12 indices, Diagenode #C05010012,). a kit with 48 indices is also available: MicroPlex Library Preparation Kit v2 (48 indexes, Diagenode #C05010014).

index number	index ID	index sequence
	1 iPCRtagT1	ATCACGTT
	2 iPCRtagT2	CGATGTTT
	3 iPCRtagT3	TTAGGCAT
	4 iPCRtagT4	TGACCACT
	5 iPCRtagT5	ACAGTGGT
	6 iPCRtagT6	GCCAATGT
	7 iPCRtagT7	CAGATCTG

8	iPCRtagT8	ACTTGATG
9	iPCRtagT9	GATCAGCG
10	iPCRtagT10	TAGCTTGT
11	iPCRtagT11	GGCTACAG
12	iPCRtagT12	CTTGTA

110 gently mix by pipetting, use sealing foil and  **1300 x g, 4°C, 00:02:00**


111 run on real time quantitative PCR machine (Roche Lightcycler 480) to monitor library amplification

LightCycler® 480 Instrument II
real-time quantitative PCR machine
Roche 05015278001

	temperature	time	ramp rate
extension	72 °C	3 min	3° C/sec
cleavage	85 °C	2 min	3° C/sec
denaturation	95 °C	2 min	3° C/sec
addition of indices	98 °C	20 sec	3° C/sec
	67 °C	20 sec	2.2 °C/sec
	72 °C	40 sec	3 °C/sec
<i>repeat steps 6 to 8 four times</i>			
library amplification	98 °C	20 sec	3 °C/sec
	72 °C	50 sec: record fluorescence using "single acquisition"	2.2 °C/sec
<i>repeat steps 11 & 12 for x* number of cycles</i>			
cool - hold	37 °C **	1 h	2.2 °C/sec

* monitor fluorescence after each cycle: the optimal phase is reached when Fluorescence (465-510) linearly increases to 3.5 – 4.5. at this point stop library amplification and move to step 15 - cooling. we find it takes approx 8 - 13 cycles to amplify libraries sufficiently. the exact number of cycles will depend on how much chromatin your antibody pulls down.

** 37°C is the lowest temperature the Roche Lightcycler 480 will cool to: transfer plate containing amplified library to ice 1 min after 37°C is reached

112 use  **1 µl** of amplified library to quantify the amount of DNA with **Qubit dsDNA HS assay Kit** according to manufacturers instruction

 **Qubit dsDNA HS Assay Kit Thermo**

Fisher Catalog #Q32851

Qubit Fluorometer

Fluorometer

Invitrogen

Q33238



5 - 15 ng/ μ l

if concentration is significantly lower return sample to real time PCR machine for extra amplification

113 dilute 1 μ l of amplified library in 4 μ l TE to assess DNA integrity and size distribution using Bioanalyzer **High Sensitivity DNA Kit** according to manufacturers instructions

[Agilent High Sensitivity DNA Kit](#) **Agilent**

Technologies Catalog #5067-4626

2100 Bioanalyzer Instrument

Sizing, quantification, and sample quality control of DNA, RNA, and proteins on a single platform

Agilent Technologies

G2939BA



bell curve 200 - 2000 bp, average size approx. 400 bp

[amplified_library.pdf](#)

114

pool libraries with different indices at equal molarities in 1.5 ml DNA LoBind tube (aim to reach a volume just over

100 μ l)


molecular mass of dsDNA = 660 [g/mol/bp]

calculation:

concentration [ng/ μ l] * 10^6 * 1/660 * 1/average size [bp] = molarity [nM]

115 

transfer library from plate to 1.5 ml DNA LoBind tube (keep  **On ice**)

NOTE: libraries (individual or pooled) can be stored at  **-20 °C** for up to 2 weeks

Library purification using AMPure beads



4h

116 NOTE: work in area designated for library preparation (distinct from low input DNA area), use designated library preparation pipettes with sterile RNase and DNase free tips, clean area and pipettes as well as all other equipment with 1 % Distel, 70 % Ethanol and DNA-ZAP before starting.

117 bring pooled libraries, **AMPure beads** and freshly prepared 80 % Ethanol to  **Room temperature**

 [Agencourt AMPure XP magnetic beads](#) **Beckman**

Coulter Catalog #A63880

118 resuspend AMPure beads until homogenous solution and add  **100 µl** to  **100 µl** of pooled library (1:1 ratio); mix until homogenous

119 incubate at  **Room temperature** for  **00:05:00**

120 spin briefly ( **00:00:03**) to collect all liquid in the bottom of the tube

121 place tube in DynaMag- 2 magnet, wait for  **00:02:00** until all beads are bound (solution clear)















122 discard supernatant

123 add  **300 µl** 80 % Ethanol

124 rotate the tube clockwise by 90 °, wait for  **00:00:10** and repeat 3 more times

10s

125 discard supernatant

- 126 add  300 µl 80 % Ethanol
- 127 rotate the tube clockwise by 90 °, wait for  00:00:10 and repeat 3 more times 10s
- 128 discard supernatant
- 129 spin briefly ( 00:00:03) to collect all liquid in the bottom of the tube
- 130 place tube in DynaMag- 2 magnet, wait for  00:02:00
- 131 remove all residual Ethanol
- 132 remove tube from magnet and dry the beads with lid open for max  00:02:00 in ThermoMixer ( 37 °C)
- 133 resuspend the beads in  50 µl TE, spin briefly ( 00:00:03) to collect all liquid in the bottom of the tube
- 134 place tube in magnet, wait for  00:02:00
- 135 
carefully transfer the eluted DNA to a new 1.5 ml DNA LoBind tube
- 136 use  1 µl of pooled purified library to quantify the amount of DNA with **Qubit dsDNA HS assay Kit** according to manufacturers instruction (see step 112 for details)
-  expect 10 ng/ µl (i.e. a total of 500 ng in 50 µl)
- 137 dilute  1 µl of pooled purified library in  4 µl TE to asses integrity and size distribution of libraries using Bioanalyzer **High Sensitivity DNA Kit** according to manufacturers instructions (see step 113 for details)



bell curve 200 - 2000 bp (all small fragments removed by AMPure XP bead purification), average size approx. 400 bp

[pooled_purified_library.pdf](#)

138 pooled, purified libraries can be stored at -20°C before sequencing for up to 2 months

our samples were sequenced by Edinburgh Genomics <https://genomics.ed.ac.uk> on the Illumina NovaSeq S1 yielding approx 750×10^6 100 bp paired end reads per lane. we aimed for a depth of 70×10^6 paired end reads for each sample. ChIPed samples (for all different histone modifications) and matched input sample should be sequenced on the same lane. we used the motif discovery software HOMER for data analysis (<http://homer.ucsd.edu/homer/>). our ChIPseq data is publicly available: GEO accession number [GSE150478](#).

NOTE: we subscribe to the notion that ChIPseq is *qualitative* (it can reveal the presence or absence of a histone modification at a particular genomic location) - not *quantitative* (it does not reveal biologically meaningful differences in peak height, which are often influenced by the efficiency of immunoprecipitation). please consider this when analysing your results.

Ma Z, Wang H, Cai Y, Wang H, Niu K, Wu X, Ma H, Yang Y, Tong W, Liu F, Liu Z, Zhang Y, Liu R, Zhu ZJ, Liu N (2018). Epigenetic drift of H3K27me3 in aging links glycolysis to healthy longevity in Drosophila.. eLife.
<https://doi.org/pii:e35368.10.7554/eLife.35368>

Orlando DA, Chen MW, Brown VE, Solanki S, Choi YJ, Olson ER, Fritz CC, Bradner JE, Guenther MG (2014). Quantitative ChIP-Seq normalization reveals global modulation of the epigenome.. Cell reports.
<https://doi.org/10.1016/j.celrep.2014.10.018>