



Dec 07, 2020

Fluorescence analysis using CF imager

Forked from [Fluorescence analysis using CF imager](#)Steven J Burgess¹¹University of Illinois at Urbana-Champaign*In Development*

This protocol is published without a DOI.

Ag SynBio Lab UIUC



Benjamin Haas

ABSTRACT

Measurement of photosynthesis parameters using the Chlorophyll Fluorescence Imager from Technologica
<http://www.technologica.co.uk/products/cfimager/index.html>

PROTOCOL CITATION

Steven J Burgess 2020. Fluorescence analysis using CF imager. **protocols.io**
<https://protocols.io/view/fluorescence-analysis-using-cf-imager-sbmeak6>

FORK NOTE

FORK FROM

Forked from [Fluorescence analysis using CF imager](#), Steven Burgess

KEYWORDS

Chlorophyll fluorescence analysis, photosynthesis

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 02, 2018

LAST MODIFIED

Dec 07, 2020

PROTOCOL INTEGER ID

14413

ABSTRACT

Measurement of photosynthesis parameters using the Chlorophyll Fluorescence Imager from Technologica
<http://www.technologica.co.uk/products/cfimager/index.html>

Preparing plants

- 1 Dark adapt plants for at least 20 minutes prior to taking measurements.


This is done to ensure the photosynthetic electron transport chain is fully oxidized and reaction centres are open. In an ideal situation plants are allowed to dark adapt overnight prior to measurement.

A properly adapted, healthy plant should give a Fv/Fm value of ~0.8. This has been shown to be highly stable between species. Significant deviation from this (e.g. <0.7) either suggests incomplete adaptation or stressed plant material.


- 2 Turn on the cf imager and open the FluorImager software.



- 3 Start the FluorImager software



- 4 The surface of the leaf should be 140mm from the base of the imaging chamber, and can be adjusted by lowering or raising the plant under analysis. Position plant/leaf in the chamber




Set focus

- 5 Set the focus by adjusting the dial above the chamber, and lock in position by turning the screw on the side.


Set exposure

- 6 Manually adjust the aperture as shown on the right to allow an optimal amount of light into the imager so as not to overexpose measurements




- 7 Isolate the plant or leaf of interest


right click
on image




After applying isolation background areas will be masked out in blue as shown below.



- 8 Remove residual noise from image
Sometimes the software picks up background noise as real signals, shown below as black dots on the blue background. It is advisable to mask these, otherwise they will be counted as a separate 'colony' during analysis and measurements will be recorded for each of these spots in the final data sheet. Noise can be masked by moving the cursor over the dot, pressing CTRL+left click simultaneously.



9 Load protocol



10

11



- 12 Once you are happy with the scheduled program click on the protocol icon in the toolbar (black arrow below) to start the run



- 13 Export data

