

3

Jun 21, 2022

Direct nuclear tagmentation and RNA-sequencing (DNTR-seq) V.3

Vasilios Zachariadis¹, Huaitao Cheng¹, Nathanael Andrews¹, Martin Enge¹

¹Karolinska Institutet

1



protocol .

Enge lab

Vasilios Zachariadis

Understanding how genetic variation alters gene expression - how genotype affects phenotype - is a central challenge in biology. To address this question in complex cell mixtures, we developed Direct Nuclear Tagmentation and RNA-sequencing (DNTR-seq), which enables whole genome and mRNA sequencing jointly in single cells.

<https://www.biorxiv.org/content/10.1101/2020.03.04.976530v1.full>

Vasilios Zachariadis, Huaitao Cheng, Nathanael Andrews, Martin Enge 2022. Direct nuclear tagmentation and RNA-sequencing (DNTR-seq). **protocols.io**

<https://protocols.io/view/direct-nuclear-tagmentation-and-rna-sequencing-dnt-b65yrg7w>
Vasilios Zachariadis

protocol

A highly scalable method for joint whole genome sequencing and gene expression profiling of single cells bioRxiv 2020.03.04.976530; doi: <https://doi.org/10.1101/2020.03.04.976530>

single-cell, scRNA-seq, scWGS

protocol ,

Apr 05, 2022

Jun 21, 2022

60312

☒ UltraPure™ DNase/RNase-Free Distilled Water Thermo

Fisher Catalog #10977035

☒ ERCC RNA Spike-In Mix Thermo

Fisher Catalog #4456740

☒ USB Dithiothreitol (DTT), 0.1M Solution Thermo

Fisher Catalog #707265ML

☒ Sera-Mag Speed Beads Ge

Healthcare Catalog #65152105050250

☒ RNase

Inhibitor Takara Catalog #2313A

☒ Hard-Shell® 384-Well PCR Plates thin wall skirted BioRad

Sciences Catalog #HSP3801

Bleach clean environment - to avoid DNA contamination. And RNase away or similar to avoid degraded RNAs. Prepare solutions in a strictly pre-PCR environment. Keep plates and reagents on ice unless otherwise noted.

Prepare lysis buffer plates for cell sorting

1 Prepare lysis buffer mix

NOTE: Reagents are prepared on ice, working quickly. ERCC is stored in single-use aliquots at **-80 °C**, thawed on ice and added last.

A	B	C	D
Reagent	Reaction conc.	μL per reaction	384w plate (400x)
Nuclease free H ₂ O	-	1.965	786
RNase Inhibitor (40u/μL)	1 unit/μL	0.075	30
ERCC (1:1 200 000)	-	0.075	30
Triton-X100 (10% solution)	0.2%	0.06	24
dNTP (10mM each)	2.5mM/each	0.75	300
Oligo-dT (100μM)	2.5μM	0.075	30
To dispense		3	1200

Add **3 μL lysis buffer mix** to each well. Cover with appropriate lids. Spin down.

Snap freeze on **dry ice**. Store until use at **-80 °C**

Sort single-cells

2 Sort single cells into **3 μL lysis buffer mix**

Immediately seal with appropriate seals (approved for -80C > 100C) and centrifuge at

2000 x g, 4°C, 00:05:00

Snap freeze on **dry ice**. Store until use at **-80 °C**

Separation of nuclear and cytosolic fractions

- 3 Thaw plate on ice.
Centrifuge at **500 x g, 4°C, 00:05:00**.
Keep on ice.
- 4 Transfer **2 µL** from each well of the sorted plate into an empty 384-well plate. Use a low flow rate (**2mm/s**) and an aspiration height of **0.9mm** above the bottom.

NOTE: We use the Eppendorf EpMotion 5073m benchtop liquid handler. We have successfully used other solutions, including the Hamilton STARlet, a semi-manual Gilson Platemaster 96-well pipette, and even manual 8-channel pipettes.

- 5 Spin down and freeze nuclear fraction at **-20 °C** to aid complete lysis.

If proceeding with **cDNA protocol** → step 12.

If proceeding with **DNA protocol** (step 6): spin down and snap freeze cytosolic fraction on **dry ice** and store at **-80 °C**

NOTE: We will typically proceed with cDNA synthesis, unless experimental design dictates otherwise, to avoid an additional freeze-thaw cycle for mRNAs in the cytosolic fraction.

Single-cell genomic libraries

- 6 Using plate with nuclear fraction, with remaining volume 1µl/well.

Proteinase K treatment

1. Dilute Proteinase K (stock 20mg/ml) to 0.2mg/ml by 30mM Tris-HCl pH8.0
2. Add **2 µL** diluted Proteinase K (0.2mg/ml) to each well. Makes 0.13mg/ml reaction concentration.
3. Incubate in thermocycler at:
 - **50 °C 01:00:00**
 - **80 °C 00:30:00**
 - **4 °C hold**

- 7 **Tn5 digestion**

Tn5 is produced from psfTn5 (Addgene #79107), purified to ~3mg/ml and assembled with Illumina Tn5 adapters (see **oligos**) as in *Picelli et al, 2014*.


Picelli S, Björklund AK, Reinius B, Sagasser S, Winberg G, Sandberg R (2014). Tn5 transposase and tagmentation procedures for massively scaled sequencing projects.. Genome research. <https://doi.org/10.1101/gr.177881.114>

- 7.1 Prepare 2X Tn5 Buffer. Keep assembled Tn5 enzyme (Picelli *et al*, 2014) on ice block and add last.

A	B	C	D
Reagent	Reaction conc	µL per reaction	384w plate (420x)
5X TAPS-PEG (50mM TAPS, 25mM MgCl ₂ , 40% PEG-8000)	10mM TAPS 5mM MgCl ₂ 8% PEG-8000	1.6	672
psfTn5, loaded with 50µM MEDS-A/B		0.1	42
Nuclease free H ₂ O		3.3	1386
<i>To dispense</i>		5	2100

Add  5 µL per well. Vortex and spin down plate.

NOTE: Buffer contains PEG, which is viscous. 5X TAPS-PEG buffer should be allowed to assume room temperature before dispensing to allow proper mixing.


- 7.2 Incubate in thermocycler: 🔥 55 °C ⌚ 00:10:00
Remove immediately and stop reaction by adding  2 µL per well of 0.2% SDS.
Vortex, spin down and incubate ⌚ 00:10:00 at 🔥 55 °C

8 PCR amplification and barcoding

1. Prepare PCR master-mix

A	B	C	D
Reagent	Reaction conc.	μl per reaction	384w plate
Nuclease free H ₂ O	-	2.9	1218
KAPA HiFi Buffer (5X)	1X	3.9	1638
dNTP (10mM/each)	0.3mM/each	0.6	252
KAPA enzyme (1u/μl)	0.02u/μl	0.4	168
Tween-20 (10%)	0.1%	0.2	84
To dispense		8	3360

2. Dispense  **8 μL per well**

3. Add primers/barcodes  **1.5 μL per well** (from 384-well index plates, with 3.75μM/each forward/reverse primers; see **oligos**). Total reaction volume is now 19.5μl (10μl sample + 9.5μl PCR mix and primers).

4. Vortex plate, spin down and incubate in thermocycler with the following program:

Step	Temperature	Time	Cycles
Gap fill	72°C	3 min	1x
First denature	95°C	30 sec	1x
Denature	95°C	15 sec	18x
Anneal	67°C	30 sec	
Extend	72°C	45 sec	
Final extension	72°C	4 min	1x
	4-10°C	hold	





9 Pool  **3 μL from each well** into a 1.5mL Eppendorf tube.

10 Library cleanup

We prepare SPRI-beads in 20% PEG-8000 solution as in:

https://openwetware.org/wiki/SPRI_bead_mix#Ingredients_for_50_mL_2


(optional) Take an aliquot of your pool (300μl)

1. Add 0.9X SPRI-beads in 20% PEG solution. Incubate for  **00:05:00**  **Room temperature**
2. Place on magnetic rack  **00:03:00**
3. Remove supernatant
4. Add 1 volume 80% EtOH (fresh). Incubate for  **00:00:30**
5. Remove supernatant
6. Repeat EtOH wash

7. Air dry for ☰00:10:00 - ☰00:15:00
8. Re-suspend beads thoroughly in 📄100 µL EB or TE buffer
9. Repeat cleanup (from step 1-7) and elute in 📄30 µL EB or TE buffer

11 (optional) Quality control of DNA libraries

Using Agilent HS 5000 DNA chips (or equivalent)



Pooled (and diluted) DNA-library from 384-well plate.

This library was sequenced on a NextSeq 550 loading 2.5pM based on a peak of 334bp. Sequencing was paired-end 37bp, 8bp dual index.

Reverse transcription and cDNA amplification

- 12 Following step 4, cytosolic/RNA fraction plate contains 2µl solution per well.

Primer annealing

Thaw plate. Spin down. Incubate in thermocycler at 🔥72 °C for ☰00:03:00. Remove to ice immediately.

- 13 Prepare RT master-mix

A	B	C	D
Reagent	Reaction conc.	µl per reaction	384w plate (420x)
Maxima H Minus RT (200/µl)	2u/µl	0.05	21
RNase Inhibitor (40u/µl)	1.66u/µl	0.125	52.5
5X First Strand buffer	1X	1	420
DTT (100mM)	8.33mM	0.25	105
Betaine (5M)	1.66M	1	420
MgCl ₂ (1M)	10mM	0.03	12.6
TSO (100uM)	1.66µM	0.05	21
Nuclease free H ₂ O	-	0.495	207.9
Total		3	1260

Dispense **3 µL per well** . Total reaction volume will be 5µl.

Cover plate with new film and spin down.

14 Incubate in thermocycler

⚙ 42 °C ⌚ 01:30:00

⚙ 70 °C ⌚ 00:05:00

⚙ 4 °C hold

15 cDNA preamplification

A	B	C	D
	Reaction conc.	µl per reaction	384w plate (420x)
Nuclease free H ₂ O	-	0.82	345
Kapa HiFi HotStart ReadyMix (2X)	1X	6	2520
IS_PCR primer (10µM)	0.1µM	0.12	50.4
Lambda Exonuclease (10u/µl)	0.05 units	0.06	25.2
Total		7	2940

Dispense **7 µL per well** . Total reaction volume will be 12µl.

16 Vortex, spin down. Cover with new lid. Incubate in thermocycler with the following program:

A	B	C	D
Step	Temperature	Time	Cycles
Lambda exonuclease	37°C	30 min	1x
Initial denaturation	95°C	3 min	1x
Denaturation	98°C	20 sec	18-24x
Annealing	67°C	15 sec	
Elongation	72°C	4 min	
Final elongation	72°C	5 min	
	4C	Hold	

NOTE: The number of cycles of pre-amplification will be different for different cell types. We suggest running a pilot (ideally qPCR-monitored to determine inflection point, for example by using 1X dsGreen to the reaction above)

17 cDNA cleanup

31m 30s

Using 20% SPRI-bead solution (as in step10 for DNA library cleanup).

1. Add 0.7X volume of SPRI beads per well. Mix well by pipetting
2. Incubate 🕒 00:05:00 🌡 Room temperature
3. Place on magnetic stand for 🕒 00:03:00
4. Carefully remove supernatant
5. Add 📏 40 µL 80% EtOH and incubate 🕒 00:00:30
6. Remove EtOH (without disturbing the beads)
7. Wash again with EtOH. Make sure to remove well.
8. Allow beads to air-dry for 🕒 00:05:00 - 🕒 00:10:00 . Take care not to over-dry the beads.
9. Remove plate from magnetic stand
10. Elute beads in 📏 12 µL EB or TE buffer Mix well by pipetting
11. Incubate 🕒 00:05:00 🌡 Room temperature
12. Place on magnetic plate for 🕒 00:03:00
13. *Optional: Carefully remove supernatant to the elution plate. cDNA plates can also be stored at -20C with beads.*

18 cDNA quantification

Option 1: Measure concentration of random wells using Qubit HS dsDNA, adapted to a 96-well plate reader.

1. Add 📏 98.5 µL of 1X Qubit HS dsDNA solution (or mix dye and buffer separately) to a flat-bottom, black plate


2. Add **1.5 µL** of cDNA sample
3. Add Standards (NOTE: We make a 8-step ladder from 0ng/µl → 10ng/µl Qubit Standard DNA in TE buffer)
4. Read in plate reader using 485nm excitation/528nm emission
5. Calculate cDNA concentration from linear model of Standards ladder

Option 2: Measure full plate using Qubit HS dsDNA in black, flat-bottom 384-well plate

1. Add **20 µL 1X Qubit HS dsDNA solution**
2. Add **1 µL cDNA sample**
- 3-5 as above

19 (optional) cDNA quality control

Using Agilent HS 5000 DNA chips (or equivalent)



Example of a single immune (=small) cell cDNA profile (cytosolic fraction from DNTR protocol)

20 Make cDNA dilution plate

Dilute cDNA based on average concentration from Qubit measurements.

Target concentration **150 pg per µl** in **15 µL** (optionally in same plate)

Optional: if using a 384w-plate reader, one can normalize each well to 150pg/µl with variable water addition.

cDNA tagmentation

21 Prepare Tn5 master mix

Let TAPS-PEG equilibrate at 37°C and mix well before use.

A	B	C	D
Reagent	Reaction conc.	µl per reaction	384w plate (420x)
Nuclease free H ₂ O	-	0.750	315
TAPS-PEG (50mM TAPS, 25mM MgCl ₂ , 40% PEG-8000)	10mM TAPS 5mM MgCl ₂ 8% PEG-8000	0.500	210
psfTn5, loaded with 50µM MEDS-A/B		0.250	105
Total		1.5	630

Dispense **1.5 µL per well** in a new plate (**tagmentation plate**)

- 22 Add **1 µL** cDNA (normalized to **150pg/µl**)
Mix well by vortexing plate. Cover with new lid and spin down.

- 23 Incubate in thermocycler at **55 °C** **00:10:00**
Remove immediately and stop reaction by adding **1 µL per well** of 0.1% SDS.
Vortex, spin down and incubate **00:10:00** at **55 °C**

20m

cDNA library PCR and barcoding

24 Make PCR master-mix

A	B	C	D
Reagent	Reaction conc.	µl per reaction	384w plate (420x)
H ₂ O	-	4.85	2037
KAPA HiFi Buffer (5X)	1X	2.5	1050
dNTP (10mM/each)	0.3mM/each	0.3	126
KAPA enzyme (1u/µl)	0.02u/µl	0.2	84
Tween-20 (10%)	0.12%	0.15	63
Total		8	

Dispense **8 µL per well** to **tagmentation plate** (containing 3.5µl sample after step 23)

- 25 Add primers/barcodes **1 µL per well** (from 384-well index plates, with 3.75µM/each forward/reverse primers; see **oligos**; final primer concentration 0.3µM per primer and reaction).

Total reaction volume is 12.5µl (3.5µl sample + 9µl PCR mix and primers).

26 Vortex. Spin down and cover. Incubate in thermocycler as below:

Step	Temperature	Time	Cycles
Gap fill	72°C	3 min	1x
First denature	95°C	30 sec	1x
Denature	95°C	15 sec	12x
Anneal	67°C	30 sec	
Extend	72°C	45 sec	
Final extension	72°C	4 min	1x
	4-10°C	hold	









cDNA library pooling and clean-up

27 Pool  2.5 µL from each well to an 1.5ml Eppendorf tube


28 Library cleanup (as for DNA libraries)

We prepare SPRI-beads in 20% PEG-8000 solution as in:

https://openwetware.org/wiki/SPRI_bead_mix#Ingredients_for_50_mL_2

1. Add 0.9X SPRI-beads in 20% PEG solution. Incubate for  00:05:00  Room temperature
2. Place on magnetic rack  00:03:00
3. Remove supernatant
4. Add 1 volume 80% EtOH (fresh). Incubate for  00:00:30
5. Remove supernatant
6. Repeat EtOH wash
7. Air dry for  00:10:00 -  00:15:00
8. Re-suspend beads thoroughly in  100 µL EB or TE buffer
9. Repeat cleanup (from step 1-7) and elute in  30 µL EB or TE buffer

29 Pooled library QC



Pooled cDNA library of 784 cells on HS D5000 Agilent tapestation