

Aug 29, 2024

## Purification of NIX-GST

DOI

[dx.doi.org/10.17504/protocols.io.q26g711dkgwz/v1](https://dx.doi.org/10.17504/protocols.io.q26g711dkgwz/v1)

Elias Adriaenssens<sup>1</sup>

<sup>1</sup>Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

---

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.q26g711dkgwz/v1](https://dx.doi.org/10.17504/protocols.io.q26g711dkgwz/v1)

**Protocol Citation:** Elias Adriaenssens 2024. Purification of NIX-GST. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.q26g711dkgwz/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 23, 2024

**Last Modified:** August 29, 2024

**Protocol Integer ID:** 101107

**Keywords:** ASAPCRN

**Funders Acknowledgement:**

**Aligning Science Across  
Parkinson's (ASAP)**

**Grant ID:** ASAP-000350

**Marie Skłodowska-Curie  
MSCA Postdoctoral  
fellowship**

**Grant ID:** 101062916

## Abstract

This protocol details the purification of NIX-GST.

## Materials

### Lysis buffer:

A	B
Tris-HCl	50 mM
pH	7.4
NaCl	300 mM
Triton X-100	1%
glycerol	5%
MgCl <sub>2</sub>	2 mM
DTT	1 mM
β-mercaptoethanol	2mM
cOmplete EDTA-free protease inhibitors (Roche)	
CIP protease inhibitor (Sigma)	
DNase (Sigma)	

### Wash buffer:

Tris-HCl	50 mM
pH	7.4
NaCl	300 mM
DTT	1 mM

### High salt wash buffer:


Tris-HCl	50 mM
pH	7.4
NaCl	700 mM
DTT	1 mM


### SEC buffer:


Tris-HCl	25 mM
pH	7.4
NaCl	300 mM

DTT	1 mM
-----	------

## Materials:

- pET-DUET1 vector (available from Addgene).  pETDuet-1 TIM9,10 **addgene Catalog #170280**
- NIX E72A/L75A/D77A/E81A (4A; ΔWIPI2) (available from Addgene)
- NIX W35A/L38A (ΔLIR) (available from Addgene).
- Rosetta pLysS cells (Novagen Cat# 70956-4)
 






 Rosetta™(DE3)pLysS Competent Cells - Novagen **Merck Catalog #70956-4**
- 10 kDa cut-off Amicon filter (Merck Millipore)
 

 Amicon® Ultra Centrifugal Filter, 10 kDa MWCO **Merck MilliporeSigma (Sigma-Aldrich) Catalog #UFC801008**



## Purification



16h

- 1 To purify NIX-GST, fuse the cytosol-exposed domain of NIX (1-182aa) to a C-terminal GST-tag through cloning into a pET-DUET1 vector (available from Addgene).
- 2 Introduce the point mutants in vitro mutagenesis to generate NIX E72A/L75A/D77A/E81A (4A;  $\Delta$ WIPI2) (available from Addgene), and NIX W35A/L38A ( $\Delta$ LIR) (available from Addgene).
- 3 After the transformation of the pET-DUET1 vector encoding NIX-GST wild-type or mutants in *E. coli* Rosetta pLysS cells (Novagen Cat# 70956-4), grow the cells in 2x Tryptone Yeast extract (TY) medium at  37 °C until an OD<sub>600</sub> of 0.4 and then continue at  18 °C .
- 4 Once the cells reaches an OD<sub>600</sub> of 0.8, induce the protein expression with  100 micromolar ( $\mu$ M) isopropyl  $\beta$ -D-1-thiogalactopyranoside (IPTG) for  16:00:00 at  18 °C .
- 5 Collect the cells centrifugation and resuspend in lysis buffer.

16h

Lysis buffer:

A	B
Tris-HCl pH 7.4	50 mM
NaCl	300 mM
Triton X-100	1%
Glycerol	5%
MgCl <sub>2</sub>	2 mM
DTT	1 mM
$\beta$ -mercaptoethanol	2mM
cOmplete EDTA-free protease inhibitors (Roche)	
CIP protease inhibitor (Sigma)	
DNase (Sigma)	

- 6 Sonicate the cell lysates twice for 30 s and clears by centrifugation at  18.000 rpm, 4°C, 00:45:00 in a SORVAL RC6+ centrifuge with an F21S-8x50Y rotor (Thermo Scientific).
- 6.1 Sonicate the cell lysates for  00:00:30 (1/2).

45m

30s



6.2 Sonicate the cell lysates for 00:00:30 (2/2).

30s



7 Collect the supernatant and incubate with pre-equilibrated Glutathione Sepharose 4B beads (GE Healthcare) for 02:00:00 at 4 °C with gentle shaking to bind NIX-GST.

2h



8 Centrifuge the samples to pellet the beads and remove the unbound lysate.



9 Wash the beads twice with wash buffer, once with high salt wash buffer, and two more times with wash buffer.

#### Wash buffer:

A	B
Tris-HCl pH 7.4	50 mM
NaCl	300 mM
DTT	1 mM

#### High salt wash buffer:

A	B
Tris-HCl pH 7.4	50 mM
NaCl	700 mM
DTT	1 mM

10 Incubate the beads Overnight with 4 mL of [M] 50 millimolar (mM) reduced glutathione dissolved in wash buffer at 4 °C, to elute NIX-GST from the beads.

2h



#### Wash buffer:

A	B
Tris-HCl pH 7.4	50 mM
NaCl	300 mM



A	B
DTT	1 mM

11 To collect the supernatant, collect the beads by centrifugation.



12 Wash the beads twice with  4 mL of wash buffer, and collect the supernatant.

13 Pool the supernatant fractions, filter through a 0.45 µm syringe filter, concentrate with 10 kDa cut-off Amicon filter (Merck Millipore), and load onto a pre-equilibrated Superdex 200 Increase 10/300 GL column (Cytiva). Elute the proteins with SEC buffer. Analyze the fractions by SDS-PAGE and Coomassie staining. Pool fractions containing purified NIX-GST.


SEC buffer:

A	B
Tris-HCl pH 7.4	25 mM
NaCl	300 mM
DTT	1 mM

14 After concentrating the purified protein, aliquote the protein and snap-frozen in liquid nitrogen.



#### Note

Store the proteins at  -80 °C .