


Jul 29, 2025

# NanoTag

DOI

[dx.doi.org/10.17504/protocols.io.kxygx34zkg8j/v1](https://dx.doi.org/10.17504/protocols.io.kxygx34zkg8j/v1)



Maria Dimitriu<sup>1</sup>, Leonard Steg<sup>1</sup>, Rodrigo Arzate-Mejia<sup>1</sup>, Isabelle Mansuy<sup>1</sup>

<sup>1</sup>ETH Zurich

Epigenomics



Maria Dimitriu

ETH Zurich

OPEN ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.kxygx34zkg8j/v1](https://dx.doi.org/10.17504/protocols.io.kxygx34zkg8j/v1)

**Protocol Citation:** Maria Dimitriu, Leonard Steg, Rodrigo Arzate-Mejia, Isabelle Mansuy 2025. NanoTag. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kxygx34zkg8j/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** December 13, 2023

**Last Modified:** July 29, 2025

**Protocol Integer ID:** 92256

**Keywords:** histone mark h3k4me3 genome, tagged protein, nanotag result, nanotag, ctcf in mouse embryonic stem cell, wide profiling of dna, chromatin, most current methods for genome, nanotag background, proteins such as chip, mouse embryonic stem cell, gene regulation, bound protein, binding protein, protein interactions in cell

**Funders Acknowledgements:**

ETH Zurich

Grant ID: ETH-10 15-2

## Abstract

**Background:** Genome-wide profiling of DNA-protein interactions in cells can provide important information about mechanisms of gene regulation. Most current methods for genome-wide profiling of DNA-bound proteins such as ChIP-seq and CUT&Tag use conventional IgG antibodies to bind the target protein(s). This limits their applicability to targets with available high affinity and specificity antibodies and prevents their use for other targets. Here we describe NanoTag, an IgG-free method derived from CUT&Tag to profile DNA-protein interactions. NanoTag is based on a fusion between an anti-GFP nanobody and Tn5 transposase that can map GFP-tagged proteins associated with chromatin in a fast, cost-effective and animal-free manner.

**Results:** We used NanoTag to indirectly profile the histone mark H3K4me3 genome-wide via its binding partner TATA box-binding protein-associated factor 3 (TAF3) and the transcription factors Nanog and CTCF in mouse embryonic stem cells (mESCs). NanoTag results show high inter-replicate reproducibility, high signal-to-noise ratio and strong correlation with CUT&Tag datasets, validating its accuracy and reliability.

**Conclusions:** NanoTag provides a novel, flexible and cost-effective IgG-free method to generate high resolution DNA-binding profiles in cells and tissues.

## Materials

This protocol requires having prepared GFP nanobody-Tn5 transposomes with a known concentration. For details on preparation of the protein and transposomes, see the last section of the protocol.

Materials required for the protein purification:

### HEGX buffer

20 mM HEPES (pH 7.5), 800 mM NaCl, 1 mM EDTA, 0.1% Triton X-100, 10% Glycerol, protease inhibitors

### Tn5 dialysis buffer

100 mM HEPES (pH 7.2), 200 mM NaCl, 0.2 mM EDTA, 2mMDTT, 0.2% (v/v) Triton X-100, 20% (v/v) glycerol

### Oligonucleotides

Tn5ME-A: TCGTCGGCAGCGTCAGATGTGTATAAGAGACAG

Tn5ME-B: GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAG

Tn5-ME-rev: /5Phos/CTGTCTCTTATACACATCT

Equipment needed for the purification of GFP nanobody-Tn5:

42°C water bath

Shaker incubator

5-liter Erlenmeyer flasks

500-ml centrifuge tubes

Refrigerated centrifuge with swinging-bucket rotor

Laboratory balance

Probe sonicator with 10-mm tip (e.g., Branson Sonifier)

15- and 50-ml conical centrifuge tubes

Econo-Pac Chromatography Columns (BioRad, cat. no. 7321010)

Membrane Dialysis, 12 to 14 kD (Spectrum Labs, cat. no. 132706)

PCR tubes

NanoTag buffers:

### Wash buffer

20 mM HEPES pH 7.5, 150 mM NaCl, 0.5 mM Spermidine, protease inhibitors

### Binding buffer

20 mM HEPES-KOH pH 7.9, 10 mM KCl, 1 mM  $\text{CaCl}_2$ , 1 mM  $\text{MnCl}_2$

### Dig-wash buffer

0.05% Digitonin, 2 mM EDTA in Wash buffer (nanobody-Tn5 will be diluted in this buffer and added to samples)

### Dig-med buffer

20 mM HEPES pH 7.5, 300 mM NaCl, 0.5 mM Spermidine, 0.01% Digitonin, protease inhibitors



### **Tagmentation buffer**

10 mM MgCl<sub>2</sub> in Dig-med buffer

### **Elution buffer**

1 mM Tris-HCl pH 8.0, 0.1 mM EDTA

Equipment needed for performing the NanoTag protocol and library preparation:

Centrifuge

Nutator

End-over-end rotator

Thermocycler

Lightcycler

Thermomixer/water bath



## Cell collection

17m

1 Remove media from culture dish and wash with D-PBS.



2m

2 Dissociate cells using Accutase and count cells.

10m

### Note

Use of trypsin to dissociate cells will disrupt the binding of the cells to the concanavalin A beads in subsequent steps, so use of Accutase is crucial.

3 Wash cells with  1 mL Wash buffer and resuspend cells at 1 million cells/mL of Wash buffer.  400  $\mu$ L per sample. Transfer in a 2 mL tube. Use 400,000 cells per NanoTag sample.



5m

4 Take  400  $\mu$ L of cell suspension (equivalent to 400k cells) per sample.


## Binding cells to beads

15m


5 Prepare concanavalin A beads:

5.1 Withdraw 10  $\mu$ L of bead slurry per sample (i.e. use  170  $\mu$ L for 16 sample), transfer it into  1.6 mL Binding buffer and mix by pipetting.




1m

5.2 Place tube on magnet until solution is clear. Remove the liquid, then remove tube from magnet and resuspend in  1.5 mL Binding buffer. Mix by pipetting.

2m











5.3 Place tube on magnet until solution is clear, then remove the liquid and resuspend the slurry in the original volume (  170  $\mu$ L for 16 samples) of Binding buffer. Keep the beads at room temperature until cells/nuclei are ready.



- 6 Add the correct amount of beads to the corresponding amount of cell suspension (  10  $\mu\text{L}$  of bead slurry and  400  $\mu\text{L}$  of cell suspension) per sample while gently vortexing.
- 7 Mix by pipetting and incubate on an end-over-end rotator for  00:15:00 . During the incubation, continue with step 8. 15m

## Tagmentation




1h 5m



- 8 During incubation of step 6 mix the correct amount of GFP nanobody-Tn5 transposomes (  2  $\mu\text{L}$  per sample for [M] 23 micromolar ( $\mu\text{M}$ ) protein) and Dig-wash buffer (  100  $\mu\text{L}$  per sample).
- 9 After the incubation from step 6 aliquot the bead-cell suspension to the corresponding amount of tubes.
- 10 Place tubes on the magnet until solution is clear. Remove the liquid, remove tubes from magnet and resuspend beads in  100  $\mu\text{L}$  of cold GFP nanobody-Tn5-containing Dig-wash buffer by gently pipetting.
- 11 Incubate on a nutator  Overnight at  4  $^{\circ}\text{C}$  to allow the GFP nanobody-Tn5 to bind to its target.
- 12 The next day place tubes on the magnet until solution is clear and remove the liquid. Wash beads with  1 mL of Dig-med buffer. Incubate on nutator for  00:05:00 . 5m
- 13 Repeat step 11 once to remove any anti-GFP nanobody-Tn5 protein that is unspecifically bound.
- 14 Resuspend beads in  300  $\mu\text{L}$  of Tagmentation buffer and incubate at  37  $^{\circ}\text{C}$  for  01:00:00 . 1h

## Tagmentation stop and DNA extraction


1h 14m




15 Add  10  $\mu\text{L}$  0.5M EDTA,  3  $\mu\text{L}$  10% SDS and  2.5  $\mu\text{L}$  of 20 mg/mL Proteinase K to each sample. Vortex tubes on full speed for 2-3 seconds.



16 Incubate at  55 °C for  01:00:00 .


1h

17 Add  300  $\mu\text{L}$  Phenol:Chloroform:Isoamyl Alcohol and full speed vortex for 2s.



18 Transfer solution to a phase-lock tube and centrifuge  16000 x g, 00:03:00 .

3m

19 Transfer the aqueous (top) layer to a fresh 1.5 mL tube containing  750  $\mu\text{L}$  100% ethanol, pipet up and down to mix and keep tubes  On ice .

20 Centrifuge  16000 x g, 4°C, 00:10:00 . Pour off the liquid and quickly drain on a paper towel.

10m

21 Add  1 mL 100% ethanol and centrifuge  16000 x g, 4°C, 00:01:00 . Pour off the liquid and let tubes dry on a paper towel for about 10 minutes.

1m

22 When tubes are dry add  25  $\mu\text{L}$  Elution buffer and full-speed vortex to dissolve the DNA.






23 Store tagmented DNA at  -20 °C .



II

## Library preparation

5m

24 To identify the optimal amount of cycles to perform during the library amplification, perform a qPCR:

24.1 Create a master mix by adding  2  $\mu\text{L}$  of i5 primer and  2  $\mu\text{L}$  of i7 primer,  0.15  $\mu\text{L}$  100X SYBR Green dye,  0.15  $\mu\text{L}$  water and  0.85  $\mu\text{L}$  NEBNext PCR Master mix per sample.

- 24.2 Add  12.5  $\mu\text{L}$  of master mix to each well of a qPCR plate. Add  2.5  $\mu\text{L}$  of tagged DNA to each well and run the qPCR using the following program:

Cycle 1: 72 °C for 5 min (gap filling)

Cycle 2: 98 °C for 30 sec

Cycle 3: 98 °C for 10 sec

Cycle 4: 63 °C for 30 sec


Cycle 5: 72 °C for 60 sec

Repeat Cycles 3-4-5 20 times

- 24.3 Check the amount of cycles necessary to achieve 1/3 of the maximum fluorescence for each sample and perform library amplification using that amount of cycles.

#### Note

The amount of cycles varies by target but usually between 18 and 21 cycles are necessary to ensure a minimum library concentration of 2 nM after clean-up.










20231221\_3edit.jpg

- 25 Once the amount of cycles has been determined according to step 23 amplify libraries:





- 25.1 Transfer  21  $\mu\text{L}$  of tagmented DNA from each sample into a PCR tube. Add  2  $\mu\text{L}$  of i5 index,  2  $\mu\text{L}$  of i7 index (use different i5 and i7 combinations for each library) and  25  $\mu\text{L}$  of NEBNext PCR master mix, mix and briefly spin down the sample then perform the following PCR program:
- Cycle 1: 72 °C for 5 min (gap filling)  
Cycle 2: 98 °C for 30 sec  
Cycle 3: 98 °C for 10 sec  
Cycle 4: 63 °C for 10 sec  
Repeat Cycles 3-4 for the amount of times determined in step 23  
72°C for 1 min and hold at 8 °C
- 26 Add  60  $\mu\text{L}$  AMPure beads to each library. Pipet up and down 10x to mix. Briefly spin down and incubate for 10 min. Place on magnet and allow solution to clear.
- 27 Remove the liquid and wash beads 2x with  200  $\mu\text{L}$  of 80% ethanol with the tubes still on the magnet. After the second wash remove the liquid using 2 pipetting steps.
- 28 Add  20  $\mu\text{L}$  Tris-HCl pH 8.0 and full-speed vortex. Incubate for 5 min. Place on magnet and allow solution to clear. Transfer liquid to fresh tubes.
- 29 Check the library quality using a fragment analyzer (e.g. Tapestation).

## Note




Example of Tapestation profile of a NanoTag library targeting TAF3 in mESC.

## Sequencing

- 30 NanoTag libraries can be sequenced using Illumina platforms. Before sequencing, serially dilute the libraries to equal molarity and pool to the desired molarity according to the sequencing platform specifications. We recommend adding 10% PhiX spike-in to the pool and performing paired-end sequencing using 50bp reads.

## Preparing GFP nanobody-Tn5 transposomes

- 31 Transform LysY/Iq competent *E. coli* with pMD1 plasmid encoding GFP nanobody-Tn5 with a intein-CBD tag.
- 32 Inoculate 2L of LB (+100ug/mL Ampicilin) with 1 colony and induce overexpression of GFP nanobody-Tn5 with 0.25 mM IPTG overnight at  18 °C .
- 33 The next day collect the bacterial medium, centrifuge and freeze the pellet.
- 34 Resuspend the pellet in 40 mL HEGX buffer and sonicate for 2.5 min using 10s cycles at 70% amplitude.



- 35 Centrifuge for 30 min at 20,000xg.
- 36 Collect the supernatant and add 1.1 mL of 10% PEI (pH 7.5) to 40 mL of supernatant.
- 37 Centrifuge for 20 min at 20,000xg.
- 38 Bind the supernatant to chitin beads for 2h at 4 °C in a disposable column. Wash thoroughly with 10 column volumes of HEGX buffer.
- 39 Induce cleavage of the GFP nanobody-Tn5 from its intein tag by incubating the column overnight at 4 °C in HEGX buffer + 100 mM DTT.
- 40 Collect the eluate. Prepare the oligonucleotides and add the oligonucleotides to the eluate and incubate for 1h at 30 °C .
- 40.1 Resuspend oligonucleotides at 200 uM in water.
- 40.2 Incubate for 5 min at 95 °C in a thermocycler and cool down over
- 41 Dialyze overnight in Tn5 dialysis buffer.
- 42 Concentrate the protein using 30kDa Amicon centrifugation filter. Store the concentrated protein in Tn5 storage buffer, prepare aliquots, flash freeze in liquid nitrogen and store at -80 °C .