

Aug 15, 2024

Immunohistochemical staining of wholemount major pelvic ganglia (MPG) for analysis of myelinated bladder afferents

Forked from a private protocol

DOI

dx.doi.org/10.17504/protocols.io.yxmvme3nbg3p/v1

Janet R Keast¹, Peregrine B Osborne¹, Nicole Wiedmann¹

¹University of Melbourne

SPARC

Tech. support email: info@neuinfo.org



John-Paul Fuller-Jackson

University of Melbourne

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.yxmvme3nbg3p/v1

Protocol Citation: Janet R Keast, Peregrine B Osborne, Nicole Wiedmann 2024. Immunohistochemical staining of wholemount major pelvic ganglia (MPG) for analysis of myelinated bladder afferents. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.yxmvme3nbg3p/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 23, 2024

Last Modified: August 15, 2024

Protocol Integer ID: 100336

**Funders Acknowledgement:****NIH SPARC****Grant ID: 3OT2OD023872**

Abstract

This protocol describes immunohistochemical procedures applied to wholemount major pelvic ganglia (MPG) for the visualization of myelinated axons filled with cholera toxin subunit B (CTB) and/or neurons transfected with adeno-associated virus (AAV) expressing TdTomato. In this protocol, samples were obtained from rats in which CTB was microinjected into the bladder, and AAV-PHP.S was intravenously administered to preferentially and sparsely label peripheral neurons. In this context, neurofascin (paranode marker) identifies myelinated axons.

Materials

- ⊗ Horse serum **Sigma Aldrich Catalog #12449C**
- ⊗ Triton X-100 **Sigma Aldrich Catalog #T8787-50ML**
- ⊗ Anti-RFP antibody (guinea pig) **Synaptic Systems Catalog #390004**
- ⊗ Anti-cholera toxin subunit B antibody (goat) **List Labs Catalog #703**
- ⊗ Anti-neurofascin antibody (rabbit) **Alomone Labs Catalog #AIP-025**
- ⊗ Cy3 Donkey anti-guinea pig IgG **Jackson ImmunoResearch Laboratories, Inc. Catalog #706-165-148**
- ⊗ AF488 Donkey anti-goat IgG **Jackson ImmunoResearch Laboratories, Inc. Catalog #705-545-147**
- ⊗ AF647 Donkey anti-rabbit antibody **Invitrogen Catalog #A32795**

Solutions:

- PBS: phosphate-buffered saline, 0.1 M, pH 7.2
- PBS containing 0.1% sodium azide
- PB: phosphate-buffer, 0.1M, pH7.2
- Blocking solution: PBS containing 10% normal horse serum and 0.5% triton X-100
- PBS containing 0.1% sodium azide, 2% normal horse serum and 0.5% triton X-100

Primary Antibodies:

A	B	C	D	E
Abbreviation	Synonym	RRID	Host Species	Dilution
RFP	Red fluorescent protein	AB_2737052	Guinea pig	1:500
CTB	Cholera toxin subunit B	AB_10013220	Goat	1:10,000
Neurofascin	NF155	AB_2756657	Rabbit	1:1000

Secondary Antibodies:

A	B	C	D
Tag-antibody	Host Species	RRID	Dilution
Cy3 anti-guinea pig	Donkey	AB_2340460	1:2000
AF488 anti-goat	Donkey	AB_2336933	1:1000
AF647 anti-rabbit	Donkey	AB_2762835	1:1000



Immunohistochemistry

- 1 Wash whole MPGs in phosphate buffer (PB; 0.1M; pH 7.2) (3 x 30 min)
- 2 Incubate sections in blocking solution (PB; 10% horse serum; and 0.5% Triton-X) at room temperature for 2 h
- 3 Incubate sections in appropriate dilutions of primary antibodies (or combinations of primary antibodies) for 72h. Antibodies are diluted in PBS containing 0.1% sodium azide, 2% horse serum, and 0.5% triton-X.
- 4 Wash tissue in PBS (3 x 30 min)
- 5 Incubate sections in appropriate dilutions of secondary antibodies (or combinations of secondary antibodies) 24 h in the dark. Antibodies are diluted in PBS containing 2% horse serum, and 0.5% triton-X.
- 6 Wash tissue in PBS (3 x 30 min)
- 7 Mount tissue onto glass slides and coverslip in preferred anti-fade mountant.