

May 10, 2024

RNA extraction from milk for HPAI surveillance

DOI

dx.doi.org/10.17504/protocols.io.n2bvjn6obgk5/v1

Andrew Lail¹, Will Vuyk¹, Isla Emmen¹, David O'Connor¹

¹University of Wisconsin-Madison



Andrew Lail

University of Wisconsin-Madison

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.n2bvjn6obgk5/v1

Protocol Citation: Andrew Lail, Will Vuyk, Isla Emmen, David O'Connor 2024. RNA extraction from milk for HPAI surveillance. protocols.io <https://dx.doi.org/10.17504/protocols.io.n2bvjn6obgk5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 01, 2024

Last Modified: May 10, 2024

Protocol Integer ID: 99087

Keywords: Dairy, Milk, Influenza, Viral RNA, RNA Extraction



Abstract

This is an RNA extraction method that isolates HPAI RNA segments from commercially-available milk samples. The protocol involves processing with a MagMax RNA extraction kit on a Kingfisher instrument. Then, the samples are concentrated using a Zymo Clean and Concentrate kit. The total RNA yields should be suitable for performing qPCR. After trying other methods involving centrifugation, filtration, and trizol, we determined this method to be the easiest and most effective.

Materials

MagMax Wastewater Ultra x96 kit

Zymo RNA Clean and Concentrator-5 kit

ThermoFisher TurboDNase kit

Heat block

Lab plasticware

Safety warnings

- ⚠ Performing this protocol with unpasteurized milk from an infected cow might expose the user to active flu virus. Please use pasteurized samples or take appropriate precautions.



RNA extraction

- 1 Use the MagMax Wastewater Ultra x96 kit¹ (or comparable magnetic bead isolation kit²⁻⁵) on the ThermoFisher Kingfisher Apex (or comparable instrument). For each sample, load 400 uL of pasteurized milk or heavy whipping cream into the binding plate.
- 2 Elute in 50 uL of elution solution or nuclease-free water.

Optional Sample Concentration and Cleaning

- 3 Pool the eluate from 6 extraction replicates of the same sample.
- 4 Take 200 uL forward into the DNase treatment. Using ThermoFisher TurboDNase, combine 200 uL of pooled eluate, 22 uL 10X TurboDNase buffer, and 2 uL TurboDNase.
- 5 Incubate at 37 C for 30 minutes.
- 6 Use a Zymo RNA Clean and Concentrator-5 kit⁷ to clean and concentrate the sample. Elute this into 30 uL final volume of nuclease-free water for downstream analysis.

Protocol references

1. MagMAXTM Wastewater Ultra Nucleic Acid Isolation Kit [Internet]. [cited 2024 May 10]. Available from: <https://www.thermofisher.com/order/catalog/product/A52606>
2. Guidelines and minimum requirements for diagnosis of H5N1 HPAIV infection in cattle [Internet]. Istituto Zooprofilattico Sperimentale delle Venezie; 2024 Apr [cited 2024 May 10]. Report No.: Version 2024.04.10. Available from: <https://www.izsvenezie.com/documents//reference-laboratories/avian-influenza/diagnostic-protocols/guidelines-diagnosis-h5n1-cattle.pdf>
3. MagMAXTM-96 Viral RNA Isolation Kit [Internet]. [cited 2024 May 10]. Available from: <https://www.thermofisher.com/order/catalog/product/AM1836>
4. MagMAXTM Viral/Pathogen Ultra Nucleic Acid Isolation Kit [Internet]. [cited 2024 May 10]. Available from: <https://www.thermofisher.com/order/catalog/product/A42356>
5. MagMAXTM Microbiome Ultra Nucleic Acid Isolation Kit, with bead tubes [Internet]. [cited 2024 May 10]. Available from: <https://www.thermofisher.com/order/catalog/product/A42358>
6. TURBOTM DNase (2 U/ μ L) [Internet]. [cited 2024 May 10]. Available from: <https://www.thermofisher.com/order/catalog/product/AM2238>
7. ZYMO RESEARCH [Internet]. [cited 2024 May 10]. RNA Clean & Concentrator-5. Available from: <https://www.zymoresearch.com/products/rna-clean-concentrator-5>