

VERSION 2
NOV 06, 2023

OPEN ACCESS



DOI:
dx.doi.org/10.17504/protocols.io.81wgb6z9olpk/v2

Protocol Citation: Sophia Ferchiou, France Caza, Yves St-Pierre 2023. Liquid biopsy in sentinel mussels .


protocols.io
<https://dx.doi.org/10.17504/protocols.io.81wgb6z9olpk/v2>
ersion created by [Sophia Ferchiou](#)

MANUSCRIPT CITATION:
Caza, F. et al. Liquid biopsies for omics-based analysis in sentinel mussels. PLoS One. 14, e0223525 (2019)

🌐 Liquid biopsy in sentinel mussels V.2

Sophia Ferchiou¹, France Caza¹, Yves St-Pierre¹

¹Institut National de la Recherche Scientifique, Université du Québec



Sophia Ferchiou

ABSTRACT

This protocol has been optimized for sampling hemolymph in sentinel mussels in remote areas, such as polar regions. A detailed protocol can be found in [Caza *et al.* \(2019\)](#).

MATERIALS

A. Sampling Section

- 1) Syringe (from 3 to 5 cc)
- 2) 25 gauge needle
- 3) Sterile 1.5 mL Eppendorf tubes
- 4) Multi Spin Battery-Powered Mini-Centrifuge (3000 x *g*) or equivalent.
- 5) Whatman 903™ FTA cards
- 6) Plastic ziplock bags
- 7) Silica gel desiccants (1gr/bag)
- 8) Pipette P200
- 9) Sterilized filter tips
- 10) Knife
- 11) Ethanol 70%

B. Extraction of circulating cell-free DNA Section

- 1) QIAamp DNA Investigator Kit (Qiagen - Cat #56504)
- 2) Ethanol 100%
- 3) Sterile 1.5 mL Eppendorf tubes
- 4) Pipettes (P1000, P100, P20, P2)
- 5) Sterilized filter tips

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Nov 04, 2023

Last Modified: Nov 06, 2023

PROTOCOL integer ID: 90415

Keywords: Mussels, liquid biopsy, remote area, circulating cell-free DNA

BEFORE START INSTRUCTIONS

A. Sampling Section (Protocol optimized for freshly collected mussels)

Be sure to gather all materials before collecting mussels.

B. Extraction of circulating cell-free DNA Section

- 1) Equilibrate samples to room temperature (15–25°C).
- 2) Equilibrate Buffer ATE to room temperature.
- 3) Set a thermomixer or heated orbital incubator to 56°C for use in step #13.
- 4) Ensure that buffers AW1 and AW2 have been prepared :
 - a) Add 25 mL ethanol (100%) to the bottle containing 19 mL Buffer AW1 concentrate.
 - b) Add 30 mL ethanol (100%) to the bottle containing 13 mL Buffer AW2 concentrate.
- 5) If Buffer AL or Buffer ATL contains precipitates, dissolve by heating to 70°C with gentle agitation.
- 6) Reconstitute carrier RNA by adding 310 µL of buffer ATE.

Sampling		50m
1	Collect mussels (<i>Mytilus spp.</i>) with a shell length range between 50 to 70 mm.	10m
2	Remove intervalvar liquid by opening gently the valves with the tip of a knife.	5m
3	Withdraw hemolymph from the adductor muscle using a syringe fitted with a 25G needle.	10m
4	Transfer immediately the hemolymph to a sterile 1.5 mL Eppendorf tube (Optional: You can pool 3 or 4 hemolymph samples to obtain a total volume of 1.5 mL per tube).	1m

- 5 Centrifuge for 3 minutes at maximum speed (approx. 3000 $\times g$) at room temperature using a battery-powered mini-centrifuge (TOMY, Japan) or equivalent. 3m
- 6 After centrifugation, transfer the supernatant (clear solution above pellet at bottom of tube) in a new 1.5 mL Eppendorf tube and gently resuspend the pellet. The pellets can be frozen (-20°C) or spotted on FTA cards until used. Supernatants can also be frozen (-20°C) or spotted on FTA cards (using 70 μ L aliquots) until use. We routinely use Whatman 903™ FTA cards with 13 mm discs to spot pellets or supernatants. 5m
- Note :** Once spotted, FTA cards are dried for 15 minutes at room temperature and stored individually in small ziplock plastic bags containing one silica gel desiccant moisture absorber.
- Note :** Unless otherwise indicated, hemolymph should be collected and processed within one hour after sampling.
- 7 Bring the mussels back to their respective mussel bed. 5m

Extraction of circulating cell-free DNA

1h 13m

- 8 Unless otherwise indicated, you can use the [QIAamp DNA Investigator Kit](#) (Qiagen) for extraction of the hemolymphatic circulating cell-free DNA.
- 9 First, thaw frozen supernatants at room temperature and clarify by centrifugation for 10 minutes at 4500 $\times g$ at room temperature. 10m
- 10 Transfer 70 μ L of the clarified supernatant into a new sterile 1.5 mL Eppendorf tube. 2m
- 11 Add 30 μ L of buffer ATL and 10 μ L of proteinase K (>600 mAU/mL) to the sample. 5m

- 12** Add 100 μL of buffer AL and 1 μL of reconstituted RNA Carrier (1 $\mu\text{g}/\mu\text{L}$). Mix by pulse-vortexing for 15 seconds or until it is a homogeneous solution. 2m
- 13** Incubate at 56°C for 10 minutes on a lab orbital shaker set at 900 rpm. 10m
- 14** Add 50 μL of ethanol (100%) and mix thoroughly by pulse-vortexing for 15 seconds. Then, incubate for 3 minutes at room temperature. 5m
- 15** Transfer the entire lysate from step #14 to the QIAamp MinElute column. 5m
- 16** Centrifuge at 6000 $\times g$ for 1 minute. Place the QIAamp MinElute column in a clean 2 mL collection tube and discard the collection tube containing the flow-through. 1m
- 17** Add 500 μL of buffer AW1 to the QIAamp MinElute column. 1m
- 18** Repeat Step #16. 1m
- 19** Add 700 μL of buffer AW2 to the QIAamp MinElute column. 1m
- 20** Repeat Step #16. 1m

- 21** Add 700 μL of ethanol (100%) to the QIAamp MinElute column. 1m
- 22** Repeat Step #16. 1m
- 23** Centrifuge at full speed (20,000 x g) for 3 minutes at room temperature. 3m
- 24** Place the QIAamp MinElute column in a clean 1.5 mL Eppendorf tube and discard the collection tube containing the flow-through. 1m
- 25** Open the lid and incubate at room temperature for 10 minutes. 10m
- 26** Apply 60 μL of buffer ATE to the center of the membrane. 1m
- 27** Close the lid and incubate at room temperature for 10 minutes. 10m
- 28** Centrifuge at full speed (20,000 x g) for 1 minute. 1m

29 Discard the QIAamp MinElute column and keep the eluate.

1m