



Dec 07, 2020

Protein Concentration Determination using Qubit

In 1 collection

Steven J Burgess¹¹University of Illinois at Urbana-Champaign

In Development

This protocol is published without a DOI.

Ag SynBio Lab UIUC



Steven Burgess

University of Illinois at Urbana-Champaign

ABSTRACT

Procedure for quantification of protein concentration. This version is adapted for use with samples that have been extracted in protein extraction buffer as part of processing soybean/cowpea samples, but the general outline is similar for other samples. See the manual for further details: https://assets.thermofisher.com/TFS-Assets/LSG/manuals/Qubit_Protein_Assay_UG.pdf

PROTOCOL CITATION

Steven J Burgess 2020. Protein Concentration Determination using Qubit. **protocols.io**
<https://protocols.io/view/protein-concentration-determination-using-qubit-bqhkm4w>

COLLECTIONS ⓘ

 **Immunoblot Analysis of Leaf Tissue**

LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 07, 2020

LAST MODIFIED

Dec 07, 2020

PROTOCOL INTEGER ID

45324

PARENT PROTOCOLS

Part of collection

[Immunoblot Analysis of Leaf Tissue](#)

MATERIALS TEXT

- Qubit® protein assay kit (Life Technologies; [Q33211](#))
- Qubit® assay tubes (Life Technologies; [Q32856](#))
- Qubit® 4 Fluorometer (Life Technologies; [Q33238](#))

ABSTRACT

Procedure for quantification of protein concentration. This version is adapted for use with samples that have been extracted in protein extraction buffer as part of processing soybean/cowpea samples, but the general outline is similar for other samples. See the manual for further details: https://assets.thermofisher.com/TFS-Assets/LSG/manuals/Qubit_Protein_Assay_UG.pdf

Create Working Solution for Analysis

Create a working solution of Qubit assay buffer by diluting the reagent 1:200 in the provided buffer.

1


The final volume in each tube must be 200 μL . Each standard tube requires 190 μL of Qubit® working solution, and each sample tube requires anywhere from 180–199 μL . Prepare sufficient Qubit® working solution to accommodate all standards and samples.

Create Sample Dilution for Analysis

15m 6s

2 Dilute sample 1:50 ( 196 μL dH_2O +  4 μL sample)

After taking into consideration the Qubit dilution factor (1:20; 10 μL sample + 190 μL qubit assay buffer) the sample being analyzed has been diluted 1:1000. This should give a value in the linear range for qubit (1.25-25 $\mu\text{g}/\text{mL}$). This should also reduce the impact of SDS on the quantification, which must be >0.2% in the 10 μL sample added (conc SDS in PEB is 2 %, so diluted 1:50 yields 0.04 %)

3 Add  190 μL of Qubit working solution to a fresh Qubit assay tube, one for each sample to be analyzed (including the three protein standards)4 Add  10 μL of protein standard to the appropriate tube and mix by vortexing  00:00:03 s.

3s



5 Add  10 μL of diluted sample to the appropriate tube and mix by vortexing  00:00:03 s.

3s

6 Allow samples to incubate at room temperature for  00:15:00

15m

7 Measure protein sample concentration using the Qubit, following the instructions on the machine (i.e. start by measuring the standard curve).

Remember to adjust measured values according to the dilution factor applied. In this example, after taking into consideration the Qubit dilution factor (1:20;  10 μL sample +  190 μL Qubit assay buffer) the sample has been diluted 1:1000

Qubit™ 4 Fluorometer, with WiFi
Fluorometer

Invitrogen

Q33238

