



OCT 16, 2023

Expression and purification of PI3KC3-C1 complex

Minghao Chen¹, Xuefeng Ren¹

¹Department of Molecular and Cell Biology, University of California, Berkeley, Berkeley, CA 94720, USA

ASAP Collaborative Research Network



Minghao Chen

OPEN  ACCESS



ABSTRACT

Expression and purification of PI3KC3-C1 complex

DOI:

dx.doi.org/10.17504/protocols.io.5qpvorjezv4o/v1

Protocol Citation: Minghao Chen, Xuefeng Ren 2023. Expression and purification of PI3KC3-C1 complex.

protocols.io

<https://dx.doi.org/10.17504/protocols.io.5qpvorjezv4o/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working







Created: Jun 05, 2023


Last Modified: Oct 16, 2023

PROTOCOL integer ID:
82916

Keywords: ASAPCRN

Funders
Acknowledgement:
Aligning Science Across
Parkinson’s (ASAP) initiative
Grant ID: ASAP-000350

| Expression | | 2d 0h 40m |
|------------|---|-----------|
| 1 | Transfect HEK GNTi cells at concentration of 2×10^6 cells/ml | |
| 2 | Dilute PEI with Warm Hybridoma-SFM(1X) | |
| 3 | In a separate tube, dilute DNA with Hybridoma-SFM(1X) | |
| 4 | Add PEI to DNA dilution. Incubate mixture for  00:30:00 at  37 °C | 30m |
| 5 | Add mixture to cells. Let cells grow for  48:00:00 | 2d |
| 6 | Harvest Cells  500 rpm ,  4 °C ,  00:10:00 | 10m |

7 Wash pellet with cold PBS. Store pellet at  -80 °C until purification.



Purification

2h 50m

8 Pellets were homogenized 20 times by Pyrex douncer (Corning) in lysis buffer (25 mM HEPES pH 7.5, 200 mM NaCl, 2 mM MgCl₂, 10% Glycerol) with 25mM TCEP/proteinase inhibitors (Thermo Scientific), and then add 10% Triton X-100 stock to final 1% concentration.



9 Rocking at  4 °C for  01:00:00

1h

10 Clarify lysate for  17000 rpm for  00:40:00 at  4 °C .

40m


11 Wash strep-tactin resin (IBA Lifesciences, Germany) into lysis buffer (without Triton). Load clarified lysate onto resin

12 Rock supernatant with equilibrated Strep-Tactin Sepharose resin for  Overnight at  4 °C

1h

13 Wash with 5CV lysis buffer (25 mM HEPES pH 7.5, 200 mM NaCl, 2 mM MgCl₂, 25 mM TCEP, 10% Glycerol)

14 Elute with lysis buffer plus 4 mM desthiobiotin for STREP resin

- 15** Concentrate elution and inject onto pre-equilibrated Superose 6 Increase 10/300 GL column (Cytiva)
(25 mM HEPES pH 7.5, 300 mM NaCl, 1 mM MgCl₂, 1 mM TCEP)
- 16** Pool peak fractions, concentrate, snap freeze, and store at  -80 °C