



BG11 hypersaline medium V.2

Tanja Bosak¹

¹Massachusetts Institute of Technology

Version 2

Aug 27, 2020

1

Works for me

dx.doi.org/10.17504/protocols.io.bkcmksu6

Bosak Lab



Emilie J Skoog

Massachusetts Institute of Technology

DOI

dx.doi.org/10.17504/protocols.io.bkcmksu6

PROTOCOL CITATION

Tanja Bosak 2020. BG11 hypersaline medium. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bkcmksu6>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 26, 2020

LAST MODIFIED

Aug 27, 2020

PROTOCOL INTEGER ID

41069

- 1 Add the following to a beaker and dissolve using a stir bar on a magnetic stir plate.

BG11 hypersaline medium	amount	unit
initial Milli-Q water	968	ml
NaCl	49.8	g
NaNO3	1.5	g
Na2CO3	0.02	g
KCl	1.3867	g
MgCl2*6H2O	8.569	g
MgSO4*7H2O	6.4998	g
CaCl2*2H2O	2.042	g
Stock A	10	ml
Stock B	10	ml
Stock C	10	ml
Stock 5 (trace metals)	1	ml
Vitamins	1	ml

- 2 Add Milli-Q water to bring total volume to 1 L.

3 Adjust pH to ~7 - 7.4

4 Filter sterilize near Bunsen burner or in sterile biohood.

5 STOCK SOLUTIONS

Stock A	Concentration (g/L)
Na2MG EDTA	0.1
Ferric ammonium citrate	0.6
Citric acid*1H2O	0.6
CaCl2*2H2O	3.6
Add up to 1L of dH2O	

Stock B	Concentration (g/L)
MgSO4*7H2O	7.5
Add up to 1L of dH2O	

Stock C	Concentration (g/L)
K2HPO4*3H2O	4
Add up to 1L of dH2O	

Stock 5	Concentration (g/L)
H3BO3	2.86
MnCl2*4H2O	1.81
ZnSO4*7H2O	0.222
CuSO4*5H2O	0.079
CoCl2*6H2O	0.05
NaMoO4*2H2O	0.391
Add up to 1L of dH2O	