



Oct 16, 2021

# Transformation of *E. coli* with plasmid

sanjana<sup>1</sup><sup>1</sup>3rd year BS-MS, IISER Pune

1



dx.doi.org/10.17504/protocols.io.byz3px8n

iGEM IISER Pune India 2021



sanjana

This protocol can be used to transform electrocompetent *E. coli* cells with a plasmid.

DOI

dx.doi.org/10.17504/protocols.io.byz3px8n

sanjana 2021. Transformation of *E. coli* with plasmid. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.byz3px8n>



\_\_\_\_\_ protocol ,

Oct 13, 2021

Oct 16, 2021

54043

## Preparation

30m

- 1 Wash cuvette with **100 % (v/v)** ethanol and UV in the hood for **00:20:00**

20m

- 2 Pre-thaw electrocompetent cells and plasmid on wet ice for **00:10:00**

10m

## Electroporation



20m 40s

- 3 Add **1 µL** of plasmid DNA to **50 µL** electrocompetent cells

4 Transfer cell-DNA mixture into the cuvette and add equal amount (  **50 µL** ) of MilliQ

5 Let the cuvette sit in wet ice for  **00:20:00** 20m



6 Place cuvette in BioRad pulser and pulse at 2.5 keV

7 Within  **00:00:40** add  **100 µL** of 2X LB to the cuvette 40s

Incubation 30m

8 Place cuvette in an incubator at  **37 °C** for  **00:30:00** 30m

9 Then use a bent pipette tip to spread the cuvette's contents onto a plate made with LB agar and the appropriate antibiotics. Spread until plate is dry

10 Place the plate in an incubator at  **37 °C**  **Overnight**

11 Make glycerol stocks of transformed cells and store for future use