



Aug 14, 2020

# Universal sandwich ELISA for investigating the binding of Protein-LA (SpLA) to avian immunoglobulins using a peroxidase-labeled -anti-IgY conjugate.

Angel A Justiz-Vaillant<sup>1</sup>, Monica F. Smikle<sup>2</sup>

<sup>1</sup>University of the West Indies St. Augustine; <sup>2</sup>University of the West Indies. Mona Campus

**1** Works for me [dx.doi.org/10.17504/protocols.io.bjq2kmye](https://dx.doi.org/10.17504/protocols.io.bjq2kmye)

University of the West Indies [angel.vaillant@sta.uwi.edu](mailto:angel.vaillant@sta.uwi.edu)

Angel Justiz-Vaillant  
University of the West Indies St. Augustine

DOI

[dx.doi.org/10.17504/protocols.io.bjq2kmye](https://dx.doi.org/10.17504/protocols.io.bjq2kmye)

PROTOCOL CITATION

Angel A Justiz-Vaillant, Monica F. Smikle 2020. Universal sandwich ELISA for investigating the binding of Protein-LA (SpLA) to avian immunoglobulins using a peroxidase-labeled -anti-IgY conjugate.. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bjq2kmye>

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 14, 2020

LAST MODIFIED

Aug 14, 2020

PROTOCOL INTEGER ID

40442

MATERIALS

NAME	CATALOG #	VENDOR
Anti-Chicken IgY, HRP Conjugate, 300ul	G1351	Promega
Nunc&trade; 96-Well Polystyrene Round Bottom Microwell Plates, V 96 well plate, Non-Treated, clear, without lid, Sterile	260210	Thermo Fisher
Staphylococcal Protein-A		Sigma Aldrich
Protein-L from P. Magnus		

- 1 This ELISA is used to study the interaction of protein-LA (SpLA) with diverse avian immunoglobulins.
- 2 The 96 well microtitre plate is coated overnight at 4°C with 2 µg/µl per well of recombinant SpLA or a mixture of SpA with SpL in carbonate-bicarbonate buffer pH 9.6.

- 3 Then plate is treated with bovine serum albumin solution and washed 4X with PBS-Tween.
- 4 50 µl of avian egg yolk or egg white (1 mg/ml) is added and incubated for 1.30h at RT and the microplate is then rewashed 4X with PBS-Tween.
- 5 Then 50 µl of peroxidase-labeled-anti-IgY conjugate diluted 1:15000 in PBS-non-fat milk is added to each well and incubated for 1.30h at RT. After that the plate is washed 4X with PBS-Tween.
- 6 Pipette 50 µl of 3,3',5,5' - tetramethylbenzidine (TMB; Sigma-Aldrich) to each well.
- 7 The reaction is stopped with 50 µl of 3M H<sub>2</sub>SO<sub>4</sub> solution.
- 8 The plate is visually assessed for the development of colour and read in a microplate reader at 450 nm.
- 9 A cut-off point can be calculated as the mean of the optical density of negative controls x 3. The higher the OD value the higher will be the binding affinity of SpLA to avian immunoglobulins.