



Aug 12, 2021

## SaCas9 protein purification

Fatma Gomaa<sup>1</sup>, Zhu-Hong Li<sup>2</sup>, Roberto Docampo<sup>3</sup>, Virginia Edgcomb<sup>4</sup><sup>1</sup>Harvard University; <sup>2</sup>The University of Georgia; <sup>3</sup>University of Georgia; <sup>4</sup>Woods Hole Oceanographic Institution

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.bxbupinw](https://dx.doi.org/10.17504/protocols.io.bxbupinw)

### Protist Research to Optimize Tools in Genetics (PROT-G)

Fatma Gomaa  
Harvard University

#### ABSTRACT

Developing transfection protocol for *Bodo saltans*, using SaCas9/sgRNA ribonucleoprotein (RNP) complex in conjunction with DNA repair template to disrupt the Paraflagellar rod 2 gene (*BsPFR2*) and increase the efficiency of targeted homologous recombination when a repair template DNA is provided. The exogenous repair template is double stranded DNA and it consists of *eGFP* fused with the drug selection gene *nptII/neo* and flanked by 500 bp of the untranslated regions (UTRs) upstream and downstream of the targeted *BsPFR2* as homologous repair arms.

#### DOI

[dx.doi.org/10.17504/protocols.io.bxbupinw](https://dx.doi.org/10.17504/protocols.io.bxbupinw)

#### PROTOCOL CITATION

Fatma Gomaa, Zhu-Hong Li, Roberto Docampo, Virginia Edgcomb 2021. SaCas9 protein purification.  
**protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bxbupinw>

#### KEYWORDS

Bodo saltans transfection, SaCas9/RNP complex, homologous recombination

#### LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

#### CREATED

Aug 12, 2021

#### LAST MODIFIED

Aug 12, 2021

#### PROTOCOL INTEGER ID

52308

*Protocol 1: for SaCas9 protein purification*

SaCas9 protein was prepared following methods described in Medeiros et al. (Medeiros et al., 2017) with modifications. Briefly, the bacterial expression vector p6XHis-NLS-SaCas9 (Addgene #101086) was transformed into *Escherichia coli* Rosetta 2(DE3) competent cells (Novagen) and grown as an overnight preculture with shaking at 37°C. Subsequently, 1 ml of preculture was used to inoculate 100 ml of terrific broth media which was incubated with shaking at 37°C until it reached an optical density at of 0.6 at 600 nm. Protein expression was induced by addition of isopropyl-β-D-thiogalactopyranoside (IPTG) to a final concentration of 200 μM and cells were grown overnight with vigorous shaking at 18°C. Following incubation, cells were collected by centrifugation, and the pellet was resuspended in 2 ml / 100 mg cell weight of xTractor Buffer (Takara Bio) plus 5 mg of lysozyme. Cells were lysed by sonication, then 1 μl / 2ml of DNase I was added. The soluble fraction was purified using the His60 Ni Superflow Resin & Gravity Columns (Takara Bio) following the manufacturer's recommendations. The concentration of SaCas9 containing fractions were determined by optical density at 280 nm and then pooled. Pooled fractions were concentrated and desalted into Cytomix (120 mM KCl, 0.15 mM CaCl<sub>2</sub>, 10mM KH<sub>2</sub>PO<sub>4</sub>, 25 mM HEPES, 2 mM EGTA, 5 mM MgCl<sub>2</sub>, pH 7.6) using Amicon Ultra-15 centrifugal filter devices (Millipore Sigma). The protein concentration was determined using a BCA Protein Assay Kit (Pierce).

## SaCas9 protein purification

1 *Protocol 1: for SaCas9 protein purification*

SaCas9 protein was prepared following methods described in Medeiros et al. (Medeiros et al., 2017) with modifications. Briefly, the bacterial expression vector p6XHis-NLS-SaCas9 (Addgene #101086) was transformed into *Escherichia coli* Rosetta 2(DE3) competent cells (Novagen) and grown as an overnight preculture with shaking at 37°C. Subsequently, 1 ml of preculture was used to inoculate 100 ml of terrific broth media which was incubated with shaking at 37°C until it reached an optical density at of 0.6 at 600 nm. Protein expression was induced by addition of isopropyl-β-D-thiogalactopyranoside (IPTG) to a final concentration of 200 μM and cells were grown overnight with vigorous shaking at 18°C. Following incubation, cells were collected by centrifugation, and the pellet was resuspended in 2 ml / 100 mg cell weight of xTractor Buffer (Takara Bio) plus 5 mg of lysozyme. Cells were lysed by sonication, then 1 μl / 2ml of DNase I was added. The soluble fraction was purified using the His60 Ni Superflow Resin & Gravity Columns (Takara Bio) following the manufacturer's recommendations. The concentration of SaCas9 containing fractions were determined by optical density at 280 nm and then pooled. Pooled fractions were concentrated and desalted into Cytomix (120 mM KCl, 0.15 mM CaCl<sub>2</sub>, 10mM KH<sub>2</sub>PO<sub>4</sub>, 25 mM HEPES, 2 mM EGTA, 5 mM MgCl<sub>2</sub>, pH 7.6) using Amicon Ultra-15 centrifugal filter devices (Millipore Sigma). The protein concentration was determined using a BCA Protein Assay Kit (Pierce).