



Jan 07, 2021

Measurement of heat production (thermogenesis) in cells using ERthermAC dye_microplate method

Chih-Hao Wang¹, Yu-Hua Tseng¹, Yu-Hua Tseng²¹Graduate Institute of Biomedical Sciences, China Medical University, Taichung 40402, Taiwan;²Section on Integrative Physiology and Metabolism, Joslin Diabetes Center, Harvard Medical School, Boston, MA, USA**1** Works for me dx.doi.org/10.17504/protocols.io.bqprvm6

Tseng Lab, Joslin Diabetes Center

Chih-Hao Wang

DOI

dx.doi.org/10.17504/protocols.io.bqprvm6

PROTOCOL CITATION

Chih-Hao Wang, Yu-Hua Tseng, Yu-Hua Tseng 2021. Measurement of heat production (thermogenesis) in cells using ERthermAC dye_microplate method. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bqprvm6>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 14, 2020

LAST MODIFIED

Jan 07, 2021

PROTOCOL INTEGER ID

45521

MATERIALS TEXT

IBioTracker ERthermAC Temperature-Sensitive Live Cell Dye (SCT057, Merck)

https://www.merckmillipore.com/TW/zh/product/BioTracker-ERthermAC-Temperature-Sensitive-Live-Cell-Dye,MM_NF-SCT057?[ReferrerURL=https%3A%2F%2Ffacebook.com%2F&fbclid=IwAR1xNHOYdLL4iTg8IRPS0fHiwa88ZmXFFCu2vYnZTVg4ivkjz3_YKd1AFb8](https://www.facebook.com/?fbclid=IwAR1xNHOYdLL4iTg8IRPS0fHiwa88ZmXFFCu2vYnZTVg4ivkjz3_YKd1AFb8)

IPrepare stock solution of ERthermAC dye (1 mM in DMSO)

IBlack 96-well plate with clear bottom

BEFORE STARTING

IDifferentiate brown/white adipocytes in a black 96-well plate with clear bottom

IPrepare staining medium: 250 nM ERthermAC dye in DMEM-H serum-free medium

1 Wash differentiated adipocytes with PBS twice

- 2 Add 100 μ L of staining solution per well of a 96-well plate
- 3 Incubate 37°C for 30~45 min
- 4 Wash with PBS twice
- 5 Add 90 μ L of DMEM-H serum-free and phenol red-free medium each well
- 6 Equilibrate plate at 25 °C for 15 min
- 7 Measure 2-3 points of basal red fluorescence (ex/em, 543/590) in a microplate reader
- 8 Add 10 μ L of medium or 10 mM CL-316,243 (CL)/forskolin (Fsk) each well (10x dilution, final 1 μ M)
- 9 Measure the red fluorescence (ex/em, 543/590) every 5 min up to 1.5 ~ 2 hr
- 10 Calculate the decrease of fluorescence over time after drug treatment