

Nov 20, 2024

Minibulk v2 (modified Prime-seq)

Forked from [prime-seq](#)

DOI

[dx.doi.org/10.17504/protocols.io.kxygx34qog8j/v1](https://doi.org/10.17504/protocols.io.kxygx34qog8j/v1)


Daniel V Brown¹

¹Walter and Eliza Hall Institute



Daniel V Brown

Walter and Eliza Hall Institute



OPEN ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.kxygx34qog8j/v1](https://doi.org/10.17504/protocols.io.kxygx34qog8j/v1)

Protocol Citation: Daniel V Brown 2024. Minibulk v2 (modified Prime-seq). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.kxygx34qog8j/v1>

Manuscript citation:

This protocol is a fork of Prime-Seq:

Janjic, A., Wange, L.E., Bagnoli, J.W. et al. Prime-seq, efficient and powerful bulk RNA sequencing. *Genome Biol* 23, 88 (2022).

<https://doi.org/10.1186/s13059-022-02660-8>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 19, 2023

Last Modified: November 20, 2024

Protocol Integer ID: 92532

Keywords: bulk RNA sequencing, RNA-seq, transcriptomics, SCRB-seq



Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

Cost-efficient library generation by early barcoding has been central in propelling single-cell RNA sequencing. Here, we optimize and validate prime-seq, an early barcoding bulk RNA-seq method. We show that it performs equivalently to TruSeq, a standard bulk RNA-seq method, but is fourfold more cost-efficient due to almost 50-fold cheaper library costs. We conclude that prime-seq is currently one of the best options to set up an early barcoding bulk RNA-seq protocol from which many labs would profit.

WEHI advanced Genomics facility amendments from Prime-seq:

The plasticware has been changed from 96w plates to 384w and volumes have been reduced 1/4. The reverse transcription of Prime-Seq does not use RNase inhibitors, which I have included.

The primers used for full length cDNA amplification have been taken from 10x Genomics chemistry. This adds compatibility with ONT sequencing (FLT-Seq) developed by Jafar and compatibility with 10x Genomics indexing primers set TT.

Guidelines

- All reagents and plastic-ware can be found in the 'Materials' section.
- Use only RNase free supplies and clean all surfaces and tools with RNase Away prior to working
- Make sure all steps involving cell lysate and RNA before reverse transcription are carried out swiftly and on ice.
- All primer sequences are listed below:

A	B	C	D	E	F
Oligo	Vendor	Purification	Working Co nc.	Sequence	Notes
MB2_KF_DS RT pr imer	IDT	Standard Des alting	100 µM	ACACGACGCTCTTCCGATCTN NNNNNNNNYRNNNNYRNNN[10n t_BC]TTTTTTTTTTTTTTTTTT TTTTTTTTTTTTTVN	
Biotin RAP-TSO	Sigma	RNase-Free H PLC	100 µM	/5Biosg/AAGCAGTGGTATCAA CGCAGAGTACrGrGrG	Biotin red uces conc atemers
RAP_PCR-1_Fwd	IDT	Standard Des alting	10 µM	CTACACGACGCTCTTCCGAT* C*T	PCR1 pri mer
10uM RAP-FL_Re v	IDT	Standard Des alting	10 µM	AAGCAGTGGTATCAACGCAG* A*G	PCR1 pri mer
tenX_adaptTop	IDT	Standard Des alting	1.5 µM	5Phos/GATCGGAAGAGCACAC GTCTGAACTCCAGTCAC	Orded as pre-duple xed DNA
tenX_adaptBot	IDT	Standard Des alting	1.5 µM	GCTCTTCCGATC*T	Orded as pre-duple xed DNA

Specific barcoded oligodT (MB2_Kx_DS_RT primer):

Base design is based off Karst et al., <https://doi.org/10.1038/s41592-020-01041-y>

- MB2 = protocol
- KF = Karst design forward UMI = NNNNNNNNNYRNNNNYRNNN
- KR = Karst design reverse UMI = NNNNNNNNNRYNNNNRYNNN
- Forward and reverse designs to improve complexity in Illumina sequencing
- DS = David Squire designed the well barcodes and thermodynamic properties
- RT = reverse transcription



mb_v2_384w_v1.xlsx



Materials

MATERIALS

- ✕ DNase I Reaction Buffer - 6.0 ml **New England Biolabs Catalog #B0303S**
- ✕ DNase I (RNase-free) - 1,000 units **New England Biolabs Catalog #M0303S**
- ✕ SPRIselect reagent kit **Beckman Coulter Catalog #B23317**
- ✕ Deoxynucleotide Solution Mix - 40 umol of each **New England Biolabs Catalog #N0447L**
- ✕ Exonuclease I (E.coli) - 3,000 units **New England Biolabs Catalog #M0293S**
- ✕ 5 M Sodium chloride (NaCl) **Sigma Aldrich Catalog #S5150-1L**
- ✕ Maxima H Minus Reverse Transcriptase (200 U/uL) **Thermo Fisher Scientific Catalog #EP0752**
- ✕ SUPERaseIN RNase Inhibitor **Thermo Fisher Scientific Catalog #AM2696**
- ✕ NEBNext Ultra II FS DNA Library Prep with Sample Purification Beads - 24 rxns **New England Biolabs Catalog #E6177S**
- ✕ EDTA **Sigma Aldrich Catalog #E7889**
- ✕ Ethanol absolute **Carl Roth Catalog #9065.4**
- ✕ KAPA HiFi 2x RM **Kapa Biosystems Catalog #KR0370**
- ✕ UltraPure DNase/RNase Free Distilled Water **Catalog #10977-049**
- ✕ Buffer EB **Qiagen Catalog #19086**
- ✕ Trizma hydrochloride solution **Sigma Aldrich Catalog #T2694**
- ✕ Aluminium seals for cold storage **Catalog #391-1275**
- ✕ PCR Seals **Thermo Scientific Catalog #AB0558**
- ✕ Dual Index Kit TT Set A **10x Genomics Catalog #1000215**

Protocol materials

☒ Deoxynucleotide Solution Mix - 40 umol of each **New England Biolabs Catalog #N0447L** Materials

☒ NEBNext Ultra II FS DNA Library Prep with Sample Purification Beads - 24 rxns **New England Biolabs Catalog #E6177S**

Materials

☒ SPRIselect reagent kit **Beckman Coulter Catalog #B23317** Materials

☒ EDTA **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7889** Materials

☒ Aluminium seals for cold storage **Catalog #391-1275** Materials

☒ Buffer EB **Qiagen Catalog #19086** Materials

☒ Trizma hydrochloride solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T2694** Materials

☒ PCR Seals **Thermo Scientific Catalog #AB0558** Materials

☒ 5 M Sodium chloride (NaCl) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S5150-1L** Materials

☒ Exonuclease I (E.coli) - 3,000 units **New England Biolabs Catalog #M0293S** Materials

☒ SUPERaseIN RNase Inhibitor **Thermo Fisher Scientific Catalog #AM2696** Materials

☒ Ethanol absolute **Carl Roth Catalog #9065.4** Materials

☒ UltraPure DNase/RNase Free Distilled Water **Catalog #10977-049** Materials

☒ Dual Index Kit TT Set A **10x Genomics Catalog #1000215** Materials, Step 64

☒ DNase I (RNase-free) - 1,000 units **New England Biolabs Catalog #M0303S** Materials

☒ KAPA HiFi 2x RM **Kapa Biosystems Catalog #KR0370** Materials

☒ Maxima H Minus Reverse Transcriptase (200 U/uL) **Thermo Fisher Scientific Catalog #EP0752** Materials

☒ DNase I Reaction Buffer - 6.0 ml **New England Biolabs Catalog #B0303S** Materials

Safety warnings

- ⚠ Please follow all Manufacturer safety warnings and recommendations.
Keep a separation of pre and post PCR steps

Before start

Wipe bench surfaces with PCR clean wipes and keep working environment clean.



Preparation

12m

- 1 Clean all surfaces and pipettes with PCR clean wipes 5m
- 2 Turn on a lab oven to 50 °C and another to 65 °C
- 3 Thaw frozen buffers and primers on ice 10m
- 4 Prepare fresh 80% EtOH 2m
- 5 Prepare 900 µL diluted AMPure RNA XP beads: 1:4 in SPRI buffer per tube
 225 µL RNA beads to 675 µL SPRI buffer
- 6 Input to minibulk v2 is extracted RNA or sorted cells. It is essential, however, that the samples either have the same input or that they are normalized after the RNA is extracted, otherwise sequencing depth per sample will be impacted.

STEP CASE

Extracted RNA 76 steps

Follow this case if you are testing samples that have already been **RNA extracted**, were DNase treated, and have been normalized to the same input. The step here will guide you in preparing RNA-seq libraries and sequencing.

Example: 1 ng of total RNA from a cell line

Sample Preparation


- 7 Normalize the samples so they are all the same concentration.

Note

Minimum: 100 pg
Maximum: 4 ng



8 Transfer  1.2 μL of the **normalized RNA** to a new plate. Return the unnormalized RNA to the freezer.

9 You may store the sample plates at  $-80\text{ }^{\circ}\text{C}$ at this point



Reverse Transcription

30m

10 Prepare **Reverse Transcription Mix**

Store on ice until use

A	B	C	D
Reagent	Well μL	384w plate + 20%	Conc
UltraPure Water	0.5625	285.1	
Maxima RT Buffer (5x)	0.5	253.4	1x
100mM dNTPs	0.0125	6.3	500nM
100uM Biotin RAP-TSO	0.025	12.7	1uM
SuperaseIN	0.125	63.4	1U
Maxima H Minus RT	0.075	38.0	1U
Total (incl RNA)	2.5	1267.2	

11 Remove 384w plates from -80°C freezer, thaw on ice and quick spin

12 Incubate  $65\text{ }^{\circ}\text{C}$ for  00:05:00

After this step increase the oven to  $80\text{ }^{\circ}\text{C}$ for later use

13 Add  1.3 μL **Reverse Transcription Mix** to each well with the FlexDrop

14 Incubate for  01:00:00 at  $50\text{ }^{\circ}\text{C}$

1h






15 Incubate for  00:10:00 at  $80\text{ }^{\circ}\text{C}$

10m




cDNA Pooling & Purification

40m

- 16 Pool all wells of each plate by inversion into a SBS reservoir using a centrifuge
 400 x g, 4°C, 00:02:00
collect into a 5mL LoBind tube.
- 17 Measure the approximate volume in each cDNA pooled tube by measuring by pipette
Should be less than  870 µL
- 18 Perform 1x SPRI ratio cleanup with SPRI beads **diluted** 1:4 in SPRI buffer
- 19 Incubate for  00:05:00 at  Room temperature to allow binding of the cDNA onto beads
- 20 Wash with  2 mL of **80% EtOH** while the tube is on the magnet. Discard the supernatant

Note

Volume of EtOH should be adjusted depending on the number of samples. More samples will require more EtOH to cover the beads completely.





- 21 Repeat wash step once more
- 22 Remove tubes from magnet and quick spin in centrifuge.
Use a large centrifuge for 5mL tubes if necessary
- 23 Air dry beads for  00:02:00

2m

Note

Daniel prefers 2 minutes to 5 minutes to reduce risk of overdrying. It is important to regularly check the beads and avoid over-drying.








- 24 Elute the beads in  45 μL of **UltraPure Water**
- 25 Incubate for  00:05:00 at RT off the magnet
Incubate for  00:02:00 on the magnet and transfer  43 μL to a new PCR tube or plate

Exonuclease I Treatment





- 26 Add 2 μL of ExoI Buffer (10x) and 1 μL of Exonuclease I. Incubate as follows:

A	B	C
Step	Temp	Time
Digest	37	20
Heat inactivation	80	10
Storage	4	Infinite

- 27 Perform 0.8x SPRI bead size selection with  40 μL undiluted SPRI beads.
- 28 Incubate for  00:05:00 at Room temperature to allow binding of the cDNA onto beads
- 29 Place the tube on the magnet stand for  00:05:00 and discard supernatant
- 30 Wash with  180 μL of 80% EtOH while the tube is on the magnet, discard the supernatant
- 31 Repeat wash step once more
- 32 Air dry beads for  00:02:00

**Note**

Depending on temperature and humidity, the beads may dry faster. Therefore it is important to regularly check the beads and avoid over-drying.

- 33 Elute the beads in  22 μ L of UltraPure Water
- 34 Incubate for  00:05:00 at RT off magnet
 00:02:00 on magnet then transfer  20 μ L to a new PCR tube or plate

Full length cDNA Amplification

30m

- 35 Prepare **Pre Amplification Mix**

A	B
Reagent	1x
KAPA HiFi 2x RM	25 μ L
10uM RAP_PCR-1_Fwd	2 μ L
10uM RAP-FL_Rev	2 μ L
UltraPure Water	1 μ L
Total	30 μL

- 36 Add  30 μ L **Pre Amplification Mix** to  20 μ L cDNA sample

- 37 Incubate the Pre Amplification PCR as follows:



A	B	C	D
Step	Temperature	Time	Cycles
Initial Denaturation	98 C	3 min	1 cycle
Denaturation	98 C	15 sec	10 - 18 cycle s*
Annealing	65 C	30 sec	



A	B	C	D
Elongation	72 C	4 min	
Final Elongation	72 C	10 min	1 cycle
Storage	4 C	∞	

Note

Adjust the number of cycles based on input (sample number, cell number, or concentration).





As a general guide we recommend:

Daniel's experience is it is better to err on more PCR cycles for this full length PCR and less for the 2nd library PCR.


A	B
Total RNA Input	Cycle s
1 ng	18
10 ng	16
50 ng	14
100 ng	12
>500 ng	10

cDNA Bead Purification

30m

- 38 Mix sample with  40 µL **SPRIselect** for a 0.8x SPRI to sample ratio
- 39 Incubate for  00:05:00 at  20 °C (Room Temp)
- 40 Place the tube on the magnet stand until clear (~3 min) and discard supernatant
- 41 Wash with  180 µL of **80% EtOH** while the tube is on the magnet, discard the supernatant
- 42 Repeat wash step once more





43 Air dry beads for  00:02:00

2m

Note


Depending on temperature and humidity, the beads may dry faster. Therefore it is important to regularly check the beads and avoid over-drying.

44 Elute cDNA in  12 μL **UltraPure Water**

45 Incubate for  00:05:00 at RT and transfer  10 μL to a new PCR tube or plate

II

Note

Stopping Point. cDNA samples can be safely stored at  -20 °C and protocol can be continued at a later date.

cDNA Quantification and Quality Check

20m

46 Quantify the cDNA using Qubit DNA HS following the manufacturer's protocol. Use 2 μL of clean cDNA for quantification.

*


47 Quality check the cDNA using the Tapestation D5000 tape.

Note

Passing the cDNA quality check does not guarantee that the data will be of high quality, however, if the cDNA fails the quality check it will usually not yield good libraries and will therefore generate lower quality data.



Expected result



Example tapestation D5000 HS

Library Preparation

15m

- 48 The input to library preparation is amplified full length cDNA from each pooled 384w plate.

Note

Before starting, read the library preparation section carefully as there are a few steps that are very time sensitive.

- 49 Prepare **Fragmentation Mix**

Take note of the total cDNA input in ng as it will be relevant in later PCR step


A	B
Reagent	1x
Ultra II FS Reaction Buffer	1.4 μ L
Ultra II FS Enzyme Mix	0.4 μ L
Amplified cDNA	2.5 μ L
Water	1.7 μ L
Total	6 μL

Note

Ensure that the Ultra II FS Reaction Buffer is completely thawed. If a precipitate is seen in the buffer, pipette up and down several times to break it up, and quickly vortex to mix. Place on ice until use.

Note

Vortex the Ultra II FS Enzyme Mix for 5-8 seconds prior to adding to master mix for optimal performance.

50 Vortex the **Fragmentation Mix** for  00:00:05 and incubate in thermocycler immediately

51 Incubate the Fragmentation reaction in a thermocycler as follows:



A	B	C
Step	Temperature	Time
Pre-Cool	4 C	∞
Fragmentation	37 C	5 min
A Tailing and Phosphorylation	65 C	30 min
Storage	4 C	∞

Note

Set heated lid to 75° C. Make sure the lid is at the correct temperature before you start the reaction.
After placing your samples, skip the 4C step to begin the 37C incubation once you have added your samples.

Adapter Ligation



45m


52 Prepare **Adapter Ligation Mix**



A	B
Reagent	1x
NEBNext Ultra II Ligation Master Mix	6 µL
NEBNext Ligation Enhancer	0.2 µL
1.5uM IDT custom adapter (pre-annealed tenX_top and tenX_bottom)	0.5 µL
Total	6.7 µL

53 Add  6.7 µL **Adapter Ligation Mix** to each fragmentation sample



54 Incubate for  00:15:00 at  20 °C in a thermocycler
Set heated lid to off

55 Add  7.7 µL Buffer EB to Samples


56 Mix samples with  16 µL **SPRI select beads** for a 0.8x SPRI ratio

Note

The volume of SPRI select beads used during library size selection can be adjusted based on desired library size. Optimization for your samples may be required.

57 Incubate for  00:05:00 at  20 °C (Room Temp)

58 Place the plate on the magnet stand until clear and discard supernatant.

59 Wash with  180 µL of **80% EtOH** while the plate is on the magnet. Discard the supernatant

60 Repeat wash step once more

61 Air dry beads for 00:02:00

2m

Note

Depending on temperature and humidity, the beads may dry faster. Therefore it is important to regularly check the beads and avoid over-drying.

62 Elute samples in 22 μ L water for 00:05:00 off magnet

7m

Incubate 00:02:00 on magnet

Transfer 20 μ L sample to a new tube

Library PCR

30m

63 Prepare **Library PCR Mix** by adding 25 μ L Kappa HiFi polymerase to the 0.2mL tube containing purified adapter ligated sample.

64 Add 5 μ L of Dual Index Kit TT Set A 10x Genomics Catalog #1000215 to each sample.



Take care to change tips and avoid cross contamination

Note

Ensure a different index for each samples as this is the unique index that will be used for demultiplexing libraries.

65 Incubate the **Library PCR** reaction as follows:

A	B	C	D
Step	Temperature	Time	Cycles
Initial Denaturation	98 C	30 sec	1 cycle
Denaturation	98 C	20 sec	8 cycles *
Annealing	54 C	30sec	
Elongation	72 C	20 sec	



A	B	C	D
Final Elongation	72 C	2 min	1 cycle
Storage	4 C	∞	

Note

Adjust the number of cycles based on cDNA input.
As a general guide we recommend:





A	B
cDNA Input	Cycle s
20 ng	8
10 ng	9
5 ng	10

Note

Daniel found the more PCR cycles used the more the library distribution skews to the left.
Better to err on the fewer PCR cycles side.

Library Double Size Selection

1h

- 66 Mix Index PCR with  25 µL 0.5x ratio **SPRI select beads**
- 67 Incubate for  00:05:00 at  20 °C (Room Temp)
- 68 Place the plate on the magnet stand until clear and **transfer**  75 µL **supernatant to clean well.**





**Note**

Be careful **NOT** to discard! This is your library.


69 Mix supernatant with  15 µL total 0.8x ratio **SPRI select beads**

Note


The volume of SPRI select beads used during library size selection can be adjusted based on desired library size. Optimization for your samples may be required.

70 Incubate for  00:05:00 at  20 °C (Room Temp)

71 Place the plate on the magnet stand until clear and discard supernatant.

72 Wash with  180 µL of **80% EtOH** while the plate is on the magnet, discard the supernatant

73 Repeat wash step once more

74 Air dry beads for  00:02:00

2m

Note


Depending on temperature and humidity, the beads may dry faster. Therefore it is important to regularly check the beads and avoid over-drying.

75 Elute in  22 µL **UltraPure Water**.

76 Incubate for  00:05:00 and then place on magnet until clear.

Transfer  20 µL eluted library to new tube.

**Note**

Stopping point. The libraries can be safely stored at  -20 °C until they will be QCed and sequenced.


Library QC and quantification

30m

77 Quantify and quality control the library using the Tapestation D1000 kits.

Expected result

Libraries will typically exceed 1-5 ng/μl concentration



Sample on D1000 tape after 8 PCR cycles

Sequencing

78 Samples should be submitted according to your Sequencing Facility specifications.

At least **8 cycles** are required for the Index Read 1 (**i7**)

**28 cycles for Read 1 (BC+UMI).**

Read 2 (DNA) should be adjusted based on the quality of the genome being mapped to, but for human and mouse 50 cycles is sufficient.

Some potential sequencing options:

A	B	C	D	E	F
Sequencer	Read 1	Read 2	Index Read (i7)	Index Read (i5)	Kit
NovaSeq	28	90	10	10	NovaSeq SP v1.5 100 cycle
NextSeq 2000	28	52	10	10	NextSeq 2000 P3 50 cycles
HiSeq	28	110	10	0	HiSeq 3000/4000 150 cycles

Bioinformatic preprocessing

6h

- 79 The fastqs may be preprocessed with your pipeline of choice. I prefer zUMIs where the parameters are:

```
file1:
  name: path_to_read1.fastq.gz
  base_definition:
    - BC(19-28)
    - UMI(1-18)
file2:
  name: path_to_read2.fastq.gz
  base_definition:
    - cDNA(1-90)
```

- 79.1 When I have multiple plates in a run I have the bclconvert script write the index reads to file and concatenate fastqs from all plates into a single file:

```
cat *I1_001.fastq.gz > combined_S1_I1_001.fastq.gz
cat *R1_001.fastq.gz > combined_S1_R1_001.fastq.gz
cat *R2_001.fastq.gz > combined_S1_R2_001.fastq.gz
```

- 79.2 I then use the index read as the first part of the well barcode:



```
sequence_files:
  file1:
    name: combined_S1_I1_001.fastq.gz
    base_definition:
      - BC(1-8)
  file2:
    name: combined_S1_R1_001.fastq.gz
    base_definition:
      - BC(19-28)
      - UMI(1-18)
  file3:
    name: combined_S1_R2_001.fastq.gz
    base_definition:
      - cDNA(1-90)
```

- 79.3 You will need to create a well barcode whitelist by concatenating the index reads used to all well barcodes:



```
# Read index read 1 barcodes
plateBC <- read.csv(here::here(
  "i7_only_Kit_TT_Set_A.csv"
))

# Read well barcodes
wellBC <- read.csv(here::here(
  "barcodeOrder_v1.csv"
))

# Define the vectors to be concatenated
vector1 <- plateBC$i7_8nt
vector2 <- wellBC$Well_BC

# Get all combinations
combinations <- expand.grid(vector1, vector2)

# Concatenate the combinations
combinations <- as.data.frame(paste(combinations$Var1,
  combinations$Var2, sep=""))

# Write the result to a CSV file
write.csv(combinations, file = here::here("all_sample_BCs.csv"),
  row.names = FALSE, quote = F, col.names = FALSE)
```

all_sample_BCs.csv is the barcode_file in the zUMIs yaml