



JAN 31, 2024

# 🌐 Juxtacellular Recordings in Ventral tegmental area and Locus coeruleus

Cristian González-Cabrera<sup>1</sup>, Kaushik S. More<sup>1</sup>, Matthias Prigge<sup>1</sup>

<sup>1</sup>Neuromodulatory Network Group, Leibniz Institute for Neurobiology, Magdeburg

ASAP Collaborative Research Network

TeamPrigge



priggelab

OPEN  ACCESS



## ABSTRACT

This protocol describes juxtacellular recordings in the VTA and LC of adult mice, involving anesthesia, precise brain targeting, and neuron labeling. Using specialized electrodes and a dual-channel amplifier, neuronal activity is recorded and modulated with current injections. Following recordings, brains are perfused, fixed, and sectioned for detailed morphological analysis.

## SAFETY WARNINGS

### DOI:

[dx.doi.org/10.17504/protocols.io.3byl4q4novo5/v1](https://dx.doi.org/10.17504/protocols.io.3byl4q4novo5/v1)

**Protocol Citation:** Cristian González-Cabrera, Kaushik S. More, Matthias Prigge 2024. Juxtacellular Recordings in Ventral tegmental area and Locus coeruleus. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.3byl4q4novo5/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working



- Ensure all solutions and equipment are prepared and calibrated according to your lab standards.
- Follow all ethical guidelines and protocols for animal care and use.
- Maintain a sterile environment throughout the procedure to prevent contamination.

Created: Jan 30, 2024

Last Modified: Jan 31, 2024

PROTOCOL integer ID: 94389

Funders Acknowledgement:

ASAP

Grant ID: 020505

## Preparation of Adult Mice for juxtacellular recordings

- 1
  - Use adult mice weighing 25-30 g.
  - Anesthetize the animal with urethane at a dosage of 1.25 g/kg.
  - Place the mouse in the stereotaxic frame and align the head angle according to Bregma and Lambda, with a tolerance of a maximum of 50 microns between them.
  - Expose the skull by shaving the mouse's head and making a sagittal incision on the skin.
  - Find the coordinates and drill the bone. Once the dura is exposed, carefully remove it with the forceps (n°).

**VTA: A/P 3.25, L/M 0.5**

**LC: A/P 5.34, L/M 1.1**

- Drill a small hole around Lambda and install the reference screw.



## Juxtacellular set-up

- 2
  - Employ an intracellular recording amplifier (NeuroData Dual channel IR-283), an analog-to-digital converter (CED Micro-1401-3), and the Spike2 software..
  - Use 5-16 M $\Omega$  (40-45 M $\Omega$  for intracellular) glass electrodes with a tip diameter of less than 1  $\mu$ m.
  - Fill the electrodes with a solution containing:
  - 250 mM k-gluconate, 5 mM KCl, 1 mM MgCl<sub>2</sub>, 2 mM EGTA, 5 mM HEPES, 2 mM MgATP, and Tetramethyl-Rhodamin Biocytin tracer or Neurobiotin (2% w/v); pH 7.2.

## Juxtacellular/Intracellular Recordings and Labelling

- 3
  - After penetrating the brain, ensure the electrode impedance is correct by injecting 0.5 nA pulses and adjusting the bridge/balance module.
  - Move down the electrode at a constant low speed until reaching the dorsal VTA ( $z=-4.0$ ).
  - Start the exploration by moving down not faster than a micron per second.
  - Isolate spontaneously active neurons until the extracellular spike amplitude is approximately 0.7 to 1.0 mV (for intracellular penetration about 2.0-2.5 mV.)
  - After baseline recordings, label neurons using the juxtacellular (Pinault, 1996) or intracellular methods.

- 4
  - For juxtacellular labeling:
  - Label neurons by injecting a small current (1-4 nA for 5-10 minutes, 250 ms ON-250 ms OFF). Make sure that the neuron is being modulated by the current injection.
- 5
  - For intracellular labeling:
  - Apply AC pulses to gain intracellular access. If unsuccessful, adjust by 2 micrometers and reapply AC pulses.
  - Once intracellular access is achieved, stabilize neurons with negative current.
  - Label neurons by injecting a small current (1-2 nA for 10-15 minutes, 250 ms ON-250 ms OFF).
- 6 Wait for at least 2 hours before perfusion.
- 7

### Post-Recording and Labeling Process:

- 8
  - Perfuse the animals via the ascending aorta with 0.01 M PBS followed by 4% paraformaldehyde in PBS.
  - Post-fix brains between 6-12 hours in 4% PFA.
  - Cryoprotect the brain in 30% sucrose until it sinks.
  - Section the brains at 40  $\mu$ m using a freezing microtome.