

JAN 25, 2024

OPEN ACCESS



DOI:
dx.doi.org/10.17504/protocols.io.q26g7pkd8gwz/v1

External link:
<https://clellandlab.ucsf.edu/protocols>

Protocol Citation: Katie Jing Kay Lam, Olubankole Arogundade, Claire D Clelland 2024. Copy number variation analysis by ddPCR. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q26g7pkd8gwz/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Created: Jan 25, 2024

Copy number variation analysis by ddPCR

Katie Jing Kay Lam¹, Olubankole Arogundade¹, Claire D Clelland¹

¹University of California, San Francisco



Katie Jing Kay Lam
University of California, San Francisco

ABSTRACT

This protocol describes copy number variation analysis using Bio-Rad Droplet Digital PCR ddPCR and QuantaSoft Software (with modifications).

MATERIALS

A	B	C
Reagent	Manufacturer	Catalog No.
ddPCR Copy Number Assay: RPP30, Human, Homo sapiens	Bio-Rad	dHsaCP1000485
ddPCR Supermix for Probes (No dUTP)	Bio-Rad	1863024
Droplet Generation Oil for Probes	Bio-Rad	1863005

A	B	C
Equipment/Consumable	Manufacturer	Catalog No.
C1000 Touch Thermal Cycler with 96-Deep Well Reaction Module	Bio-Rad	1851197
ddPCR 96-Well Plates	Bio-Rad	12001925
DG8 Cartridges	Bio-Rad	1864008
DG8 Cartridge Holder	Bio-Rad	1863051
DG8 Gaskets	Bio-Rad	1863009
PCR Plate Heat Seal, foil, pierceable	Bio-Rad	1814040
PX1 PCR Plate Sealer	Bio-Rad	1814000
QX200 Droplet Digital PCR System	Bio-Rad	

Last Modified: Jan 25, 2024

PROTOCOL integer ID: 94179

Keywords: Bio-Rad, ddPCR

ddPCR workflow

1




Figure 1. ddPCR workflow ^[1]

ddPCR reaction set up

- 2 Thaw all components to room temperature

3 Prepare 20X target primers/probe

A	B
Forward primer (100 µM)	18 µL
Reverse primer (100 µM)	18 µL
Target probe (FAM) *light sensitive	5 µL
RNase/DNase-free H ₂ O	59 µL

Note: Store at -20°C

4 Prepare Master Mix – Number of reactions + 1 (as extra)

A	B
2X ddPCR Supermix for Probes (No dUTPs)	11 µL
20X target primers/probe (FAM)	1.1 µL
20X reference primers/probe (HEX/VIC)	1.1 µL
RNase/DNase-free H ₂ O	Up to 22 µL (including DNA sample)

Note: RRP30 is used as the reference primer/probe in this protocol

5 For each PCR tube: DNA sample (Up to 350 ng) + Master Mix = 22 µL

Droplet generation

6 Insert DG8™ Cartridge into cartridge holder




Figure 2. Loaded DG8 Cartridge [2]

7 FIRST - Load >20µL sample into wells

Note:

- Add a 50:50 mix of Supermix & H₂O in empty wells
- All 8 wells must be loaded with samples or Supermix + H₂O
- Make sure there are no bubbles

8 Load 70 µL of generation oil into all 8 wells

Note:

- Make sure samples are loaded before generation oil
- Generation oil is one-time use/can only be used within a day once opened

9 Apply a DG8 Gasket on cartridge

10 Load cartridge into QX200™ Droplet Generator



Figure 3. Loaded GD8 Cartridge in QX200™ Droplet Generator [3]

Transfer Droplets

- 11 Droplets should be transferred to 96-well plate within 5 mins of generation



- 12 Transfer 40 μ L droplets to a ddPCR 96-Well Plate
Note: Pipette **SLOWLY** for droplet suction and dispensing

- 13 Cover 96w plate with a foil seal sheet
Note: The side with red line should be facing up




Figure 4. 96-well plate covered with foil seal sheet ^[1]

14 Seal plate in a PX1 PCR Plate Sealer at 180°C for 5 secs



Figure 5. Bio-Rad PX1 PCR Plate Sealer ^[4]

ddPCR cycle

15 ddPCR cycle should be started within 40 mins after sealing



16 Place plate into Bio-Rad Thermal Cycler with 96-Deep Well C1000 block



A	B	C	D
Step	Temperature	Time	
Enzyme Activation	94°C	10 mins	Ramp rate
Denaturation	94°C	30 sec	2°C/sec
Annealing/Extension	*50-65°C	1 min	
Repeat from step 2 (39x)			
Enzyme Deactivation	98°C	10 mins	
Hold	4°C	∞	

Note: Optimized annealing temperature when using primers/probe set for the first time by using gradient temperatures

Reading the plate

17 Set plate in room temperature for 1-2mins

Note: The plate can be stored at 4°C for ≤3 days before reading

18 Place plate in Bio-Rad QX200™ Droplet Reader



Figure 6. Bio-Rad QX200™ Droplet Reader [5]

19 Open QuantaSoft Software & set up

19.1 Select a well

19.2 E.g. Sample
 - Enter sample name
 - Experiment: **CNV2**
 - Supermix: **ddPCR Supermix for Probes (no dUTP)**

19.3 E.g. Target 1 (FAM)
 - Enter name
 - Type: **Ch 1 Unknown**

- 19.4 E.g. Target 2 (HEX/VIC)
- Enter name (RPP30)
- Type: **Ch 2 Reference**

- 19.5 Repeat (copy & paste) for every well used in the plate

- 20 Save template & RUN

Data analysis

- 21 Check if there are >10k events for each sample

- 22 Set threshold manually (pink line) - if necessary






Figure 7. ddPCR 2-D plot ^[1,6]

Note: Make sure the positive droplets (blue/green) are separated from the negative droplets (grey)