



VERSION 2

APR 05, 2024

OPEN ACCESS



DOI:

[dx.doi.org/10.17504/protocols.io.yxvm3726l3p/v2](https://dx.doi.org/10.17504/protocols.io.yxvm3726l3p/v2)

**Protocol Citation:** Nimalka Weerasuriya 2024. Baiting *Pythium myriotylum* from Infested Soil . **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.yxvm3726l3p/v2>Version created by [Nimalka Weerasuriya](#)

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

## Baiting *Pythium myriotylum* from Infested Soil V.2

Nimalka Weerasuriya<sup>1</sup>

<sup>1</sup>USDA, Oklahoma State University

Peanut and Small Grain Research Unit, Stillwater, OK



Nimalka Weerasuriya

USDA, Oklahoma State University

### ABSTRACT

Baiting *Pythium* from seedlings in infested soil or infected hulls.

### PROTOCOL REFERENCES

1. Singleton, L. L., Mihail, J. D., & Rush, C. M. (Charles M. (1992). *Methods for research on soilborne phytopathogenic fungi*. APS Press.
2. Stanghellini, M. E., & Kronland, W. C. (1985). Detrimental effect of surface sterilization on isolation of *Pythium* spp. from feeder roots (Abstr.). *Phytopathology*, 75(11), 1333–1334. [\[link\]](#)
3. Garren, K. H. (1966). Peanut (Groundnut) Microfloras and Pathogenesis in Peanut Pod Rot. *Journal of Phytopathology*, 55(4), 359–367. <https://doi.org/10.1111/j.1439-0434.1966.tb02238.x>
4. Jeffers, S. N., & Martin, S. B. (1986). Comparison of Two Media Selective for *Phytophthora* and *Pythium* Species. *Plant Disease*, 70(11), 1038–1043. <https://doi.org/10.1094/PD-70-1038>

Created: Apr 05, 2024

Last Modified: Apr 05, 2024

PROTOCOL integer ID: 97857

## MATERIALS

P<sub>5</sub>ARP plates

CMA/PDA+amp plates

Soil with *Pythium* infestation signs

Seeds of susceptible peanut cultivar; we used PI 378012 from the USDA Germplasm Collection

Fume hood and tools for sterile plating technique

Sieve and beakers for surface sterilization

RO water

Forceps and scalpel/blade

95% ethanol

Flame

For seedling baiting:

4" pots with catch tray

greenhouse or adequate outdoor conditions

Mature banker plants (Alyssum and Peppers) for IPM

## Preparation

- 1 Prep P<sub>5</sub>ARP Plates 1-2 days before plating.  
Prepare working culture plates (CMA, PDA + amp) up to 1 week before plating.

## Seed Pod or Seedling Baiting



8m




- 2 Select whether you want to bait from diseased hulls in soil or from seedlings grown in diseased soil.

## STEP CASE

### Peanut Seedlings 14 steps

To bait from soil infested with *Pythium* using susceptible germinated seedlings.

3 Germinate susceptible seedlings by wrapping in RO-dampened paper towels and incubating for 2-4 days at  28 °C in the dark.  Check daily to re-wet paper towels as necessary.

4 Plant germinated seeds in infested soil, in triplicate, in 4" diameter pots with a drip tray. Maintain greenhouse at  70 % humidity with temperatures between  24 °C and  30 °C .

4.1 Water to maintain a level of 1-2" in drip tray for the first 2 days, and then to maintain soil moisture without excessive water until seedlings reach 4-6" in height.

5 At 4 DAP, trim to 1 seed/pot, or carefully repot in separate pots.

6 Grow for 15-20 days and check for rot symptoms. Brown rot will start at ~4-6 DAP, wilt at ~6-18 DAP. Those exposed to higher temperatures may be more susceptible to rot.

## Harvest and Surface-Sterilization

7 Harvest seedlings. Rinse plants thoroughly in RO water to remove soil particles.

8 Select appropriate samples with black rotted segments.



Rotted segments of peanut seedling suitable for plating

**8.1** Cut up segments into 1/4 inch (1-2 cm) pieces.

**9** Try with or without [M] 95 % (v/v) ethanol soak for ⌚ 00:04:00 or rinse in sterile RO water. Use a sieve 4m in a small beaker.

**Do not use bleach** due to *Pythium* sensitivity. [2]

**9.1** Blot dry on sterile paper towels.

10 Plate on P<sub>5</sub>ARP or other media with bacterial inhibition (CMA or PDA + amp) [4].

## Transfer Culture

11 Check plates after 24-48 h.



12 Hyphal tip transfer using sterile pin to clean plates of CMA (V8 or PDA).

## Oospore Check

13 Check for "gold coin" oospores at plate edges to indicate *Pythium myriotylum*.



- Use CMA for oospore production
- Use full strength PDA to increase hyphal growth for DNA extraction
- Use V8 to increase hyphal growth for DNA extraction (untested)