




Version 2 ▼

Jun 26, 2022

# *In vitro* germination of *Austropuccinia psidii* urediniospores V.2

Alyssa M Martino<sup>1</sup>, Rebecca M Degnan<sup>2</sup><sup>1</sup>School of Life and Environmental Sciences, The University of Sydney, Camperdown NSW 2006, Australia;<sup>2</sup>School of Chemistry and Molecular Biosciences, The University of Queensland, Brisbane City QLD 4072, Australia

1 Works for me

 Share[dx.doi.org/10.17504/protocols.io.3byl4b73jvo5/v2](https://dx.doi.org/10.17504/protocols.io.3byl4b73jvo5/v2) Alyssa Martino

## ABSTRACT

Optimisation of *Austropuccinia psidii* urediniospore germination for use in RNA extraction and cytogenetics.

## DOI

[dx.doi.org/10.17504/protocols.io.3byl4b73jvo5/v2](https://dx.doi.org/10.17504/protocols.io.3byl4b73jvo5/v2)

## PROTOCOL CITATION

Alyssa M Martino, Rebecca M Degnan 2022. In vitro germination of *Austropuccinia psidii* urediniospores . **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.3byl4b73jvo5/v2>  
Version created by [Alyssa Martino](#)

## LICENSE



\_\_\_\_\_ This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Jun 26, 2022

## LAST MODIFIED

Jun 26, 2022

## PROTOCOL INTEGER ID

65406

## Spore collection, desiccation, and storage

2d

- 1 Harvest fresh spores from heavily infected leaves by shaking into a paper bag.
- 2 Move fresh spores to a glass petri dish with no lid. Transfer petri dish to a desiccator with silica gel beads for 24 - 48 hours to dry the spores.
- 3 Sieve spores to remove any plant matter or dirt.
- 4 Aliquot spores into Nunc CryoTubes or similar for short or long term storage.

Spores can be frozen at -80°C or in liquid nitrogen for long-term storage but may reduce germination rates. Use of spores immediately after desiccation is recommended.

## Inocula preparation and plating

2h

- 5 Create 2% water agar plates.
- 6 In a 15 mL centrifuge tube, make up inocula to a concentration of approximately 1 mg/mL in sterile distilled water (SDW) with 5ul/mL Tween 20 (0.05%). Mix thoroughly by gently inverting the tube several times.

A haemocytometer can be used to count and confirm spore concentration.

- 7 Rehydrate spores in inocula for 30 minutes prior to use.

Pipette 2 uL aliquots of inocula onto agar plate, covering the entire plate with distinct 2 uL

8    droplets and seal the plates with parafilm.

Incubation                      8h

9    Germinate in incubator for a minimum of 8 hours in the dark at 18°C and 75% relative humidity.

Determine Germination Percentage                      1h

10    For each aliquot, count the number of germinated spores under a dissecting microscope. A germinated spore is considered to be a spore with a germination tube at least double the length of the urediniospore.

Germination rates of approximately 70 - 80% are expected when using fresh spores.