



Sep 25, 2024

Meat identification protocol

This protocol is a draft, published without a DOI.

openbioscience Adrian¹

¹OBS



openbioscience Adrian

OBS

OPEN  ACCESS



Protocol Citation: openbioscience Adrian 2024. Meat identification protocol. [protocols.io](https://protocols.io/view/meat-identification-protocol-dm3q48mw) <https://protocols.io/view/meat-identification-protocol-dm3q48mw>

Manuscript citation:

Based on <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3550954/>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: September 25, 2024

Last Modified: September 25, 2024

Protocol Integer ID: 108368

Keywords: alkaline, LAMP, pork, identification, DNA

Disclaimer

This protocol is work in progress for now.



Abstract

Identifying species by identifying unique DNA sequences is common practice.

In this case we will identify the presence or absence of pork DNA in a meat sample.

We will use the alkaline method for DNA extraction and LAMP method for DNA amplification and detection.

This is not a new protocol but rather the detailing of a previously published protocol listed in the manuscript citation.

Guidelines

Follow the country specific guidelines.

Materials

Consumables: 3 x 1.5 ml Centrifuge tubes or equivalent, 1 PCR tube that matches the device well size to store the ready to be tested sample (aliquot).

Devices: Pestle and mortar, LAMP/PCR/Water bath and thermometer, scale, pH meter, 1- 10 µl micropipette, 100- 1000 µl micropipette. Optional: Vortex.

Reagents: NaOH, Tris, HCl, Dnase free water, meat to be tested

Safety warnings



Use appropriate PPE and follow Safe Work Practices.

NaOH (sodium hydroxide or caustic soda) needs to be handled with gloves because it's a highly corrosive alkali.

Seek and follow Safety Data Sheets for NaOH and HCl

Before start

Pork meat identification protocol. Meat Preparation protocol based on

<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3550954/>

"About 500 mg of meat was triturated in pestle & mortar with four volumes (4 ml) of 0.2N NaOH. 5 µl of the extract was again mixed with eight volumes (40 µl) of 0.2N NaOH in a sterile micro centrifuge tube. Mix was heated at 75 °C in water bath for 20 min and then added eight volumes (360 µl) of 0.04 M Tris HCl (pH 7.75) for neutralization of pH. 1 µl of final mix containing about 100 ng DNA was used for."



Meat Preparation

- 1 Measure 500 mg pork meat without fat. This is our sample (aliquot).
- 2 Grind (triturate) the meat in a pestle and mortar. This will break the structure and results in a paste/liquid.
- 3 Measure 1 ml of the liquid meat and add it to a sterile centrifuge tube labeled "raw sample". This is done with a 100- 1000 μ l micropipette.
- 4 Prepare 5 ml of NaOH 0.2 mol/L (sodium hydroxide or caustic soda). This will be used twice, once 4 mL and once 40 μ L.
- 4.1 In a new centrifuge tube labelled "NaOH 0.2 mol/L" add 5 mL DNAfree water and 40 mg or NaOH. Mix to dissolve manually or using a vortex.

*For general molarity calculations see **Molarity calculation***

Molar concentration is the number of moles contained in one Liter of a substance. The concentration of a solution is expressed in mol/L.

$C=(m/M)/V$ where V is the volume in liters m is mass in grams and M is molar weight? in grams/mol

<https://www.alloprof.qc.ca/en/students/vl/sciences/calculating-molar-concentration-s1622>

*Molecular Weight of sodium hydroxide (NaOH) =40 g/mol: Na=22.989g/mol O=15.999g/mol H=1.008g/mol (<http://chembook.org/page-nonav.php?chnum=1§=5>) so:
22.989g/mol+15.999g/mol+1.008g/mol=39.996g/mol.*

*so for us: $m=C*M*V=0.2\text{mol/L}*40\text{g/mol}*5\text{mL}= 8\text{ g/L}*5\text{ mL}=40\text{ mg of NaOH}$*

- 5 Add 4 ml from the tube labelled "NaOH 0.2 mol/L" to the "raw sample" tube. *This will break down the cell membranes freeing the DNA and also break down (denature) some of the proteins.*
- 6 Measure 5 μ L from previous step and add it to a new sterile centrifuge tube labeled "prepared sample". This is done with a 1- 10 μ l micropipette.
- 7 Add in the "prepared sample" new tube 40 μ L from the tube labelled "NaOH 0.2 mol/L".
- 8 Heat the resulting mix at 75 Celsius for 20 minutes in a water bath or using the LAMP device or a PCR device. *This will further break down the cell membranes freeing the DNA.*



- 9 Prepare (360 µl) of 0.04 mol/L Tris HCl (pH 7.75).
This will be done by mixing Tris with HCl

MW for TRIS Hydrochloride ($\text{C}_4\text{H}_{11}\text{NO}_3 \cdot \text{HCl}$) 157.6 mol/Liter

0.04 M TRIS TRIS Hydrochloride: $m = CMV = 0.04 \text{ mol/L} * 157.6 \text{ g/mol} * 360 \text{ uL} = 2.26 \text{ mg}$ of TRIS Hydrochloride

Normally we make half the solution and titrate with HCl

https://www.reddit.com/r/labrats/comments/193ei6n/trying_to_prepare_trishcl_buffer_ph_seems/?rdt=42146

<https://www.cshlpress.com/pdf/sample/2014/LabRefV1/LabRefV1Ch1S1.pdf> has collection of reagents

pH calculations= <https://www.youtube.com/watch?v=387t-nh00Mg>

<https://www.cusabio.com/m-296.html#a21> for calculating Tris HCL pH from Tris and Hydrochloric acid

- 10 Add 360 µl of the 0.04 mol/L Tris HCl. This is done for for neutralization of pH.
- 11 Sample 1 µl in a PCR tube and store at -20.

LAMP

12

Protocol references

<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3550954/>