

May 03, 2022

Electroporation of fluorescent antisense molecules into *Bodo saltans* and live imaging

Mastaneh Ahrar¹, g.hurst¹, Ewa Chrostek¹¹Department of Evolution, Ecology and Behaviour, Institute of Infection, Veterinary and Ecological Sciences, University of Liverpool

1

dx.doi.org/10.17504/protocols.io.e6nvwk8e2vmk/v1

Symbiosis Model Systems

Bodo protocols

Ewa Chrostek

This is the protocol used in our Laboratory at the University of Liverpool to electroporate molecules into Bodo saltans.

DOI

dx.doi.org/10.17504/protocols.io.e6nvwk8e2vmk/v1

Mastaneh Ahrar, g.hurst, Ewa Chrostek 2022. Electroporation of fluorescent antisense molecules into Bodo saltans and live imaging. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.e6nvwk8e2vmk/v1>



protocol ,

May 03, 2022

May 03, 2022

61862

Overview

- 1 Antisense peptide nucleic acids (PNAs, synthesised by Panagene, South Korea) were electroporated into *Bodo saltans* cells using Neon[®] Transfection System MPK1025 (Invitrogen).

Bodo culture conditions

- 2 *Bodo saltans* was cultured in a cerophyl-based medium enriched with 3.5 mM sodium phosphate dibasic (Na₂HPO₄)¹. Cultures were incubated at 22 °C in T25 tissue culture flasks containing 20 ml of media bacterized with *Klebsiella pneumoniae subsp. Pneumoniae* (ATCC[®] 700831). 3 to 4

day old cultures were used for electroporation.

Preparing cells for electroporation

3 The steps below describe how to prepare *Bodo* cells for electroporation.

- 3.1 1-Filter the culture through 100 and 8 μ m filter.
- 3.2 2-Harvest the cells by centrifugation at $1200 \times g$ for 12 mins at 19 °C.
- 3.3 3-Wash the cells with 10 ml PBS and centrifuge as above.
- 3.4 4-Re-suspend the cells in 5 ml PBS, count the cells using hemocytometer and take the volume of cells which contains 5×10^5 cells, as recommended for Neon transfection kit for 10 μ l tip.
- 3.5 5-Centrifuge at $1200 \times g$ for 12 mins at 19 °C.
- 3.6 6-Remove the PBS and resuspend the cells in resuspension buffer for electroporation.
- 3.7 7-Add the antisense molecule at the final concentration of 50 μ M to the *Bodo* cells and mix well by pipetting.
- 3.8 8-Aspirate the mix into Neon pipette and electroporate using 1800 V, 10 ms pulse width and 1 pulse.

Preparing electroporated cells for imaging

4 The steps below describe how to prepare electroporated *Bodo* cells for imaging.

- 4.1 1-Empty the pipette of Neon system in a well of a 96-well plate after electroporation, add 5 μ l of PBS, mix, add 15 μ l of warm low melting temperature agarose (Thermo Fisher Scientific), and mix with the cells. Let is set for a few seconds.
- 4.2 2-Add 200 μ l of Hoechst 33342 (Thermo Fisher Scientific) diluted 1:2000 in PBS to the agarose-embedded *Bodo* and incubate 10 mins at room temperature (RT).
- 4.3 3-Wash the agarose with PBS (2 \times 5mins) at RT.
- 4.4 4-Remove the agarose from the well with clean forceps and place it on a microscope slide.
- 4.5 5-Add a drop of a mounting medium (eg. Vectashield, Vector Laboratories), and flatten the agarose as much as you can using the coverslip.
- 4.6 6-Proceed with either fluorescence or confocal imaging.




Figure 1. Confocal image of live *Bodo* cell electroporated with antisense peptide nucleic acid (PNA). A) TAMRA-labelled antisense molecules inside *Bodo* cell. B) DNA of bacterial cells inside *Bodo* stained with Hoechst 33342. C) Overlay of two channels.

References

- 5 Gomaa Fatma, Li ZuHong, Docampo Roberto, Girguis Peter, E. V. *Bodo* saltans culture protocol V.2. *Protocols.io* (2018). doi:10.17504/protocols.io.sh6eb9e

