



Apr 22, 2021

FL-BEEHIV HIV-1 Genotyping and Drug Resistance by Next Generation Sequencing

Brenna M McGruder Rawson¹, Jason Blanton¹¹Florida Department of Health, Bureau of Public Health Laboratories

In Development

dx.doi.org/10.17504/protocols.io.btrnm5e

Jason Blanton

Florida Department of Health, Bureau of Public Health Labora...

ABSTRACT

The Florida Department of Health's Bureau of Public Health Laboratories in Jacksonville is developing a protocol for the Next Generation Sequencing (NGS) of HIV, primarily for the purpose of drug-resistant mutation identification. This HIV-1 protocol uses amplicon based sequencing based on primers designed by the BEEHIVE Consortium (<https://www.beehive.ox.ac.uk/>). The amplified pol gene regions are used in both genotyping and drug resistance determination. Our protocol utilizes newer enzymes with higher fidelity for sequencing and Illumina sequencing technology. The NGS data generated can also be used in surveillance and outbreak monitoring, giving epidemiologist more information about circulating viral genomes. There is also the potential that this protocol can be expanded to whole genome sequencing for HIV-1.

The imminent sunset of ViroSeq (Abbott Molecular) has required many labs to look for new methods to continue identifying HIV-1 drug resistance strains for both clinical management and epidemiological study. NGS was chosen as it is more cost effective than investing in a single pathogen platform. NGS allows for one sample to produce results and data that can aid not just a patient but an entire population.

DOI

dx.doi.org/10.17504/protocols.io.btrnm5e

PROTOCOL CITATION

Brenna M McGruder Rawson, Jason Blanton 2021. FL-BEEHIV HIV-1 Genotyping and Drug Resistance by Next Generation Sequencing. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.btrnm5e>

KEYWORDS

null, HIV-1, HIV-1 genotyping, HIV-1 drug resistance, Next Generation sequencing, drug resistance mutations

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Mar 29, 2021

LAST MODIFIED

Apr 22, 2021

PROTOCOL INTEGER ID

48654

MATERIALS TEXT

RNA Extraction by Qiagen Viral RNA Mini Kit (DSP or Regular)

SSIV VILO Master Mix (Thermofisher, Cat 11756050)

Q5 Master Mix (NEB, Cat M0492S)

Primers (Gall A, et al. Journal of Clinical Microbiology. 2012; 50:12)

Set and primer	Sequence (5'–3')	Positions ^a	Product size ^a
2			
Pan-HIV-1_2F	GGG AAG TGA YAT AGC WGG AAC	1031–1051	3,574 bp
Pan-HIV-1_2R	CTG CCA TCT GTT TTC CAT ART C	4604–4583	
3			
Pan-HIV-1_3F	TTA AAA GAA AAG GGG GGA TTG GG	4329–4351	3,066 bp
Pan-HIV-1_3R	TGG CYT GTA CCG TCA GCG	7394–7377	

^aAccording to HIV-1 reference strain HXB2 (GenBank accession number NC001802).

Single/multichannel pipettes with p20/p200/p1000 tips

Thermocycler

Nuclease-free water

AMPure XP Beads (Beckman Coulter)

Magnetic stand

Tapestation or Agarose gel

Qubit or other quantitation method

Illumina Nextera XT DNA Library Prep Kit

Illumina Nextera v2 Index Kits

Illumina iSeq

Illumina iSeq 100 i1 v2 cartridge

<https://www.smartgene.com/>

Pipeline: HIV1-PR+RT+IN (2.2.0HIV1 V1.6)

DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://www.protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://www.protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

The use of reagents, equipment, or services in this protocol is not an endorsement either stated or implied by the Florida Department of Health.

RNA Extraction

1 Extract RNA using the Qiagen Viral RNA Mini Kit (DSP or RUO)

1.1 Add 5.6 µl carrier RNA per sample

1.2 Follow Kit instructions for extraction

cDNA Synthesis 25m

- 2 Master Mix
 - ▢ 4.0 µl SuperScript IV VILO Master Mix
 - ▢ 6.0 µl Nuclease Free Water
 - ▢ 10.0 µl RNA template

25m

- 3 Run the following protocol on a thermocycler
 - ⌚ 25 °C ⌚ 00:10:00
 - ⌚ 50 °C ⌚ 00:10:00
 - ⌚ 85 °C ⌚ 00:05:00

Amplicon PCR 6m

- 4 Each fragment will need to be amplified in an individual PCR reaction
 - Set 1
 - Pan-HIV-1_2F GGG AAG TGA YAT AGC WGG AAC
 - Pan-HIV-1_2R CTG CCA TCT GTT TTC CAT ART C

Set 2

Pan-HIV-1_3F TTA AAA GAA AAG GGG GGA TTG GG

Pan-HIV-1_3R TGG CYT GTA CCG TCA GCG

- 4.1 Master Mix
 - ▢ 12.5 µl 2x Q5 Master Mix
 - ▢ 0.5 µl Forward Primer [M]20 Micromolar (µM)
 - ▢ 0.5 µl Reverse Primer [M]20 Micromolar (µM)
 - ▢ 6.5 µl Nuclease Free Water
 - ▢ 5.0 µl cDNA template

6m

- 5 Run the following protocol on a thermocycler
 - ⌚ 105 °C Lid
 - ⌚ 50 °C ⌚ 00:00:30
 - ⌚ 98 °C ⌚ 00:00:30
 - 35 cycles ↻
 - ⌚ 72 °C ⌚ 00:05:00
 - ⌚ 4 °C Hold

- 6 Bead clean up using a ratio of 0.5- follow the AMPure XP bead protocol for PCR purification

7 Check fragment on Tapestation or gel

- 7.1 Band size should be
Amplicon 1- 3.5 kB
Amplicon 2- 3.0 kB

Fragment Normalizing and Pooling

8 Fragments can be pooled in eqimolar amounts or in equal concentrations.

8.1 Pool fragments

8.2 Dilute as needed to achieve 1.0 ng input concentration for library preparation

Library Prep

9 Follow Illumina Protocol for Nextera XT DNA Library Sample Prep

Pooling Samples

10 Amplicon quality can effect how many samples can be pooled onto one run. Use caution in deciding how many samples to pool.

Analysis

11 We currently use SmartGene HIV-1 pipeline (<https://www.smartgene.com/>) for analysis Pipeline Name: HIV-1 PR+RT+IN Version 2.2.0HIV1 V1.6 Noise Filter [%] 0.5 Interpretation cut off [%] 5.0

Minimum read depth and additional criteria should be determined by your institution

References

12 Gall A, Ferns B, Morris C, Watson S, Cotten M, Robinson M, Berry N, Pillay D, Kellan P. Universal Amplification, Next-Generation Sequencing, and Assembly of HIV-1 Genomes. Journal of Clinical Microbiology. 2012; 50:12. doi: 10.1128/JCM.01516-12

Cornelissen M, Gall A, Vink M, Zorkrager F, Binter S, Edwards S, Jurriaans S, Bakker M, Ong SH, Gras L, van Sighem A, Bexemer D, de Wolf F, Reiss P, Kellam P, Berkhout B, Fraser C, van der Kuyl AC, the BEEHIVE Consortium. From clinical samples to complete genome: Comparing methods for the extraction of HIV-1 RNA for high-throughput deep sequencing. Virus Research. 2017; 239:10-16. doi: [10.1016/j.virusres.2016.08.004](https://doi.org/10.1016/j.virusres.2016.08.004)