



Oct 11, 2022

Colony PCR

Georgina Elisa Diego Macías¹, Ana Belem García González¹,
Irán Alessandra Chaparro Rodríguez¹, Jair Alexis Gardea Sáenz¹

¹Tecnologico de Monterrey Campus Chihuahua

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.dm6gpjwx8gzp/v1

G Georgina Elisa Diego Macías

ABSTRACT

Protocol for a PCR from a colony

DOI

dx.doi.org/10.17504/protocols.io.dm6gpjwx8gzp/v1

PROTOCOL CITATION

Georgina Elisa Diego Macías, Ana Belem García González, Irán Alessandra Chaparro Rodríguez, Jair Alexis Gardea Sáenz 2022. Colony PCR. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.dm6gpjwx8gzp/v1>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 11, 2022

LAST MODIFIED

Oct 11, 2022

PROTOCOL INTEGER ID

71137

1 Dissolve one colony in **20 µL** of nuclease-free water and heat at 80°C for 10 minutes.

Prepare the following 50 µL reaction in a 0.2mL PCR tube on ice:

2

Component	uL Reaction
Magnesium chloride (MgCl ₂) 25 mM	4 µL
5X Colorless GoTaq Buffer	5 µL
PCR Nucleotide Mix 10 mM	0.5 µL
Forward primer 10 uM	1 µL
Reverse primer 10 uM	1 µL
GoTaq DNA polymerase (5u/uL)	0.2 µL
DNA Colony	10 µL
Nuclease-free water	3.3 µL
TOTAL	25 µL

3 Gently mix the reaction and spin down in the microcentrifuge

4 Cycling conditions for a routine PCR:

A	B	C	D
1 Cycle	Initial activation	2 minutes	94°C
35 Cycles	Denaturation	30 seconds	94°C
	Alignment	30 seconds	50°C
	Extension	2 minutes	72°C
1 Cycle			72°C
Maintenance		∞	4°C