



Sep 10, 2020

# Splitting p1 (1xT75) to p2 (2xT150)

In 1 collection

Andrea Argouarch<sup>1</sup><sup>1</sup>University of California, San Francisco

1 Works for me dx.doi.org/10.17504/protocols.io.8gnhtve

Neurodegeneration Method Development Community

Tech. support email: [ndcn-help@chanzuckerberg.com](mailto:ndcn-help@chanzuckerberg.com)

Andrea Argouarch

University of California, San Francisco

## ABSTRACT

Protocol includes splitting growing cell line from a T75 flask into two T150 flasks for expansion.

## DOI

[dx.doi.org/10.17504/protocols.io.8gnhtve](https://dx.doi.org/10.17504/protocols.io.8gnhtve)

## PROTOCOL CITATION

Andrea Argouarch 2020. Splitting p1 (1xT75) to p2 (2xT150). **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.8gnhtve>

## COLLECTIONS

**Dural Cell Isolation and Culturing - Collection**

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Oct 19, 2019

## LAST MODIFIED

Sep 10, 2020

## PROTOCOL INTEGER ID

28910

## PARENT PROTOCOLS

Part of collection

[Dural Cell Isolation and Culturing - Collection](#)

## STEPS MATERIALS

NAME	CATALOG #	VENDOR
<a href="#">DPBS, no calcium, no magnesium</a>	14190250	<a href="#">Thermo Fisher</a>
<a href="#">DMEM, high glucose, pyruvate</a>	11995073	<a href="#">Thermo Fisher</a>
<a href="#">Penicillin-Streptomycin</a>	15140122	<a href="#">Gibco - Thermo Fisher</a>
<a href="#">Fetal Bovine Serum</a>	97068-091	<a href="#">Vwr</a>
<a href="#">Trypsin-EDTA (0.05%), phenol red</a>	25300062	<a href="#">Thermo Fisher</a>

## Observations

- 1 At 90-100% confluency, split T75 flask into 2xT150 flasks

## Preparation

- 2 Turn off UV lights and clean hood with 70% ethanol
- 3 Clean items with 70% ethanol and bring into hood
  - a. DPBS -/-



**DPBS, no calcium, no magnesium**

by Thermo Fisher

Catalog #: 14190250

b. Sterile Filtered Media



**DMEM, high glucose, pyruvate**

by Thermo Fisher

Catalog #: 11995073



**Penicillin-Streptomycin**

by Gibco - Thermo Fisher


Catalog #: 15140122



**Fetal Bovine Serum**

by Vwr

Catalog #: 97068-091

- c. Label 2xT150 flasks with ID, date, and p2, add 28 mls  of media per flask
- d. 0.05% Trypsin












**Trypsin-EDTA (0.05%), phenol red**

by Thermo Fisher

Catalog #: 25300062

## Culturing

- 4 Aspirate old media
- 5 Rinse by adding 5 mls  5 mL of DPBS per flask and gently swirling
- 6 Aspirate DPBS
- 7 Add 2 mls  2 mL of trypsin per flask and place in incubator  37 °C for 2-3 mins  00:03:00 until cells are starting to detach
  - a. Can also gently tap the side of the flask to detach cells
- 8 Check under microscope
- 9 Add 10 mls  10 mL of fibroblast media per flask to inactivate trypsin
- 10 Collect cell suspension in a 15 ml conical
- 11 Spin at 600 rpm  600 rpm for 4 mins  00:04:00
- 12 Aspirate supernatant, being careful to not aspirate the pellet
- 13 Tap the pellet to resuspend
- 14 Add fibroblast media to 4 mls  4 mL
- 15 Add 2 ml  2 mL of cell suspension per flask

16 Place flask in incubator (5% CO<sub>2</sub> , 37°C  37 °C )

17 Observe and feed the next day, then feed every 2-3 days until 90-100% confluent

Clean Up

18 Throw away biohazard materials properly

19 Clean and sterilize hood with 70% ethanol and turn on UV

20 Update cell culture notes in lab notebook