



Jul 12, 2022

Small scale Lentivirus Production and Infection

Herschel Dhekne¹, Suzanne R R R Pfeffer¹¹Stanford University Department of Biochemistry1 *Works for me* Sharedx.doi.org/10.17504/protocols.io.bp2l61z2zvqe/v1 Suzanne R R Pfeffer

ABSTRACT

This protocol can be used for production and transduction of lentiviral sgRNA, shRNA and protein overexpression in conjunction with generation 2 and generation 3 lentivirus plasmids.

ATTACHMENTS

[465-974.docx](#)

DOI

dx.doi.org/10.17504/protocols.io.bp2l61z2zvqe/v1

PROTOCOL CITATION

Herschel Dhekne, Suzanne R R R Pfeffer 2022. Small scale Lentivirus Production and Infection. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bp2l61z2zvqe/v1>



KEYWORDS

Small scale Lentivirus Production , Small scale Lentivirus Infection

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 16, 2022

LAST MODIFIED

Jul 12, 2022

OWNERSHIP HISTORY

Jun 16, 2022  maria.s

Jun 28, 2022  Suzanne R R Pfeffer

PROTOCOL INTEGER ID

64696

MATERIALS TEXT

Materials:

- BSL-2+ facility cell culture lab
- Addgene plasmids - [psPAX2 addgene Catalog #12260](#) ,
[pMD2.G addgene Catalog #12259](#)), lentiviral vector
- Polyethyleneimine (PEI, Polysciences) [1 mg/mL](#) stock
- HEK 293T cells
- Polyethylene glycol (PEG) 8000
- Polybrene ([10 mg/mL](#))
- **4X lentivirus concentrator solution**6. Store at [4 °C](#) .
- Ultracentrifuge and compatible tubes

Make lentivirus

6h

1



Plate 293T cells at 40% confluency in a 6 well tissue plate submerged under [2 mL](#) medium per well.

2

After [06:00:00](#) , most cells will have attached.

6h

Day 0

3

Prepare DNA mix for transfection:

3.1



Add the following to **100 µL** Optimem per well for transfection:

A	B
1 µg	PsPAX2 1µg helper plasmid (Addgene ##12260)
0.5 µg	VSV-G / pMD2.g (Addgene #12259)
1 µg	Lentivirus vector (see below)

3.2



Add PEI (from a **1 mg/mL** stock) to this mixture solution at ratio 5:1 w/w (PEI:DNA).

Example, **12.5 µg** PEI for **2.5 µg** DNA mix.

3.3



20m

Mix DNA mix gently and incubate for **00:20:00** at **Room temperature**.

3.4



Add the mix to the cells dropwise.

Day 1 (16 hours later)

4

Check for cell viability; at this time, >70% of the cells should be transfected and virus is already being produced and is being released into the supernatant.

NOTE: Removal of residual PEI at this stage by medium change is not essential but will be present in the supernatant.

Day 2

5



2d

⌚ **48:00:00** after transfection, collect the culture supernatant in a BSL-2+ facility; centrifuge in an enclosed rotor and remove supernatant with care. This is “Day-2 virus”.

6



Carefully add an additional **2 mL** complete DMEM medium into each well without splashing or disturbing the monolayer.

7

Bleach all tips and pipettes used to collect the virus.

Day 3

8

3d

⌚ **72:00:00** after transfection, collect the culture supernatant in BSL-2+ facility as before. This is “Day-3 virus”. Day-2 and Day-3 virus are then pooled; Day-2 titre is lower than Day-3.

9



5m

The pooled virus (~ **4 mL**) is transferred into a 15ml tube and centrifuged at **250 x g** for **00:05:00**.

The pellet represents cell debris as well as 293T cells that can contaminate the target cell line to be infected with the virus; care should be taken when aspirating the virus supernatant. Filtration can decrease viral titre and is not required.

10



Prepare **0.5 mL** aliquots of the lentivirus and freeze at **-80 °C**.

Lentivirus Infection

11



Thaw a **0.5 mL** virus aliquot in a **37 °C** water bath, flicking tube gently to facilitate gentle thaw.

12 

Add **1 µL**, **10 mg/mL** Polybrene.

NOTE: Polybrene enhances infectivity but is not essential. Use at 2-8µg/ml depending on the cell type; polybrene can be toxic to cells so take care. HeLa, MEF, 3T3 and A549 cells tolerate up to 8 µg/ml.

13 

Transfer virus mixture to the medium covering 1 well of a 6 well plate containing the target cell line. Polybrene will become diluted in the cell medium to a final concentration of **4 µg/ml**.

14  **48:00:00** post infection, cells are ready for analysis or selection.

2d

Concentrating the virus

15


Rationale: To achieve 100% infection and/or if you have low titers or do not care about precise multiplicity of infection, it is beneficial to concentrate the lentivirus.


4×Lentivirus Concentrator Solution


16 Dissolve **80 g** PEG-8000 and **14.0 g** NaCl in **80 mL** MilliQ water.

17 Add **20 mL**, 10X PBS (**pH 7.4**).

18 Mix with gentle stirring, heating gently only if necessary, until the solids are dissolved then adjust pH to 7.0~7.2.

19 Adjust the final volume to  200 mL .

20 Sterilize by passage through a  0.2 µm filter.



The concentrations of PEG-8000 and NaCl in the stock solution are 40% (w/v) and  1.2 Molarity (M) , respectively.

Virus concentration protocol

21 




Carefully transfer the virus supernatant into a new 50 ml tube.

22 

Add 1 volume of concentrator solution to 3 volumes of virus supernatant (eg.  1 mL concentrator solution for  3 mL virus).

23   

4h 0m 20s

Mix by gentle shaking for ~  00:00:20 then incubate with constant rocking at least  04:00:00 at  4 °C .

 **Overnight** rotation or rocking will enhance recovery.


24 

1h


Spin down at  **1600 x g** for  **01:00:00** at  **4 °C** .

- 25 Carefully remove supernatant without disturbing the pellet.

Pellet size does not necessarily correlate with virus yield.

- 26 

Thoroughly resuspend the viral pellet in PBS or desired medium using 1/10~1/20 of the original volume by gentle pipetting using a 1ml Pipetman.

- 27 Aliquot and store at  **-80 °C** until use.


Alternative Centrifugation- based Virus concentration method

3d 1h 35m

- 28

3d

In case of low transduction efficiency, consider ultracentrifugation as follows:

 **72:00:00** after transfection, collect the virus-containing supernatant in a BSL-2+ facility (take only Day 3 supernatant).




- 29 

5m


Spin down at  **250 x g** for  **00:05:00** at  **Room temperature** .

- 30  

1h 30m

Transfer the precleared supernatant to ultracentrifuge tubes and pellet at  **90000 x g** for  **01:30:00** at  **4 °C** .

- 31 

Remove the supernatant and leave a little less than  **1 mL** in the tube. Use a 1 mL pipette to recover the remaining pellet which may be difficult to see.

32 Make aliquots of  0.2 mL concentrated virus and freeze at  -80 °C .