



APR 15, 2024

🌐 Behavior Hardware Setup Protocol

📁 In 1 collection

Sasha Burwell¹

¹Duke University

ASAP Collaborative Research Network



Sasha Burwell
Duke University

OPEN  ACCESS



ABSTRACT

This protocol details instructions for recreating the conditioning/extinction reward learning assay used in this paper.

DOI:

dx.doi.org/10.17504/protocols.io.ewov19857lr2/v1

Protocol Citation: Sasha Burwell 2024. Behavior Hardware Setup Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.ewov19857lr2/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: Apr 08, 2024

Last Modified: Apr 15, 2024

Keywords: ASAPCRN

To run the conditioning/extinction reward learning assay, any behavioral setup that incorporates head-fixation on a treadmill, a reward delivery system, a tone delivery system, and a lick detection system can be used. Below are the parts and instructions for building the behavioral setup we used. We also incorporated an infrared camera for video recordings/ visualization into the closed box, though video was not analyzed and this component is not necessary per se.

Parts List:

A	B	C	D	E
Vendor	Item	Part Number	Price	Quantity
Box to create a dark and sound attenuated environment				
Med-Associates	Expanded PVC Sound Attenuating Box	ENV-022V	\$547	1
Thor Labs	12x12 Base	MB12	\$151	1
Amazon	Speaker	Z50	\$17	1
https://www.amazon.com/gp/product/B00EZ9XLF8/ref=oh_aui_detailpage_o06_s01?ie=UTF8&psc=1				
	Auxiliary Stereo Extension audio cable	https://www.amazon.com/gp/product/B06X9SJL97/ref=oh_aui_detailpage_o06_s00?ie=UTF8&psc=1	\$10/2	1
	Egg crate Foam	https://www.amazon.com/gp/product/B076JNCQ39/ref=oh_aui_detailpage_o07_s00?ie=UTF8&psc=1	\$11.95	1
	Spray adhesive	https://www.amazon.com/gp/product/B000PCWRMC/ref=oh_aui_detailpage_o06_s00?ie=UTF8&psc=1	\$5.77	1
	Camera tripod	https://www.amazon.com/gp/product/B000PCWRMC/ref=oh_aui_detailpage_o06_s00?ie=UTF8&psc=1	\$6	1

A	B	C	D	E
		on.com/gp/		
product/B01L8QNBFA/ref=oh_aui_				
detailpage_o05_s00?ie=UTF8&psc=1				
Graftek	Basler Ace camera	aca1300-gm	\$422.75	1
	Tripod mount plate	2000029679	\$18.13	1
	Power cable	2200000167	\$52.50	1
	Lens C125-1218-5M	2000034833	\$139.65	1
B&H Photo and Video	Bosch IR Illuminator 5000 SR (850nm)	BOIIR50850SR	\$239	1
	Bosch power supply unit (35W, 1.46A)	BOPSUIR35	\$52.43	1
Reward Delivery System				
Lee Company	Solenoid	LHDA1233315H HDI-PTD-Saline-12V-30PSI	\$76.68	1
McMaster-Carr	60mL Luer Lock Syringe	7510A802	\$2.73	1
	Stainless steel tubing	5560K741	\$3.90	1
	PVC Clear Tubing	6516T11	\$0.23/ft	1
	Thumb Screw	91830A101	\$3.12	1
	On/Off Valve	7033T24	\$2.29	1
	Stainless Steel Dispensing Needle	75165A673	\$13.32	1
Thor Labs	Clamping Fork	CF175	\$10.30	1
	Studded Pedestal Base Adapter	BE1	\$9.49	1
	Post Holder (3")	PH3	\$8.27	1
	½" Optical Post (6")	TR6	\$6.91	1
	½" Optical Post (4")	TR4	\$5.87	1
	Swivel Post Clamp	SWC	\$23.56	1

A	B	C	D	E
3D printer	Lick port	10.5281/zenodo.10928257		1
Treadmill and head-fixation				
Delvie's Plastics	Plexiglas disk (8")	1/8" thickness, 8" diameter	\$5.00	1
Amazon	Silicone rubber sheet black	https://www.amazon.com/gp/product/B005584T7C/ref=oh_aui_detailpage_o06_s02?ie=UTF8&psc=1	\$20.77	1
	Bag clips	https://www.amazon.com/gp/product/B004RTBOK0/ref=oh_aui_detailpage_o08_s00?ie=UTF8&psc=1	\$6 for 7	1
U.S. Digital	Rotary Encoder	H5-100-NE-S	\$68.85	1
	5-pin connector	CA-FC5-W5-NC-10	\$12.25	1
3D printer	Treadmill to rotary encoder adapter	10.5281/zenodo.10928257		1
Thor Labs	Clamping Fork	CF175	\$10.30	1
	Studded Pedestal Base Adapter	BE1	\$9.49	1
	Post Holder (1")	PH1	\$7.03	1
	1/2" Optical Post (1")	TR1	\$4.74	1
	1/2" Optical Post (2")	TR2	\$5.19	1
	Swivel Post Clamp	SWC	\$23.56	1
McMaster-Carr	18-8 Stainless Steel Rod 7in	95412A562	\$11 for 10	2
	18- stainless steel hex nut	91841A005	\$2.61 for 100	4
	Steel wing nut	90866A029	\$10.29 for 100	2
	Nickel-plated brass washer	92917A155	\$8.56 for 100	4
3D Printer	Head post (5")	10.5281/zenodo.10928257		1
	Head post (7")		1	

A	B	C	D	E
Electronics for assay control and data collection				
National Instruments	NI USB-6351, X Series DAQ (16 AI, 24 DIO, 2 AO)	781440-01	\$1572.30	1
	Power cord, AC, 120 VAC	763000-01	\$9	1
Amazon	12V Power cable	https://www.amazon.com	\$9.39	1
/gp/product/B075VQ7NQH/ref=oh_a				
ui_detailpage_o05_s01?ie=UTF8&psc=1				
	Screwturn DC Power Jack Adapter	https://www.amazon.com	\$6.99	1
/gp/product/B0127J7NWE/ref=				
oh_aui_detailpage_o05_s00?ie=				
UTF8&psc=1				
	Male-Male jumper cables	-		
	Alligator Clips	https://www.amazon.com/gp	\$5.13	2
/product/B0002KRABU/ref=oh_				
au_detailpage_o02_s01?ie=UTF8&psc=1				
	Breadboard	https://www.amazon.com/gp/	\$10 for 3	1
product/B01EV6LJ7G/ref=oh_aui				
_detailpage_o09_s00?ie=UTF8&psc=1				
Digikey	220 and 330 Ohm Resistors			
	NPN Transistor	PN2222ATFCT-ND	\$0.18	1
	Diode	1N4007DICT-ND	\$1.30	1
	IR Emitter/Detector			

Diagram 1: NI Card connections summary:





Diagram 2: Interior of fully set up box



Behavior Hardware Setup

- 1** Use the spray adhesive to attach egg crate foam to the inside walls and ceiling of the Med-Associates box. Though the box is already sound attenuating, this will help with further sound attenuation.

- 2** Set up the NI Card:
 - 2.1** Place the NI USB-6351 card outside and adjacent to the Med-Associates box. This will allow for easy electrical connections, without potential damage from mice.

 - 2.2** Use Velcro or tape to securely attach a breadboard next to or on top of the NI card.

 - 2.3** Plug in the NI card to an outlet and the computer that will be running the MATLAB behavior code.

- 3** Build the treadmill:
 - 3.1** On the 12x12 Thor Labs base, about 4-5" in from one corner, use the clamping fork to attach the studded pedestal base adapter to the base board.

 - 3.2** Screw in the 1" post holder.

3.3 Attach the 2" optical post to the post holder.

3.4 Attach the swivel post clamp to the 2" optical post, and attach the 1" optical post to the other arm of the swivel clamp.

3.5 Screw this 1" optical post into the hole on the side of the rotary encoder.

3.6 Drill a small hole, the same size as the arm of the rotary encoder, into the direct center of the Plexiglas disk.

- Place the Plexiglas disk, such that it is tightly fitting and has no wiggle room, atop the arm of the rotary encoder.
- The rotary encoder arm should rotate when the Plexiglas disk is spun. If needed, use a 3D printed adapter to securely connect the two components.

3.7 Cut a small hole into the center of the silicon rubber sheet (to accommodate the rotary encoder arm), and place the silicon sheet on top of the Plexiglas disk to create a grippy surface for the treadmill.

3.8 Attach the 5-pin connector to the rotary encoder, with the other side leaving the box to connect to the NI card.

- Blue wire: to PFI8/P2.0
- Yellow wire: to PFI10/P2.2
- Orange wire: to +5V
- Brown wire: to ground
- Purple wire: unused

Note

The bag clips can be used as a treadmill “brake” for easier head fixation, during behavioral habituation, etc.

3.9 From these instructions, you should have a horizontal treadmill that can rotate with the rotary encoder and angle slightly in the horizontal plane, such that the mice are running slightly uphill.

4 Build the head-fixation apparatus:


4.1 In line with your treadmill, screw 2 7” stainless steel rods into the 12x12 base board on either side of the treadmill (4 total).

- Per side, place 1 rod as close to the treadmill as possible without touching it, and the other 2 peg holes away. Attach firmly to the base board with wing nuts.

4.2 Using hex nuts and washers, place the 3D printed head posts onto the stainless steel rods such that they meet over one side of the circular treadmill.

- The longer head post will cross the majority of the treadmill. Both head posts should be about 4-5” above the base board, therefore clearing the treadmill at the rough height of a mouse.


4.3 Use more washers and wing nuts to hold the head posts firmly in place. These can be loosened and tightened when head-fixing a mouse, and can be moved up and down the steel rods for height flexibility.




5 Build the infrared lick detection system:

- 5.1 On the 12x12 Thor Labs base, to the front left of the treadmill, use the clamping fork to attach the studded pedestal base adapter to the base board.
- 5.2 Screw in the 3" post holder.
- 5.3 Attach the 6" optical post to the post holder.

- 5.4 Attach the swivel post clamp to the 6" optical post, and attach the 4" optical post to the other arm of the swivel clamp.
- 5.5 Screw this 4" optical post into the small hole on the side of the 3D printed lick port.
- 5.6 The lick port has two additional holes at the front of the u-shape. Insert an IR emitter bulb into one, and an IR detector bulb into the other, such that they face each other.
- 5.7 IR detector: attach the long wire to ground, and the shorter wire to the NI card AI0 port and +5V using a 220Ω resistor (see below wiring diagram).



- 5.8** IR emitter: attach the short wire to ground, and the long wire to +5V using a 220Ω resistor (see below wiring diagram).




6 Build the reward delivery system:

- 6.1** Use glue to clog the innermost port of the 3 port solenoid.
- 6.2** Use Velcro to attach the solenoid, ports facing outwards, to the side of the Med-Associates Box.
- 6.3** Use Velcro to attach the 60mL Luer Lock syringe above the level of the lick port (this is your liquid reservoir; if it is above the lick spout, gravity will pull the liquid through to the lick spout)

whenever the solenoid is open).


- 6.4** Attach one side of the On/Off valve to the Luer Lock syringe, and the other to the stainless steel dispensing needle.
- 6.5** Using an appropriately cut length of PVC clear tubing, connect the dispensing needle of the reservoir to the outermost port of the solenoid.
- 6.6** Cut a longer length PVC clear tubing (needs to reach from solenoid to lick port) and attach the stainless steel tubing to one end. This will be the lick spout. Connect the other end of the PVC tubing to the middle port of the solenoid.
- 6.7** Put the stainless steel lick spout through the hole crossing through the middle of the “u” of the lick port, and use a thumbscrew through the top of the lick port to hold it in place.
- 6.8** Connect the solenoid to the NI card (see the below circuit diagram):
 1. Use a 12V power cable and the screwturn DC power jack adapter to create a 12V input to one of the solenoid connections.
 2. Wire the other solenoid connection to the +12V using a diode, and to the C arm of a NPN transistor.
 3. Wire the E arm of the transistor to ground.
 4. Wire the B arm of the transistor to a 330Ω resistor, and connect this to the AO1 output port on the NI card.
 5. Connect AO1 on the NI card to the AI2 input port, to read the solenoid output into the MATLAB code.



7 Set up the speaker:

- 7.1 Plug the speaker into an outlet, and place in the Med-Associates box in front of the treadmill.
- 7.2 Plug in one end of the audio cable to the speaker.
- 7.3 To have the MATLAB code control the speaker output, attach the other end of the audio cable to the NI card using alligator clips and male-male jumper cables.
 - Connect the audio cable ground to AO GND on the NI Card.

- Connect the audio cable input to both AI3 and AO0 on the NI Card.



8 Set up the Basler ace camera:

8.1 Connect the camera to the lens, tripod mount, and tripod.

8.2 Plug the camera in to an outlet and to the computer running the MATLAB code, and set in the Med-Associates box angled so that the mouse and lick port can be readily viewed.

8.3 Plug in the Bosch IR Illuminator, and place near the camera so that it illuminates the treadmill setup.

9 Connect P0.0 on the NI Card to AI9 to create an end of trial flag for the MATLAB code.

