



Jul 05, 2022

# RNAlater recipe

Yin-Tse Huang<sup>1</sup><sup>1</sup>KMU1 *Works for me* Share[dx.doi.org/10.17504/protocols.io.bp2l61w35vqe/v1](https://dx.doi.org/10.17504/protocols.io.bp2l61w35vqe/v1)

Yin-Tse Huang

## ABSTRACT

RNAlater recipe

## DOI

[dx.doi.org/10.17504/protocols.io.bp2l61w35vqe/v1](https://dx.doi.org/10.17504/protocols.io.bp2l61w35vqe/v1)

## PROTOCOL CITATION

Yin-Tse Huang 2022. RNAlater recipe. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bp2l61w35vqe/v1>



## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Jul 05, 2022

## LAST MODIFIED

Jul 05, 2022

## PROTOCOL INTEGER ID

65948

Before you start

1d

1 Submerge all your equipment and glassware in **1% bleach to get rid of nuclease**

Reagents

1d

2 **DEPC treated water**

1d

1. Add **1 mL** DEPC in **1 L** ddH<sub>2</sub>O (**0.1 % (v/v)**); **MIX it well.**
2. Let it sit for **24:00:00** at **25 °C**
3. Autoclave it to deactivate DEPC
4. Store the DEPC treated water at **4 °C**

### 3 **0.5 M Disodium dihydrate EDTA:**

1. Add **18.61 g** disodium dihydrate EDTA in **100 mL** ddH<sub>2</sub>O
2. Adjust pH to 8.0 with NaOH while stirring

### 4 **1 M tri-Sodium Citrate Dihydrate:**

Add **29.4 g** tri-Sodium Citrate Dihydrate in **100 mL** ddH<sub>2</sub>O, stir to dissolve

### 5 **1 M sulfuric acid:**

Add **5.33 mL** sulfuric acid in **100 mL** ddH<sub>2</sub>O, stir to dissolve

## Start to make RNAlater

### 6 Make a mixture of:

**40 mL** of **0.5 M EDTA** + **25 mL** **1 M tri-Sodium Citrate** + **700 g Ammonium Sulfate** + **935 mL ddH<sub>2</sub>O**; stir on a hot plate with low heat until ammonium sulfate dissolved

### 7 Allow it to cool; adjust the pH of the solution to **pH 5.2** using **1 M sulfuric acid**

To avoid contamination from the immersed pH meter, adjust pH by pouring 10% of the solution in a new beaker and add base or acid to hit the right pH; then scale it up to the needed amount for the original solution and blindly adjust its pH

### 8 The final concentration: 25 mM Sodium Citrate, 10 mM EDTA, 70 g ammonium sulfate in 100 ml ddH<sub>2</sub>O; pH 5.2.

### 9 Transfer to screw top, **nuclease-free** bottles; label the date of make; store at **4 °C**

