



Version 2

Oct 14, 2020

# BHI + v2 salts media V.2

Matthew Haines<sup>1</sup><sup>1</sup>Imperial College London

In Development

[dx.doi.org/10.17504/protocols.io.bmb8k2rw](https://dx.doi.org/10.17504/protocols.io.bmb8k2rw)**Matthew Haines**  
Imperial College London

## ABSTRACT

*Vibrio natriegens* grows exceptionally well in BHI + v2 salts media ([link](#)). Furthermore, this media is routinely used for culturing this organism (Weinstock et al., 2016). However, care must be taken during preparation not to autoclave v2 salts and BHI together. This protocol utilises a 10x v2 salts buffer to achieve this.



Weinstock MT, Hesek ED, Wilson CM, Gibson DG (2016).  
*Vibrio natriegens* as a fast-growing host for molecular  
biology.. Nature methods.  
<https://doi.org/10.1038/nmeth.3970>

## DOI

[dx.doi.org/10.17504/protocols.io.bmb8k2rw](https://dx.doi.org/10.17504/protocols.io.bmb8k2rw)

## PROTOCOL CITATION

Matthew Haines 2020. BHI + v2 salts media. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bmb8k2rw>  
Version created by [Matthew Haines](#)

## KEYWORDS

*Vibrio natriegens*, Model prokaryotes, Synthetic biology

## LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Sep 14, 2020

## LAST MODIFIED

Oct 14, 2020

## PROTOCOL INTEGER ID

42080

## MATERIALS

| NAME  | CATALOG # | VENDOR  |
|---|-----------|---------|
| <a href="#">MgCl2</a>                                 |           |         |
| <a href="#">Sodium chloride</a>                       |           |         |
| <a href="#">Potassium Chloride</a>                    |           |         |
| <a href="#">Brain Heart Infusion Broth Dry Medium</a> | B9500     | Teknova |

## ABSTRACT

*Vibrio natriegens* grows exceptionally well in BHI + v2 salts media ([link](#)). Furthermore, this media is routinely used for culturing this organism (Weinstock et al., 2016). However, care must be taken during preparation not to autoclave v2 salts and BHI together. This protocol utilises a 10x v2 salts buffer to achieve this.



Weinstock MT, Hesek ED, Wilson CM, Gibson DG (2016).  
*Vibrio natriegens* as a fast-growing host for molecular biology.. Nature methods.

<https://doi.org/10.1038/nmeth.3970>

### Prepare 10x v2 salts buffer

- 1 Dissolve the following salts in **800 mL** of ddH<sub>2</sub>O:

| Component         | Target concentration (mM) | Component (g/900 mL) | Final concentration (mM) |
|-------------------|---------------------------|----------------------|--------------------------|
| NaCl              | 2040.0                    | <b>107.30</b>        | 2040.08                  |
| MgCl <sub>2</sub> | 231.4                     | <b>19.83</b>         | 231.42                   |
| KCl               | 42.0                      | <b>2.82</b>          | 42.03                    |

- 2 Adjust the volume to **900 mL** using ddH<sub>2</sub>O.

### Prepare BHI media

- 3 Dissolve **37 g** of BHI dry medium in **900 mL** ddH<sub>2</sub>O.

### Sterilise and combine

- 4 Sterilise both BHI media and 10x v2 salts buffer by autoclaving.
- 5 Add **100 mL** 10x v2 salts buffer to the sterilised **900 mL** BHI media under sterile conditions.



The remaining 10x v2 salts buffer can be used for making further BHI + v2 salts media.