



# Standard PCR protocol V.1

Victoria LN Jackson<sup>1</sup>

<sup>1</sup>Living Systems Institute, University of Exeter

Version 1

Jun 28, 2021

1

Works for me



Share

This protocol is published without a DOI.

Victoria Jackson

Living Systems Institute, University of Exeter

## ABSTRACT

Generic PCR protocol

## PROTOCOL CITATION

Victoria LN Jackson 2021. Standard PCR protocol. **protocols.io**

<https://protocols.io/view/standard-pcr-protocol-bv54n88w>

Version created by Victoria Jackson

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Jun 28, 2021

## LAST MODIFIED

Jun 28, 2021

## PROTOCOL INTEGER ID

51100

Master mix 10m

1 For a 25 µl reaction volume use the following recipe:

14.9 µl dH<sub>2</sub>O

5 µl buffer

1.5 µl MgCl<sub>2</sub>

1 µl FWD primer (from a 10 Micromolar (µM) working solution)

1 µl REV primer (from a 10 Micromolar (µM) working solution)

0.5 µl dNTPs (from a 10 Micromolar (µM) mix)



0.1 µl Taq polymerase

2 Prepare the master mix **On ice** allowing for 10% pipetting error (eg. if carrying out 10 reactions, calculate the master mix quantities needed for 11 reactions).

Remember to include a **negative control** with no DNA template and, if possible, a **positive control** using a known DNA template.

## Preparing PCR reaction tubes

5m



- 3 Aliquot  **24 µl** of master mix into each 0.2 mL PCR tube.
- 4 Add  **1 µl** of either DNA template, dH<sub>2</sub>O (negative control) or known DNA (positive control) to the PCR tubes and place in a thermocycler.

## Thermocycling program

- 5 Standard PCR thermocycler program (30 cycles):

- |     |  |     |
|-----|--|-----|
| 5.1 | Initial denaturation at  <b>94 °C</b> for  <b>00:03:00</b>   | 3m  |
| 5.2 | Denature at  <b>94 °C</b> for  <b>00:00:30</b>   | 30s |
| 5.3 | Annealing at ~  <b>55 °C</b> (depending on primers) for  <b>00:01:00</b>   | 1m  |
| 5.4 | Extension at  <b>72 °C</b> for  <b>00:02:00</b><br> <b>x 29</b> | 2m  |
| 5.5 | Final extension at  <b>72 °C</b> for  <b>00:05:00</b>  | 5m  |

## Storage

- 6 Store PCR products at  **4 °C** , or at  **-20 °C** for long term storage.

## Analysis

- 7 Check if the PCR was successful by gel electrophoresis.







Gel electrophoresis  
by Victoria Jackson,  
Living Systems Institute, University of Exeter

PREVIEW

RUN










- 7.1 For  **200 mL** of a 1.2% agarose gel, add  **2.4 g** of agarose to  **200 mL** TAE or TBE in a conical flask and mix by swirling the flask.

- 7.2 Microwave on a high power for about  **00:03:00** , or until boiling, swirling the flask approximately every **00:00:30** during heating to ensure all the agarose dissolves. 3m 30s



Use heat proof gloves.

- 7.3 Allow the agarose to cool until you are able to comfortably touch the outside of the flask (approximately  **00:05:00** <sup>5m</sup> ).
- 7.4 Add  **10 µl** SafeView stain to the hot agarose and swirl to mix.
- 7.5 Pour the agarose into the casting tray. Pop any bubbles or move them to the edges of the mould using a pipette tip. Place the comb as desired.
- 7.6 Allow the gel to set for approximately  **00:40:00** . 40m
- 7.7 Place the gel (still in the casting tray) gently into the tank, ensuring there is enough TAE or TBE running buffer to cover it.
- 7.8 If the buffer used for PCR already contains loading dye, skip to step 11.
- 7.9 On a piece of parafilm, dot  **1 µl** of 6x loading dye using a pipette.
- 7.10 Pipette  **5 µl** ladder or  **6 µl** PCR product directly onto the loading dye and mix by pipetting up and down several times.
- 7.11 Load  **6 µl** into each well.
- 7.12 Place the lid on the tank, making sure that the wells are at the negative electrode end of the tank. The samples should run from black to red (negative to positive).
- 7.13 Set the power supply to between **80-120V**, depending on the expected amplicon size and desired duration of electrophoresis.

**7.14** Run the electrophoresis until the visible dye bands are about 75% of the way down the gel, approximately **1-1.5 hours**.



Once complete, ensure you disconnect the power supply from the gel tank, before removing the lid.