

Apr 22, 2024

Self-made chrome alum gelatin coated slides

DOI

dx.doi.org/10.17504/protocols.io.3byl49kkogo5/v1

Sonja Fritzsche¹

¹Max-Delbrueck-Center



Sonja

Sonja Fritzsche

Max-Delbrueck-Center

OPEN  ACCESS



DOI: **dx.doi.org/10.17504/protocols.io.3byl49kkogo5/v1**

Protocol Citation: Sonja Fritzsche 2024. Self-made chrome alum gelatin coated slides. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.3byl49kkogo5/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 22, 2024

Last Modified: April 22, 2024



Protocol Integer ID: 98589

Abstract

This is step-wise protocol to coating microscopy slides with chrom alum gelatin for better tissue adherence.








Materials

-  Chromium(III) potassium sulphate dodecahydrate (#3535.3 Carl Roth) **Carl Roth Catalog #3535.3**
-  Gelatine from porcine skin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #48722**







Prepare chrome alum gelatin solution

- 1 Heat up  200 mL ddH₂O to  50 °C on hot plate using a magnetic stirrer.
- 2 Add  1 g gelatin gradually until completely dissolved to prevent formation of clumps.
Solution has to be clear.
- 3 Let solution cool down to room temperature and add  0.1 g Chrome potassium sulfate until completely dissolved.
Either proceed immediately or store the solution for a few days @  4 °C . If stored, let the solution warm-up to room temperature before use.

Slide coating

10m 10s

- 4 Pour prepared chrome alum gelatin solution in a staining dish preheated at  60 °C .
- 5 Incubate to-be-treated slides in [M] 70 % (v/v) Ethanol for  00:10:00 to clean them. 10m
- 6 Rinse in distilled water, dip at least 3 times to remove alcohol.
- 7 Dip the slides 5 times in the solution,  00:00:10 each. 10s
- 8 Leave the slides to drain for a few minutes onto a filter paper.
- 9 Dry the treated slides in an oven at  37 °C over night.
- 10 Store the slides in a dry box and protect from dust.

Protocol references

The protocol is adapted based on the following published protocols:

<https://www.stainsfile.com/adhesives/chrome-alum-gelatin/>