

Aug 28, 2024

## Fluorescent Immunolabelling for endogenous ChAT in mice striatum

DOI

[dx.doi.org/10.17504/protocols.io.8epv59wk6g1b/v1](https://dx.doi.org/10.17504/protocols.io.8epv59wk6g1b/v1)

Roberta Marongiu<sup>1,2</sup>

<sup>1</sup>Department of Neurological Surgery, Weill Cornell Medical College, New York, NY 10065;

<sup>2</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA

ASAP Collaborative Rese...



Eileen Ruth Torres

Weill Cornell Medicine

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.8epv59wk6g1b/v1](https://dx.doi.org/10.17504/protocols.io.8epv59wk6g1b/v1)

**Protocol Citation:** Roberta Marongiu 2024. Fluorescent Immunolabelling for endogenous ChAT in mice striatum. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.8epv59wk6g1b/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 19, 2022

**Last Modified:** August 28, 2024

**Protocol Integer ID:** 68929

**Keywords:** ASAPCRN

**Funders Acknowledgement:**

**Aligning Science Across**

**Parkinson's Disease**

**Grant ID:** 020608



## Abstract

This immunohistochemistry protocol has been designed and adjusted specifically for ChAT immunolabelling in the mice striatum

## Materials

⊗ Goat anti-ChAT antibody **Millipore Catalog #AB144P**

⊗ Donkey anti-Goat IgG (H L) Cross-Adsorbed Secondary Antibody Alexa Fluor™ 488 **Thermo Fisher Scientific Catalog #A-11055**

⊗ DAPI (4',6-Diamidino-2-Phenylindole, Dihydrochloride) **Thermo Fisher Catalog #D1306**













⊗ Donkey serum **Sigma Aldrich Catalog #S30-100ML**

## Solution

Tris Buffer Soaline Triton TBST - 2L

- 50mM Tris pH 7.5 – 100 mL of 1M stock
- 150 mM NaCl – 17.53 g
- Distilled Water – 1.9 L
- Triton 0.2 % - 4 mL



- 1 Select floating sections of mouse striatum (30um thick).
- 2 Wash: 3 x  00:10:00 gently on a shaker, in TBST at  Room temperature . 10m
- 3 During the washes, prepare blocking solution: 10 % donkey serum (DS) in TBST.
- 4 Blocking step: Incubate the slides with blocking solution on the shaker for  01:00:00 at  Room temperature . 1h
- 5 Before the end of the blocking step, prepare the solution with the primary antibody ChAT blocking solution (TBST+10% DS) + Goat Anti-Chat Ab ♦ dilution 1:200.
- 6 Primary antibody: incubate the slides with the blocking solution/Antibody on a shaker at  4 °C for  24:00:00  48:00:00 . 3d
- 7 Washes: 3 x  00:10:00 on a shaker, in TBST. 10m
- 8 During the washes, prepare the secondary Antibody solution blocking solution (TBST+10% DS) + anti donkey anti goat Alexa Fluor 1:1000.
- 9 incubate the slides with Secondary antibody solution  01:00:00 at  Room temperature . 1h
- 10 Washes: 2 times, for  00:10:00 in TBST. 10m
- 11 *Optional: DAPI staining*  
During the washes, prepare the DAPI solution  
DAPI (1 mg/ml). Dilution 1:10.000 in TBST.
- 12 Incubation the slides  00:10:00 with the DAPI solution. 10m



13 Washes: 2 times, for  00:10:00 each in TBST.

10m

14 Mounting on Superfrost+ slides.

15 Coverslip: mounted with prolong anti-fade reagent (Invitrogen).