



Jan 12, 2021

Multiplexed RT-qPCR to screen for SARS-CoV-2 B.1.1.7 variants

Chantal Vogels¹, Joseph Fauver¹, Nathan Grubaugh¹¹Department of Epidemiology of Microbial Diseases, Yale School of Public Health*In Development* dx.doi.org/10.17504/protocols.io.brcrm2v6

Coronavirus Method Development Community

Nathan Grubaugh
Department of Epidemiology of Microbial Diseases, Yale School of Public Health

ABSTRACT

As the B.1.1.7 "UK" variant continues to spread, there is a need for labs to screen COVID-19 samples for these viruses. While sequencing is the gold standard, it cannot always be scaled or implemented in some settings. The Applied Biosystems TaqPath COVID-19 assay (ThermoFisher), a PCR test, was discovered to have a distinct signature (spike gene target failure, [SGTF]) when testing viruses containing the $\Delta 69/70$ HV deletion, like the B.1.1.7 variant. Detecting the $\Delta 69/70$ HV deletion alone is not definitive for the B.1.1.7 variant, but tracking the frequencies of SGTFs helped the UK track the B.1.1.7 variant. In the US and other countries, screening samples for the SGTF helped to identify potential B.1.1.7 variants for sequencing prioritization. ThermoFisher, however, has not released their spike probe sequence, so the assay needs to be recreated to be used more broadly. Here, we designed a multiplexed PCR assay that recreates the SGTF signature using primer/probe set that targets the $\Delta 69/70$ HV deletion, the CDC N1 set as a positive control, and the CDC RNase P set as an extraction/sample control. Samples positive for the CDC N1 set but negative for the $\Delta 69/70$ HV deletion set can be tracked and prioritized for sequencing.

DOI

dx.doi.org/10.17504/protocols.io.brcrm2v6

PROTOCOL CITATION

Chantal Vogels, Joseph Fauver, Nathan Grubaugh 2021. Multiplexed RT-qPCR to screen for SARS-CoV-2 B.1.1.7 variants. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.brcrm2v6>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jan 11, 2021

LAST MODIFIED

Jan 12, 2021

PROTOCOL INTEGER ID

46193

GUIDELINES

Disclaimer and intended use: This multiplexed protocol is still under development and is for research purposes only. It should not be used for clinical diagnosis. The intention of this assay is to screen for the probable presence of the B.1.1.7 (aka "501Y.V1 or "UK") variant lineage of SARS-CoV-2.

MATERIALS TEXT

- **NEB Luna® Universal Probe One-Step RT-qPCR Kit**

- > MM; Luna Universal Probe One-Step Reaction Mix, 2X

- > RT; Luna WarmStart® RT Enzyme Mix (20X)

- > Nuclease-free water

- **Primers/probes**

- > CDC_N1; Forward Primer (100 µM), Reverse primer (100 µM), probe (100 µM)

- > Yale_69/70del; Forward Primer (100 µM), Reverse primer (100 µM), probe (100 µM)

- > CDC_RP; Forward Primer (100 µM), Reverse primer (100 µM), probe (100 µM)

A	B	C	D	E
Set name	Nt positions	TM	Primer/probe	Sequence
CDC_N1	28,287	53.6	Forward primer	GACCCCAAAATCAGCGAAAT
	28,335	57.7	Reverse primer	TCTGGTTACTGCCAGTTGAATCTG
	28,309	63.3	Probe	FAM-ACCCCGCATTACGTTTGGTGGACC-BHQ1
Yale_69/70del	21,710-21,733	59.3	Forward primer	TCAACTCAGGACTTGTCTTACCT
	21,796-21,817	57.4	Reverse primer	TGGTAGGACAGGGTTATCAAC
	21,755-21,779	61.2	Probe (drop-out)	HEX-TTCCATGCTATACATGTCTCTGGGA-BHQ1
CDC_RP			Forward primer	AGATTTGGACCTGCGAGCG
			Reverse primer	GAGCGGCTGTCTCCACAAGT
			Probe	Cy5-TTCTGACCTGAAGGCTCTGCGCG-BHQ2

- **Positive control;** [Twist synthetic SARS-CoV-2 RNA controls](#) at 100 copies/uL

DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

As the B.1.1.7 "UK" variant continues to spread, there is a need for labs to screen COVID-19 samples for these viruses. While sequencing is the gold standard, it cannot always be scaled or implemented in some settings. The Applied Biosystems TaqPath COVID-19 assay (ThermoFisher), a PCR test, was discovered to have a distinct signature (spike gene target failure, [SGTF]) when testing viruses containing the Δ69/70 HV deletion, like the B.1.1.7 variant. Detecting the Δ69/70 HV deletion alone is not definitive for the B.1.1.7 variant, but tracking the frequencies of SGTFs helped the UK track the B.1.1.7 variant. In the US and other countries, screening samples for the SGTF helped to identify potential B.1.1.7 variants for sequencing prioritization. ThermoFisher, however, has not released their spike probe sequence, so the assay needs to be recreated to be used more broadly. Here, we designed a multiplexed PCR assay that recreates the SGTF signature using primer/probe set that targets the Δ69/70 HV deletion, the CDC N1 set as a positive control, and the CDC RNase P set as an extraction/sample control. Samples positive for the CDC N1 set but negative for the Δ69/70 HV deletion set can be tracked and prioritized for sequencing.

RT-qPCR Protocol

- 1 Briefly vortex and centrifuge reagents before use.

- 2 Prepare 20 μM working stocks of the primers and probes, by adding 20 μL of 100 μM stock to 80 μL nuclease-free water.
- 3 Use the 20 μM working stocks to prepare **primer-probe-water mix** containing the following:

A	B	C
Component	Volume (1 reaction)	Volume (100 reactions)
CDC_N1_F (400 nM/reaction)	0.4 μL	40 μL
CDC_N1_R (400 nM/reaction)	0.4 μL	40 μL
CDC_N1_P (200 nM/reaction)	0.2 μL	20 μL
Yale_69/70del_F (400 nM/reaction)	0.4 μL	40 μL
Yale_69/70del_R (400 nM/reaction)	0.4 μL	40 μL
Yale_69/70del_P (200 nM/reaction)	0.2 μL	20 μL
CDC_RP_F (200 nM/reaction)	0.2 μL	20 μL
CDC_RP_R (200 nM/reaction)	0.2 μL	20 μL
CDC_RP_P (100 nM/reaction)	0.1 μL	10 μL
Nuclease-free water	1.5 μL	150 μL

NOTE: a larger volume of primer-probe-water mix can be prepared in advance, aliquoted in LightSafe microcentrifuge tubes, and stored at -20°C .

- 4 Diagram sample, standard, and control positions on a 96-well plate map.
- 5 1. On ice, prepare a master mix containing the following (account for 10% extra lost during pipetting), except RNA:

A	B
Component	Volume in 20 μL reaction
Tube label = MM	10 μL
Tube label = RT	1 μL
Tube label = primer-probe-water mix	4 μL
Viral RNA, positive control, or negative control	5 μL (do not add to master mix)

- 6 Add 15 μL of mastermix to each well (on ice).
- 7 Add 5 μL of positive control (Twist 100 copies/ μL) and no-template control (NTC - water) to the designated wells (on ice). Mix by pipetting (avoid bubbles).
- 8 Add 5 μL of viral RNA to the designated wells (on ice). Mix by pipetting (avoid bubbles).

Cover with plate sealer. Centrifuge to remove bubbles, if present.

9

10 Set the thermocycler to read **FAM, HEX, and Cy5** fluorophores.

11 Run the following thermocycler conditions:

A	B	C
Step	Temperature	Time
1	55°C	10 min
2	95°C	1 min
3	95°C	10 sec
4	55°C	30 sec
5	Read plate	
Repeat steps 3-5 for 44 cycles.		

12 Interpreting results:

A	B	C	D
Result	CDC_N1	Yale_69/70del	CDC_RP
Potentially B.1.1.7 (or other B.1 variant)	CT ≤ 40	CT > 40 or undetected	Any value
Not B.1.1.7	CT ≤ 40	CT ≤ 40	Any value
Invalid	CT > 40 or undetected	CT > 40 or undetected	CT > 40 or undetected
Negative	CT > 40 or undetected	CT > 40 or undetected	CT < 40