


Sep 09, 2022

Preparation of LRRK2 RCKW trimer cryo-EM grids

Mariusz Matyszewski¹¹Department of Cellular and Molecular Medicine, University of California, San Diego, La Jolla, CA 92093

1 Works for me

 Share

This protocol is published without a DOI.

 Mariusz Matyszewski

ABSTRACT

This protocol has been adapted from Deniston et al (<https://doi.org/10.1038/s41586-020-2673-2>)

Original protocol by Colin Deniston. Adapted to protocols.io by Mariusz Matyszewski.

This protocol describes how to create cryo-EM grids for LRRK2 RCKW. In particular, this protocol was used to obtain a high-resolution cryo-EM structure of the LRRK2 RCKW trimer, as well as lower resolution structures of RCKW monomers and dimers.

PROTOCOL CITATION

Mariusz Matyszewski 2022. Preparation of LRRK2 RCKW trimer cryo-EM grids.
protocols.io
<https://protocols.io/view/preparation-of-lrrk2-rckw-trimer-cryo-em-grids-brypm7vn>



KEYWORDS

cryo-EM, LRRK2, structural biology, ASAPCRN

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jan 31, 2021

LAST MODIFIED

Sep 09, 2022

PROTOCOL INTEGER ID

46831

GUIDELINES

Newer protocol is available (see [updated protocol here](#)) with better results for monomeric protein. This is included mainly for archival reasons.

MATERIALS TEXT

LRRK2 Buffer:

- [M]**20 millimolar (mM)** HEPES pH 7.4
- [M]**80 millimolar (mM)** NaCl
- [M]**0.5 millimolar (mM)** TCEP
- [M]**5 % volume** Glycerol
- [M]**2.5 millimolar (mM)** MgCl₂
- [M]**20 micromolar (μM)** GDP

SAFETY WARNINGS

For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).

Take proper precautions while freezing grids.

Preparing Sample

- 1 Dialyze purified **LRRK2 RCKW** into the final LRRK2 buffer (see Materials)
- 2 Dilute the protein to the desired final concentration.

High concentrations favor dimers and trimers, while lower concentrations favor monomers. For the published trimer structure, a final concentration of 4 μM was used (see Deniston et al).

- 2.1 **(Optional)** If adding inhibitors, add them after diluting to the final concentration and incubate on ice for 🕒 **01:00:00** . If using Apo protein, then proceed immediately. 1h

Freezing Grids

20s

- 3 Glow discharge grids. 20s
We used UltrAuFoil Holey Gold 1.2/1.3 300 mesh grids and glow discharged them at 20 mA for 🕒 **00:00:20** in a K100 instrument.
- 4 Apply protein to grids and plunge freeze.
We used a Vitrobot (FEI) to blot away excess sample and plunge freeze
- 5 Store grids in liquid nitrogen until ready for imaging.