

May 22, 2024

## Digestion with NEBNext dsDNA Fragmentase (M0348)

 Forked from [Digestion with NEBNext dsDNA Fragmentase \(M0348\)](#)

DOI

[dx.doi.org/10.17504/protocols.io.j8nlk8w75l5r/v1](https://dx.doi.org/10.17504/protocols.io.j8nlk8w75l5r/v1)

New England Biolabs<sup>1</sup>

<sup>1</sup>New England Biolabs

New England Biolabs (NEB)

Tech. support phone: +1(800)632-7799 email: [info@neb.com](mailto:info@neb.com)



kassidy hebert

Rice University



OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.j8nlk8w75l5r/v1](https://dx.doi.org/10.17504/protocols.io.j8nlk8w75l5r/v1)

External link: <https://www.neb.com/protocols/0001/01/01/digestion-with-nebnext-dsdna-fragmentase-m0348>

**Protocol Citation:** New England Biolabs 2024. Digestion with NEBNext dsDNA Fragmentase (M0348). [protocols.io](https://dx.doi.org/10.17504/protocols.io.j8nlk8w75l5r/v1)  
<https://dx.doi.org/10.17504/protocols.io.j8nlk8w75l5r/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

**We are still developing and  
optimizing this protocol**

**Created:** February 09, 2020

**Last Modified:** May 22, 2024

**Protocol Integer ID:** 100201

**Keywords:** M0348, dsDNA, Fragmentation, Fragmentase, Shear, DNA, NGS,



## Abstract

NEBNext dsDNA Fragmentase is an enzyme-based reagent that shears DNA to produce fragments of the desired sizes in a time-dependent manner, for next generation sequencing library preparation protocols

- dsDNA Fragmentase provides random fragmentation, similar to mechanical methods (1,2).

## Materials

### MATERIALS


⊗ NEBNext dsDNA Fragmentase - 250 rxns **New England Biolabs Catalog #M0348L**


⊗ NEBNext dsDNA Fragmentase - 50 rxns **New England Biolabs Catalog #M0348S**

## Safety warnings

❗ Please refer to the Safety Data Sheets (SDS) for health and environmental hazards.

## Before start

Adequate mixing of NEBNext dsDNA Fragmentase is important for the success of this reaction. NEBNext dsDNA Fragmentase should be vortexed for  00:00:03 **prior to use.**



For tough digestions, add  1  $\mu\text{L}$  of [M] 200 millimolar (mM)  $\text{MgCl}_2$  to the reaction. Additional  $\text{MgCl}_2$  can be added if necessary.

The protocol listed below is for fragmentation of **5 ng–3  $\mu\text{g}$**  of DNA.

1 Vortex NEBNext dsDNA Fragmentase for  00:00:03 , quick spin and place  On ice .


2 Combine the following components in a sterile PCR tube and vortex:

A	B
Component	Amount
DNA (5 ng–3 µg)	1–16 µl
10X Fragmentase Reaction Buffer v2	2 µl
Sterile Water	variable
<b>Final Volume</b>	<b>18 µl</b>

3 Add  2 µL dsDNA Fragmentase and vortex mixture for  00:00:03 .

#### Note

Fragmentase is very viscous and should be pipetted slowly. If the enzyme has been sitting for several minutes vortex it again before adding to the sample.

4 Incubate at  37 °C for the recommended times below to generate the desired fragment size:

#### Note


If starting material is 100 ng or less, incubation times should be increased by 10 minutes.

A	B
Desired Fragment Size (bp)	Incubation Time (min)
50–200	25–35
200–1,000	15–25
1,000–2,000	10–15

#### Note

To determine the exact incubation time for a given sample type, a time course study should be performed.



5 Add  5  $\mu$ L 0.5 M EDTA to stop the reaction.



6 Clean up the fragmented DNA with column purification or using SPRI beads.

#### Note

If using SPRI beads, it is recommended to dilute the sample 1:1 with sterile water for easier handling of the sample and faster collection of the beads to the magnet.

SPRI beads are available from Beckman Coulter: A63880, A63881, A63882

#### For further analysis:

**Bioanalyzer:** Clean up the fragmented DNA prior to loading on a Bioanalyzer chip.

**End Repair:** Clean up the fragmented DNA then proceed with desired DNA end repair protocol.

**Polyacrylamide Gel Analysis:** Clean up the fragmented DNA prior to loading the samples on a PAGE gel.

**Long Term Storage:** Clean up the fragmented DNA prior to long term storage.

**Agarose Gel Size Selection/Analysis:** Samples can be loaded directly on to an agarose gel. It is not necessary to clean up the reactions prior to loading.