




Upload image

Jul 22, 2020

Autofluorescence Reduction and Imaging Marmoset NHP Tissue with a TissueFAXS Imaging System

Guoping Feng Lab¹¹Massachusetts Institute of Technology*In Development*

This protocol is published without a DOI.

 Guoping Feng Lab

PROTOCOL CITATION

Guoping Feng Lab 2020. Autofluorescence Reduction and Imaging Marmoset NHP Tissue with a TissueFAXS Imaging System . **protocols.io**
<https://protocols.io/view/autofluorescence-reduction-and-imaging-marmoset-nh-biunkeve>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 21, 2020






LAST MODIFIED

Jul 22, 2020

PROTOCOL INTEGER ID

39534

Autofluorescence Reduction 10m

- 1) Wash samples in 1X PBS for  **00:02:00** .
- 2) Stain samples free floating with  **0.05 % volume** TrueBlack in  **70 % volume** Ethanol for  **00:03:30**
- 3) Rinse samples in 1X PBS for  **00:10:00** .
- 4) Mount samples on a slide with ProLong™ Gold Antifade Mountan.

Imaging Samples 1w

- 5 Image samples with a 5x overview scan to determine the tissue regions and boarders.
- 6 Switch to 20x and image both slices with DAPI and 488 to determine the location of GFP+ cells.

- 7 Create "Defined Regions" around GFP+ cells.
- 8 Once these defined regions are created, switch to the 40x water immersion objective and "Force Focus" on multiple points across the tissue to find upper and lower bounds for focusing.
- 9 Compute focus points for each defined region -- Options, Focus (1x1) -- and set your upper and lower bounds for imaging.
- 10 Image each defined region with $\rightarrow 0.5 \mu\text{m}$ steps so that cells can be reconstructed for morphology.
- 11 Export the FOVs for morphology reconstruction. Be sure to uncheck "Use Stich" as the stitching is done later in the reconstruction pipeline. Save the Z-stack with "Original" and "16-bit" format.