

Sep 02, 2020

Agarose pads for microscopy

Joao Vitor Molino¹¹University of Zürich

1

Works for me

This protocol is published without a DOI.

Joao Vitor Molino
University of Zürich

ABSTRACT

This protocols describe the steps required for the preparation of agarose pads for live cell fluorescent microscopy.

PROTOCOL CITATION

Joao Vitor Molino 2020. Agarose pads for microscopy. **protocols.io**
<https://protocols.io/view/agarose-pads-for-microscopy-bkn8kvhw>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 02, 2020

LAST MODIFIED

Sep 02, 2020

PROTOCOL INTEGER ID

41408

GUIDELINES

All steps described in this protocol are intended to be conducted in a research laboratory. Follow aseptic procedures.

MATERIALS

NAME	CATALOG #	VENDOR
Frame-Seal™ in situ PCR and Hybridization Slide Chambers 15 x 15 mm 65 µl	SLF0601	

SAFETY WARNINGS

Follow the biosafety guidelines for the cell lines use during the experiment.

DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

BEFORE STARTING

Separate all material needed for the protocol.


Read through the protocol.

Before start

10m

1 Check all the necessary material:

- Glass slide
- Cover slip
- Agarose (*Agar should work as well, but probably with less imaging quality*)
- Media or buffer
- Chamber or well forming stick (ex: Frame-Seal™ in situ PCR and Hybridization Slide Chambers, 15 x 15 mm, 65 µl)
- Tube or flask for melting agarose
- 1 mL plastic tips



Melting agarose

5m

1. Melt the agarose in the desired concentration (**1.5 Mass / % volume** were tested) using a microwave, with interval mixing. Using the desired liquid. (e.g. *Media for keeping the cells alive*)
2. Keep the melted agarose in a plate heater at **70 °C**.



Glass slide preparation

10m

- 3
 1. Cut a 1mL tips to increase the point diameter for easier liquid handling of the agarose.
 2. Separate a glass slide to use as cover during pad preparation
 3. Stick the chamber forming material to a glass slide, maintaining the plastic barrier in the exposed side (*Alternatively you can follow the instructions at [dx.doi.org/10.17504/protocols.io.kfrctm6](https://doi.org/10.17504/protocols.io.kfrctm6) for a option without the chamber sticks*)
 4. Add a excess of melted agarose to the well, and cover it with another glass slide for shaping the agarose pad
 5. Remove the glass slide, by sliding it carefully to the side




Cutted plastic tip



Melted agarose added






Solidifying agarose




Glass slide removed after solidification

Sample addition

10m

- 4
 1. Add the desired amount of sample to the agarose pad. (5-6  1 μ l drops were tested | Be aware of cell concentration if the goal is to see individual cells)
 2. Wait the drops to dry  00:05:00 to  00:10:00 (*Specially important for cells with motility*)
 3. Carefully remove the plastic barrier exposing the adhesive face of the well forming material
 4. Add a cover slip, sealing the well. (*Avoids dehydration of the agarose pad, important for long microscopy experiments*)



Samples added to agarose pad



Glass slide ready for microscopy