

VERSION 2

FEB 12, 2024

🌐 Extraction and selection of high-molecular-weight DNA for long-read sequencing from *Chlamydomonas reinhardtii* V.2

📖 PLOS One ✓ Peer-reviewed method

Frédéric Chaux-Jukic¹, Nicolas Agier¹, Stephan Eberhard², Zhou Xu¹

¹Sorbonne Université, CNRS, UMR7238, Institut de Biologie Paris-Seine, Laboratory of Computational and Quantitative Biology, 75005 Paris, France;

²Sorbonne Université, CNRS, UMR7141, Institut de Biologie Physico-Chimique, Laboratory of Chloroplast Biology and Light-Sensing in Microalgae, 75005 Paris, France

PLOS ONE Lab Protocols

Tech. support email: plosone@plos.org

OPEN  ACCESS



zhou.xu

DISCLAIMER

The authors declare no conflict of interest.

DOI:

[dx.doi.org/10.17504/protocols.io.8epv59j9jg1b/v2](https://doi.org/10.17504/protocols.io.8epv59j9jg1b/v2)

External link:

<https://doi.org/10.1371/journal.pone.0297014>

Protocol Citation: Frédéric Chaux-Jukic, Nicolas Agier, Stephan Eberhard, Zhou Xu 2024. Extraction and selection of high-molecular-weight DNA for long-read sequencing from *Chlamydomonas reinhardtii*.

protocols.io

<https://dx.doi.org/10.17504/protocols.io.8epv59j9jg1b/v2> Version created by [zhou.xu](#)

ABSTRACT

Recent advances in long-read sequencing technologies have enabled the complete assembly of eukaryotic genomes from telomere to telomere by allowing repeated regions to be fully sequenced and assembled, thus filling the gaps left by previous short-read sequencing methods. Furthermore, long-read sequencing can also help characterizing structural variants, with applications in the fields of genome evolution or cancer genomics. For many organisms, the main bottleneck is to develop robust methods to obtain high-molecular-weight (HMW) DNA for whole genome sequencing purposes. We developed an optimized protocol to extract DNA suitable for long-read sequencing from the unicellular green alga *Chlamydomonas reinhardtii*, based on CTAB/phenol extraction followed by a size selection step for long DNA molecules. We provide validation results for the extraction protocol, as well as statistics obtained with Oxford Nanopore Technologies sequencing.

IMAGE ATTRIBUTION

Frédéric Chaux-Jukic

MANUSCRIPT CITATION:

Chaux F, Agier N, Eberhard S, Xu Z (2024) Extraction and selection of high-molecular-weight DNA for long-read sequencing from *Chlamydomonas reinhardtii*. PLOS ONE 19(2): e0297014. <https://doi.org/10.1371/journal.pone.0297014>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: May 20, 2022

Last Modified: May 20, 2022

PROTOCOL integer ID: 62922

Keywords: high-molecular weight DNA, long-read sequencing, size selection, algae DNA extraction

GUIDELINES

- This protocol aims at extracting high-molecular-weight DNA molecules with maximal removal of contaminants. Both parameters are indeed crucial to take full advantage of long-read sequencing.
- This protocol was developed for the microalga *Chlamydomonas reinhardtii*, but might be adapted to other algae and land plants.
- This protocol is relatively simple and short, however, for basic molecular biology (PCR, cloning of small fragments, etc.), simpler protocols may be preferred (e.g. www.chlamycollection.org/methods/).
- The extraction steps rely mostly on home-made solutions. Purification and size selection uses commercial reagents, the composition of which is unknown to us.

MATERIALS

- Equipment for *C. reinhardtii* culture
- Growth media: <https://www.chlamycollection.org/methods/media-recipes/tap-and-tris-minimal/>
- Ice bucket
- Water bath
- 2-ml tubes
- Wide-bore tips (or cut out the extremity of regular tips)
- 1.5-ml Lo-Bind tubes (Eppendorf)
- SpeedVac (optional)
- Magnetic rack

- Ultrapure water
- CTAB solution:
 - 50 mM Tris-HCl pH 8
 - 20 mM EDTA
 - 1.4 M NaCl
 - 2% CTAB (Hexadecyltrimethylammonium bromide; Sigma-Aldrich)
 - 1% PVP 40.000 (Polyvinylpyrrolidone; Sigma-Aldrich)
- Proteinase K: prepare stock solution at 20 mg.mL⁻¹
- RNase A: prepare stock solution at 100 mg.mL⁻¹
- Phenol:Chloroform:Isoamyl alcohol (25:24:1)
- Chloroform:Isoamyl alcohol (24:1)
- Isopropanol
- Ethanol: prepare at 70% v/v in ultrapure water

- AMPure XP beads (A63880, Beckman Coulter)
- Short Read Eliminator (SRE) kit (Circulomics)

- Optional: Qubit fluorimeter, Nanodrop
- Optional: pulsed-field gel electrophoresis (PFGE) system

SAFETY WARNINGS



Wear gloves at all steps of the protocol and apply standard safety procedure for biochemistry.

Handle phenol chloroform isoamyl alcohol and chloroform isoamyl alcohol under a chemical hood.

BEFORE START INSTRUCTIONS

On the day of extraction:

- Set water bath to 65°C (large enough for one 50-ml tube per sample)
- Chill isopropanol to -20°C (around 2.5 ml per sample)

Harvesting and storage of cells





- 1 Grow *Chlamydomonas reinhardtii* cells in 100 mL TAP medium under low light ($\sim 5 \mu\text{mol photon.m}^{-2}.\text{s}^{-1}$) with constant shaking at 100 rpm, 25°C, or in other conditions as appropriate for the specific experiment.
- 2 Harvest the cells at the end of exponential growth phase ($\sim 10^7 \text{ cells.mL}^{-1}$) in 50 mL tubes by centrifugation, 4000 x g, Room temperature, 00:05:00. Discard supernatant. 5m
- 3 Store cell pellets at -20 °C (optional).


DNA extraction

- 4 Thaw cell pellets at Room temperature (not applicable if the pellets were not frozen) and put on ice.
- 5 Centrifuge at 4000 x g, 4°C, 00:05:00 and discard any liquid left. 5m
- 6 Resuspend by gentle pipetting in 3 mL of CTAB solution preheated at 65 °C.

Note



To reduce DNA shearing, from this step forward, use wide bore tips (or cut standard tips) and pipet as slowly as possible, mix gently by inverting when needed and do not vortex.


- 7 Add  5 μL of proteinase K (stock solution: 20 mg.mL⁻¹) and  5 μL of RNase A (stock solution: 10 mg.mL⁻¹), mix by pipetting and incubate at  65 °C for  00:20:00 in a water bath. 20m

- 8 Centrifuge at  4000 x g, Room temperature, 00:05:00 and distribute 3 x 1 mL of each supernatant into 5m mL tubes. Discard pellet.

Note

Efficient lysis produces a pale green pellet, of much lighter hue than the supernatant.





- 9 Add  1 mL of phenol:chloroform:isoamyl alcohol (25:24:1) to each tube, mix gently by inverting 10 times and centrifuge at  20000 x g, 4°C, 00:05:00. Transfer the aqueous (upper) phase to new 2 mL tubes. 5m



After centrifugation, the aqueous phase should be colorless, the phenol (lower) phase should be dark green and a white precipitate may form at the interface.



Note


When pipetting the upper phase, avoid contact with the white precipitate that usually forms at the interface. To avoid traces of phenol, do not attempt to transfer more than three quarters of the aqueous phase.

- 10 Add  1 μ L of proteinase K and  1 μ L of RNase A, mix gently by inverting 3 times and incubate at  50 $^{\circ}$ C for  00:20:00 . 20m

Note




This step aims at completing RNA degradation and removal of DNA-binding proteins.





- 11 Add  1 mL of chloroform:isoamyl alcohol (24:1) to each tube, mix gently by inverting 3 times and centrifuge at  20000 x g, 4 $^{\circ}$ C, 00:05:00 . Transfer the upper phase to new 2 mL tubes. 5m

- 12 Add  700 μ L of isopropanol to each tube and mix gently by inverting 10-12 times. A visible DNA precipitates should form.

Note




If greenish/brownish contaminants precipitate together with DNA, do not use it for sequencing.

- 13 Centrifuge at  20000 x g, 4 $^{\circ}$ C, 00:05:00 , discard the supernatant carefully, add  500 μ L of ice-cold 10m
70% ethanol without disturbing the pellet and centrifuge at  20000 x g, 4 $^{\circ}$ C, 00:05:00 .



- 14 Discard the ethanol carefully by pipetting and eliminate the remaining traces using a SpeedVac for 35m
 00:05:00 at  Room temperature or leave tubes open for  00:30:00 at  Room temperature






- 15 Resuspend the pellet in  30 μL of ultrapure water overnight at  Room temperature .



Purification and selection of high-molecular-weight DNA

- 16 Pool the samples corresponding to the same initial culture in one 1.5ml Lo-Bind tube, add an equal volume  5m (e.g. 3 x  30 μL) of AMPure XP beads, mix gently by hand for  00:05:00 .




Note

The beads solution is relatively viscous, let it reach  Room temperature , then resuspend by vortexing for  00:00:10 before pipetting from the stock and ensure complete mixing with the DNA sample.

- 17 Place the tubes on a magnetic rack and let the beads aggregate on the side of the tube for  00:05:00  20m
Discard the supernatant. Remove the tubes from the magnetic rack, resuspend the beads in  60 μL of ultrapure water for at least  00:15:00 at  Room temperature .

- 18 Place the tubes on a magnetic rack and let the beads aggregate on the side of the tube for  00:05:00  5m
Transfer the supernatant to new 1.5 mL Lo-Bind tubes and discard the beads.









- 19 Measure DNA concentration and purity using a Nanodrop and/or a Qubit device.

- 20 Adjust each sample to $\sim 150 \text{ ng} \cdot \mu\text{L}^{-1}$ in  60 μL of ultra-pure water, gently mix with an equal volume of  30m
Short Read Eliminator (SRE) buffer (SRE kit, Circulomics) and centrifuge at  10000 x g, 4°C, 00:30:00 .

Pipet out the supernatant very carefully as the pellet is colorless and can easily be lost.

Note

Manufacturer's recommendation is to measure input DNA concentration using Qubit ([Short Read Eliminator Kit - Handbook v2.0](#))

- 21 Gently add  200 µL of freshly prepared 70% ethanol to the pellet. Wash by centrifuging at  10000 x g, 4°C, 00:02:00 and discarding the ethanol. 2m
- 22 Repeat Step 21.
- 23 Air dry the pellet for  00:10:00 at  Room temperature and resuspend in  50 µL of pre-warmed 30m
EB buffer (SRE kit, Circulomics) by incubating at  50 °C for  00:20:00 .
- 24 Proceed to long-read sequencing or store at  4 °C for later use (at least several weeks). Do not freeze to avoid DNA shearing.

Quality check

- 25 DNA purity can be assessed by Qubit or Nanodrop.
- 26 Absence of small fragments can be assessed by electrophoresis on agarose gel and size distribution can be assessed by pulsed-field gel electrophoresis (PFGE).

