



JAN 24, 2024

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.q26g7pk1qgwz/v1

Protocol Citation: Marta Gonzalez-Sepulveda, Thais Cuadros, Miquel Vila 2024. Protocol for preparing post-mortem tissue for intracellular neuromelanin quantification in H&E-stained sections..

protocols.io

<https://dx.doi.org/10.17504/protocols.io.q26g7pk1qgwz/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

🌐 Protocol for preparing post-mortem tissue for intracellular neuromelanin quantification in H&E-stained sections.

Marta Gonzalez-Sepulveda^{1,2}, Thais Cuadros^{1,2}, Miquel Vila^{1,3,4,5}

¹Neurodegenerative Diseases Research Group, Vall d'Hebron Research Institute (VHIR)-Network Center for Biomedical Research in Neurodegenerative Diseases (CIBERNED), 08035 Barcelona, Spain;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815;

³Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815;

⁴Institute of Neurosciences, Autonomous University of Barcelona (INc-UAB), Bellaterra, Barcelona, Spain;

⁵Catalan Institution for Research and Advanced Studies (ICREA), Barcelona, Spain

Vilalab Public



Miquel Vila

ABSTRACT

Protocol for preparing post-mortem tissue for intracellular neuromelanin quantification in H&E-stained sections

Protocol status: Working

We use this protocol and it's working

Created: Jan 24, 2024

Last Modified: Jan 24, 2024

PROTOCOL integer ID: 94079

Funders Acknowledgement:

Aligning Science Across
Parkinson's (ASAP)

Grant ID: ASAP-020505

Tissue Processing

- 1 Formalin-fixed paraffin-embedded tissue blocks from the pontine and midbrain regions were sectioned at 5µm, collected onto adhesive microscope slides and allowed to dry at room temperature for 24 hours.
- 2 Slides were incubated in the oven at 60°C for 30 minutes to melt the paraffin.
- 3 To remove the paraffin, sections were submerged in Xylene (Panreac Applichem 211769.2714) for 3 x 3 minutes, followed by rehydration in decreasing ethanol concentrations (100% ethanol for 2 x 5 minutes, 95% ethanol for 2 x 5, 70% ethanol for 2x5 minutes) and distilled H₂O for 5 minutes.
- 4 Standard hematoxylin-eosin (H&E) staining.

Image Acquisition

- 5 Sections were scanned using 20x objective (NA=0.8) with pre-set focusing and exposure parameters for optimal neuromelanin signal quality with an automated Slide Scanner Olympus (SLIDEVIEW VS200, Tokyo,

Japan).

- 6** Identical parameters were applied to each scanned section to ensure consistency in capturing neuromelanin.