



Apr 16, 2021

Single-nuclei suspensions from primary human esophagus tissue

Lucy Kimbley¹, Rachel Parker¹, Jack Harrington¹, Rob Walker¹, Annette Hayden¹, Jonathan West¹, Peter Johnson¹, Tim Underwood¹, Matthew Rose-Zerilli¹

¹Cancer Sciences Unit, University of Southampton

1 Works for me dx.doi.org/10.17504/protocols.io.t9wer7e

CZI START Project



Matthew Rose Zerilli
Cancer Sciences, Faculty of Medicine, University of Southamp...

ABSTRACT

A protocol to dissociate frozen esophagus tissue specimens for single-nuclei transcriptomics.

EXTERNAL LINK

<https://www.southampton.ac.uk/medicine/about/staff/tju.page>

DOI

dx.doi.org/10.17504/protocols.io.t9wer7e

EXTERNAL LINK

<https://www.southampton.ac.uk/medicine/about/staff/tju.page>

PROTOCOL CITATION

Lucy Kimbley, Rachel Parker, Jack Harrington, Rob Walker, Annette Hayden, Jonathan West, Peter Johnson, Tim Underwood, Matthew Rose-Zerilli 2021. Single-nuclei suspensions from primary human esophagus tissue. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.t9wer7e>

KEYWORDS

esophagus, esophageal, oesophageal, cancer, tissue, dissociation, single-nuclei, single-nucleus, RNA, sequencing

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 03, 2018

LAST MODIFIED

Apr 16, 2021

PROTOCOL INTEGER ID

16406

MATERIALS TEXT

MATERIALS

[BSA-Molecular Biology Grade - 12 mg](#) **New England**

Biolabs Catalog #B9000S

[tweezers](#) **Contributed by users**

[Proteases inhibitor](#)

Complete Roche Catalog #11 697 498

[4,6-Diamidino-2-Phenylindole, Dihydrochloride \(DAPI\)](#) **Thermo Fisher**

Scientific Catalog #D1306

[Nuclei EZ lysis buffer](#)

Sigma Catalog #EZ PREP NUC-101

[Dounce homogenizers](#)

Sigma Catalog #D8938-1SET

[Falcon 15 mL Polystyrene Conical Tube](#) **Fisher**

Scientific Catalog #352095

[NxGen RNase](#)

Inhibitor Lucigen Catalog #30281-2

[Corning® 40µm Cell](#)

Strainer Corning Catalog #431750

[DNA LoBind Tube 1.5ml](#)

Eppendorf Catalog #022431021

[DTT Sigma](#)

Aldrich Catalog #43816-10ML

[Trypan Blue Solution 0.4% Sterile-filtered](#) **Sigma**

Aldrich Catalog #T8154

[Aluminum Foil](#) **Contributed by users**

[Surgical scissors or scalpel](#) **Contributed by users**

[RNaseZap™ RNase Decontamination Solution](#) **Thermo Fisher**

Scientific Catalog #AM9780

[Dulbecco's Phosphate Buffered Saline \(PBS\)](#) **Contributed by**

users Catalog #D8537

[Polypropylene Centrifuge Tubes, 50mL Conical bottom tube; Max RCF: 15500xG; High-density polyethylene cap; No rack](#) **Thermo**

Fisher Catalog #0553868

[Sterilin™ 30 to 140mm Petri Dishes, 60mm](#) **Thermo**

Fisher Catalog #123-17

[Benchtop Dewar Flasks, 1L](#) **Thermo**

Fisher Catalog #4150-1000PK

[RNAlater™-ICE Frozen Tissue Transition Solution](#) **Thermo**

Fisher Catalog #4427575

[Ambion™ RNase Inhibitor, cloned, 40 U/µL](#) **Thermo**

Fisher Catalog #AM2684

[30 µM Cell Strainer](#) **Miltenyi**

Biotec Catalog #130-098-458

BEFORE STARTING

RNAlater-ICE pre-treatment of frozen tissue (**Important:** perform on at least the day before nuclei isolation):

1. Pre-chill a minimum of 10 volumes (compared to the sample mass/ volume) of RNAlater-ICE at 70 or °80°C in a polypropylene tube.
2. Place frozen tissue (from LN2 or °80°C), no larger than 0.5 cm in the smallest dimension, into the RNAlater-ICE. It is important that the RNAlater-ICE stays cold and the tissue remains frozen, so work quickly. Set up the tubes

- of RNAlater-ICE and tissue samples in a container of dry ice. For the frozen tissue, simply pipet at least 10 volumes of pre-chilled RNAlater-ICE into the tube. Tightly cap the tube and invert several times to mix.
3. Soak the tissue in RNAlater-ICE at -20°C for at least 16 hours. Samples transitioned to -20°C in RNAlater-ICE should be left in the RNAlater-ICE and maintain at -20°C or colder.

Reagents to prepare immediately before use:

- Prepare 6 mL of Nuclei Suspension Buffer (NSB; consisting of 1x PBS, 0.01% BSA and 0.2% v/v RNase inhibitor, 1 mM DTT). Keep cold on ice.
- Prepare 7 mL of Nuclei EZ lysis buffer (0.2% v/v RNase inhibitor, 0.2% v/v SUPERasin, 1x Protease inhibitor and 4 mM DTT). Keep cold on ice.
- Preparing the 300 nM DAPI solution:
 - i) Add 2 mL of ddH₂O to the DAPI vial to make a 14.3 mM (5 mg/mL) stock solution. Store at -20°C .
 - ii) Add 2.1 μL of the 14.3 mM DAPI stock solution to 100 μL PBS to make a 300 μM DAPI intermediate dilution.
 - iii) Dilute the 300 μM DAPI intermediate dilution 1:1000 in PBS as needed to make a 300 nM DAPI solution.

NB: Clean benchtop area with RNA-Zap and wipe down pipettes. Keep all reagents on ice at all times.

- 1 Pre-chill Dounce homogenizer with pestle B, containing 2 ml of ice-cold Nuclei EZ lysis buffer. 10m
- 2 Remove tissue from storage (-20°C in RNAlater-ICE) and hold in a petri dish on wet ice. If not utilizing all of the frozen^{1m} tissue, sever off the amount of tissue required (using a cold scalpel) and return the remainder to original storage site.
- 3 Using a cold scalpel cut tissue into very small pieces and transfer to Dounce homogeniser. Do not overload homogeniser with tissue (~30 mg per Dounce at one time). 5m
- 4 Homogenize for 30 strokes with pestle B (pre-chilled in ice-cold in Nuclei EZ lysis buffer. Tissue should be reduced to^{5m} ~90% homogenous small particles. [Perform the homogenization on ice, occasionally removing from the ice to check the homogenate].
- 5 With a 1 mL pipet, transfer homogenate to 15 mL conical centrifuge tube (on ice) and rinse the Dounce homogenizer^{6m} out with a further 1 mL of ice-cold Nuclei EZ lysis buffer. Take care not to transfer any large particulates of unhomogenized tissue. This is a potential source of activated RNAses that have undergone renaturation. Incubate on ice for 5 minutes.
- 6 Filter suspension through a 40 μM cell strainer in to a 50 ml conical centrifuge tube. 1m
- 7 Centrifuge nuclei/tissue suspension at 700 x g at 4°C for 5 minutes. 5m
- 8 Re-suspend isolated nuclei in 4 ml of ice-cold Nuclei EZ lysis buffer and incubate on ice for a further 5 minutes. 5m
- 9 Centrifuge nuclei suspension at 500 x g at 4°C for 5 minutes. 5m
- 10 Re-suspend nuclei in 3 mL NSB and filter through a 30 μM strainer in to a 50 mL conical centrifuge tube. 1m
- 11 Centrifuge nuclei suspension at 500 x g at 4°C for 5 minutes. 5m

- 12 Re-suspend nuclei in 1-2 mL NSB and hold on wet ice (4°C). 1m
- 13 Dilute 10 µL of nuclei suspension 1 in 10 with NSB in to a micro-centrifuge tube and count nuclei using 1:2 dilution with ^{5m} dye (either 10 µl Trypan blue or 10 µl DAPI stain (300 nM in PBS; incubate for 2 minutes with nuclei before counting) added to 10 µl of nuclei suspension) and a haemocytometer (Improved Neubauer).