

Jul 12, 2024

Cell counting

DOI

dx.doi.org/10.17504/protocols.io.14egn66rqI5d/v1

daniel.dautan daniel^{1,2}, Per Svenningsson^{1,2}

¹Department of Clinical Neuroscience, Karolinska Institutet, 171 76 Stockholm, Sweden;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA

ASAP Collaborative Rese...

Kaplitt Protocols



Eileen Ruth Torres

Weill Cornell Medicine

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.14egn66rqI5d/v1

Protocol Citation: daniel.dautan daniel, Per Svenningsson 2024. Cell counting. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.14egn66rqI5d/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 04, 2024

Last Modified: July 12, 2024

Protocol Integer ID: 101203

Keywords: ASAPCRN, neurons, apoe, chat

Funders Acknowledgement:

Aligning Science Across

Parkinson's

Grant ID: 020608



Abstract

Used for counting cells labeled with immunofluorescent markers in mouse brain sections. Sections should be stained, mounted, and imaged with high resolution (2048 x 2048 scanning).



- 1 Using 2-3 sections to encompass the brain region of interest, acquire high resolution images (2048 x 2048) using tile scanning and z-stack acquisition.
- 2 Import images into ImageJ.
- 3 Apply maximum average projection and color adjust as needed, making sure to apply the same settings to all images.
- 4 For each image, use the multi-section tool in ImageJ to manually quantify the positively-labeled cells of interest. Multiple channels can be viewed to assess cells labeled with multiple markers.
- 5 Measure the region's area using the selection tool. Use this value to derive the relative density of cells.