



Feb 20, 2020

# Genome editing in the choanoflagellate *Salpingoeca rosetta*

David Booth<sup>1</sup>

<sup>1</sup>University of California, San Francisco

**1** Works for me dx.doi.org/10.17504/protocols.io.89fhz3n

King Lab Protist Research to Optimize Tools in Genetics (PROT-G)

David Booth ⚡ 🌞 🌱

## ABSTRACT

This protocol details the preparation and execution of CRISPR/Cas9 genome editing in *S. rosetta*. The protocol builds on a [method to transfect macromolecules into \*S. rosetta\*](#) for delivering a purified Cas9 ribonucleoprotein from *Streptomyces pyogenes* (*SpCas9* RNP) into *S. rosetta*. Upon cleaving the *S. rosetta* genome at locations specified by the guide RNA (gRNA) of the *SpCas9* RNP, *S. rosetta* can use DNA oligonucleotides as templates to repair the double-stranded break. Those repair templates can encode foreign sequences and mutations for editing the *S. rosetta* genome, so long as DNA oligonucleotides have >30 bases of sequence that is homologous to both sides of the Cas9 cleavage site.

## GUIDELINES

Perform cell culturing and transfection procedure inside of a biosafety cabinet to maintain sterility.

## STEPS MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
crRNA	Custom Order	Integrated DNA Technologies
Duplex Buffer	11-01-03-01	Integrated DNA Technologies
tracrRNA	1072534	Integrated DNA Technologies
DNA Oligonucleotide	Custom Order	Integrated DNA Technologies
EnGen Cas9 NLS, <i>S. pyogenes</i> - 400 pmol	M0646T	New England Biolabs
<i>Salpingoeca rosetta</i> cultured with <i>Echinicola pacifica</i> (SrEpac)	PRA-390	ATCC
Falcon 525cm <sup>2</sup> Rectangular Straight Neck Cell Culture Multi-Flask, 3-layer with Vented Cap	353143	Corning
SF Cell Line 96-well Nucleofector Kit	V4SC-2096	Lonza
Papain from papaya latex	P3125-100MG	Millipore Sigma



## BEFORE STARTING

Please consult the attached file of media recipes for artificial seawater, high nutrient media, and low nutrient media.

☐ [MediaRecipes.pdf](#)

## Culture Cells

1 Seed a large culture of *S. rosetta*.

- 1.1 Two days prior to transfection, inoculate  **120 ml** of high nutrient media with a culture of *S. rosetta* feeding on *E. pacifica* to a final concentration of *S. rosetta* of  **8000 cells/ml** .



Salpingoeca rosetta cultured with  
Echinicola pacifica (SrEpac)

by ATCC

Catalog #: PRA-390

- 1.2 Grow the culture for  **48:00:00** in a 3-layer flask at  **22 °C with 60% humidity** .




Falcon 525cm<sup>2</sup> Rectangular Straight  
Neck Cell Culture Multi-Flask, 3-layer  
with Vented Cap

by Corning

Catalog #: 353143

## Prepare Gene Editing Cargo

- 2 **Prepare a guide RNA (gRNA) that binds to *SpCas9* and targets DNA by annealing CRISPR RNA (crRNA) with the trans-activating CRISPR RNA (tracrRNA) .**

- 2.1 Resuspend crRNA in duplex buffer (30 mM HEPES-KOH, pH 7.5; 100 mM potassium acetate) to a final concentration of  **200 Micromolar (μM)** .



crRNA

by Integrated DNA Technologies

Catalog #: Custom Order



Duplex Buffer

by Integrated DNA Technologies

Catalog #: 11-01-03-01

2.2 Resuspend tracrRNA in duplex buffer to a final concentration of **[M]200 Micromolar (μM)** .



**tracrRNA**

by Integrated DNA Technologies

Catalog #: 1072534

2.3

Mix equal volumes of crRNA ( [go to step #2.1](#) ) and tracrRNA ( [go to step #2.2](#) ) to have a final concentration of **[M]100 Micromolar (μM)** gRNA, which is the annealed complex of crRNA and tracrRNA.

2.4 Incubate the gRNA solution at **95 °C** in an aluminum block for **00:05:00** .

2.5 Place the aluminum block was placed at **Room temperature** to slowly cool the gRNA to **25 °C** .

2.6 Store the gRNA at **-20 °C** .

### 3 Prepare DNA oligonucleotides that serve as repair templates after *SpCas9* cleavage.

3.1 Dissolve oligonucleotides to a final concentration of **[M]250 Micromolar (μM)** in 10 mM HEPES-KOH, pH 7.5.



**DNA Oligonucleotide**

by Integrated DNA Technologies

Catalog #: Custom Order

3.2 Incubate the dissolved oligonucleotides at **55 °C** for **01:00:00** .

3.3 Store oligonucleotides at **-20 °C** .

3.4 Before starting nucleofections, ensure that the oligonucleotides are fully dissolved by incubating them at **55 °C** for **01:00:00** , which concurs with the assembly of the *SpCas9*/gRNA complex.



### 4 Assemble *SpCas9* with the gRNA to form the *SpCas9* RNP.



- 4.1 For one transfection, place  2 µl of **[M]20 Micromolar (µM) SpCas9** in the bottom of a 0.2 ml PCR tube.



EnGen Cas9 NLS, *S. pyogenes* - 400  
pmol  
by New England Biolabs  
Catalog #: M0646T

- 4.2 

Add  2 µl of **[M]100 Micromolar (µM) gRNA** (  go to step #2 ) by slowly pipetting up and down with *SpCas9* to gently mix the gRNA together. This solution is called the "*SpCas9* ribonucleoprotein (RNP)."

- 4.3 Incubate the *SpCas9* RNP at  **Room temperature** for  **01:00:00** (roughly the time to complete the preparation of *S. rosetta* for priming, see below).

## Prepare transfection reagents

### 5 Prepare SF Buffer (Lonza) for transfections.



SF Cell Line 96-well Nucleofector Kit  
by Lonza  
Catalog #: V4SC-2096

- 5.1 

Add all of buffer B (smaller volume that may also be called supplement 1) to buffer A (larger volume).

- 5.2 Store on ice until ready for use. The combined buffer can also be stored at 4°C for up to 3 months.



The combined buffer can be stored at 4°C for up to 3 months.



Because the Lonza kits can be so expensive, we recommend aliquoting large volumes of the SF components (900 µl aliquots for buffer A and 200 µl aliquots for buffer B) to prevent SF buffer from spoiling after buffers A and B have been combined.

### 6 Prepare the priming buffer.

## 6.1

Dilute papain to a final concentration of 100  $\mu$ M in dilution buffer (50 mM HEPES-KOH pH 7.5, 200 mM sodium chloride, 20% [v/v] glycerol, and 10 mM cysteine) from a stock solution of 1 mM papain (Millipore Sigma, St. Louis, MO; Cat. No. P3125-100MG), and incubate at room temperature just before priming cells for transfection.



### Papain from papaya latex

by Millipore Sigma

Catalog #: P3125-100MG



The dilution buffer [50 mM HEPES-KOH pH 7.5, 200 mM sodium chloride, 20% (v/v) glycerol and 10 mM cysteine] should be sterile filtered through a 0.22  $\mu$ m filter.

The dilution buffer may also be prepared ahead of time and stored in a -80°C freezer just before its use.

## 6.2 Make a solution of the remaining components of the priming buffer (40 mM HEPES-KOH, pH 7.5; 34 mM lithium citrate; 50 mM L-cysteine; 15% [wt/vol] PEG 8000). DO NOT combine the papain and priming buffer until just before adding the priming buffer to cells.






The priming buffer without papain should be sterile filtered through a 0.22  $\mu$ m filter.

The priming buffer without papain can also be made ahead of time and stored at -80°C until it is used. Be sure that the priming buffer is warmed to room temperature prior to use.

## Wash Cells

### 7 Prepare *S. rosetta* for transfection by washing away feeder bacteria.

7.1 The  120 ml culture of *S. rosetta* feeding on *E. pacifica* (  go to step #1 ) was homogenized by vigorous shaking and then split into  40 ml aliquots in 50 ml conical tubes.

### 7.2

Vigorously shake the aliquots and centrifuge the cells for  00:05:00 at  2000 x g and  22 °C in a swinging bucket rotor.

7.3 Use a serological pipette to gently remove from the cell pellet all but 2 ml of the supernatant, which remains cloudy with *E. pacifica* bacteria. With a fine tip transfer pipette, gently remove the remaining liquid near the pellet.

7.4 The three cell pellets were resuspended in a total volume of **50 ml** artificial seawater, combined into one conical tube, and vigorously shaken to homogenize the cells.

7.5 

For a second time, the resuspended cells were centrifuged for **00:05:00** at **2000 x g** and **22 °C** in a swinging bucket rotor.

7.6 The supernatant was removed as before ( **go to step #7.3** ).

7.7 The pellet was resuspended in **50 ml** of artificial seawater, and the cells were homogenized by vigorous shaking.

7.8 

The cells were centrifuged for a third time for **00:05:00** at **2200 x g** and **22 °C**.

7.9 Remove the supernatant as before ( **go to step #7.3** ).

7.10 Resuspend the cell pellet in **400 µl** of artificial seawater. This resuspension is called the "washed cells."

8 **Prepare 100 µl aliquots of 50000000 cells/ml .**

8.1 Dilute **2 µl** of "washed cells" ( **go to step #7.10** ) into **196 µl** of artificial seawater.

8.2 Fix the diluted cells with **2 µl** of 37.5% formaldehyde and homogenize by vortexing.

8.3 Pipet the fixed cells into a fixed chamber slide and determine the cell concentration.



Remember that concentration of diluted and fixed cells is a 100-fold dilution from the "washed cells." Be sure to factor that dilution into your concentration.



Cells can be counted on a hemacytometer (Neubauer with brightlines) or with an automated cell counter. We recommend a Luna-FL automated cell counter.



LUNA-FL  
Dual Fluorescence Cell Counter  
Logos Biosystems L20001 [↗](#)

8.4 After determining the cell concentration, dilute the "washed cells" to final concentration of **5000000 cells/ml** and split into **100 µl** aliquots.



One aliquot provides enough cells for 12 nucleofections.

## Prime Cells

9 **Prime cells for nucleofection by degrading the glycocalyx that surrounds *S. rosetta*.**

9.1 **Spin the 100 µl aliquots of washed cells ( [go to step #8.4](#) ) at 800 x g and 22 °C for 00:05:00 .**

9.2 Gently remove the supernatant from the cell pellet with a gel-loading pipette tip.







9.3 

Combine the priming buffer components ( [go to step #6](#) ) to make a final priming buffer (40 mM HEPES-KOH, pH 7.5; 34 mM lithium citrate; 50 mM L-cysteine; 15% [wt/vol] PEG 8000; and 1 µM papain)

9.4 Resuspend each cell pellet in **100 µl** of priming buffer.





9.5 

Incubate cells for **00:35:00** at **Room temperature** .


- 9.6 Add  **10 µl** of **50 mg/ml Bovine Serum Albumin** to each aliquot of primed cells for quenching proteolysis from the priming buffer.
- 9.7 Centrifuge cells at  **1250 x g** and  **22 °C** for  **00:05:00**.
- 9.8 Gently remove the supernatant from the cell pellet with a gel-loading pipette tip.
- 9.9 Resuspended each cell pellet in  **25 µl** of SF Buffer (  **go to step #5** ). This suspension of cells is called the "primed cells."
- 9.10 Store the "primed cells" on ice while preparing nucleofection reactions.







## Transfect Cells

### 10 Deliver gene editing cargo via nucleofection.


- 10.1 Add  **16 µl** of ice-cold SF Buffer (  **go to step #5** ) to the *SpCas9* RNP (  **go to step #4.3** ), which has a total volume of  **4 µl**.



For reactions that use two different gRNAs, assemble each *SpCas9* RNP separately then combine each *SpCas9* RNP at this step. After the *SpCas9* RNPs have been combined, add  **16 µl** of ice-cold SF Buffer

- 10.2 Add  **2 µl** of the repair oligonucleotide template to the PCR tube with *SpCas9* RNP and SF Buffer (  **go to step #10.1** ).
- 10.3 Add  **2 µl** of "primed cells" (from  **go to step #9.10** ) to the PCR tube with *SpCas9* RNP, SF Buffer, and the repair template (  **go to step #10.2** ). This solution, which is called the "nucleofection mix," should have a total volume of  **24 µl**.
- 10.4 Transfer the entire nucleofection mix into one well of a 96-well nucleofection plate.



At this point, prepare for the recovery step, by transferring the recovery buffer into a convenient vessel and setting the pipette to  **100 µl**.



10.5 Pulse the nucleofection plate with the CM156 pulse.



**4D-Nucleofector Core Unit**

Control system for performing nucleofection

Lonza AAF-1002B [↗](#)



**96-well Shuttle Device**

Add-on for Nucleofector 4d device to perform plate-based nucleofections

Lonza AAM-1001S [↗](#)

## Rest and Recover Cells

11 **Allow membranes to reseal by resting cells in recovery buffer before growing cells again in media.**

11.1 Immediately after transfection, add **100 µl** of ice-cold recovery buffer (10 mM HEPES-KOH, pH 7.5; 0.9 M sorbitol; 8% [wt/vol] PEG 8000) to each nucleofection transfection and gently mixed by firmly tapping the side of the plate.

11.2 Allow cells to rest in recovery buffer for **00:05:00**.

11.3 Gently mix the well in the nucleofection plate by pipetting up and down before transferring the entire volume in nucleofection well (the nucleofection mix plus the recovery buffer) into to **2 ml** of low nutrient media in one well of a 6 well plate.

11.4 Incubate at **22 °C and 60% humidity** for **00:30:00**

12 **Add *E. pacifica* food and grow transfected cells.**

12.1 Add **10 µl** of **10 mg/ml** of *E. pacifica* to the wells in the 6 well plate.



Prepare the *E. pacifica* solution by resuspending a frozen, 10 mg pellet of *E. pacifica* in 1 ml of artificial seawater.

12.2



Incubate the 6 well plate at **22 °C and 60% humidity** for **24:00:00** before using in downstream experiments.

#### (Optional) Select for Cycloheximide Resistance

13



Add **10 µl** of **1 mg/ml** of cycloheximide to the **2 ml** culture of transfected cells after allowing the cells to fully recover.



Cycloheximide is toxic. Handle carefully and properly dispose.

14



Incubate the cells in cycloheximide for **96:00:00** prior to genotyping and clonal isolation.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited