



DEC 22, 2023

DNA Barcoding Protocol for Arthropods

Valerie Warhol¹

¹Science Club



Valerie Warhol

ABSTRACT

This is a basic protocol for doing extraction and amplification of DNA barcodes (COI) for arthropods (typically small insects or spiders).

PROTOCOL MATERIALS

Wash Buffer Carolina Biological Supply Catalog #C33428	Step 16
Silica Resin Carolina Biological Supply Catalog #C33426	Step 10
Guanidine Hydrochloride 6M Carolina Biological Supply Catalog #C33427	
	Step 5
EZ PCR Master Mix 5X miniPCR	Step 36

OPEN ACCESS



Protocol Citation: Valerie Warhol 2023. DNA Barcoding Protocol for Arthropods. **protocols.io** <https://protocols.io/view/dna-barcoding-protocol-for-arthropods-c6i7zchn>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Dec 19, 2023

Last Modified: Dec 22, 2023

PROTOCOL integer ID:
92479


Prepare sample and equipment

- 1 Make sure all instruments, such as forceps and pestle, are clean and sterile.
- 2 Prepare water bath at 65°C .
- 3 Dissect sample from specimen (typically 1 leg). Return specimen to freezer.
- 4 Let sample dry for 10 minutes to remove any ethanol.
- 5 Prepare a clean 1.5 mL hinged tube by writing sample ID on it and filling with $250\ \mu\text{L}$ of
Guanidine Hydrochloride 6M Carolina Biological
Supply Catalog #C33427


Lyse cells


11m

- 6 Put sample in tube. Grind sample with pestle until broken up into tiny pieces.
- 7 Incubate sample tube in 65°C water bath for 00:10:00. 10m
- 8 Remove tube and lower temperature of water bath to 57°C .

9 Centrifuge tube for  00:01:00 at maximum speed to pellet debris.

1m

10 Remove  Silica Resin Carolina Biological
Supply Catalog #C33426 from refrigerator.

11 Label a clean 1.5 μ L tube with sample number. Transfer  150 μ L of the supernatant to the clean tube. Discard old tube containing debris.

Bind DNA

5m 30s

12 Add  3 μ L of silica resin to tube. Mix well by pipetting up and down several times.

13 Close tube and incubate for  00:05:00 in  57 °C water bath.

5m


14 Centrifuge for  00:00:30 at maximum speed to pellet the resin.



30s

15 Use a pipette with a fresh tip to remove the supernatant, being careful not to disrupt the pellet.

Wash



1m



16 Remove molecular grade water from refrigerator and
 Wash Buffer Carolina Biological
Supply Catalog #C33428 from freezer.

- 17 Add  500 µL of ice-cold wash buffer to the pellet. Mix well by pipetting up and down several times to resuspend the silica resin.
- 18 Close the tube and centrifuge for  00:00:30 at maximum speed to pellet the resin. 30s
- 19 Use a pipette with a fresh tip to remove the supernatant, being careful not to disrupt the pellet.
- 20 Again, add  500 µL of ice-cold wash buffer to the pellet. Mix well by pipetting up and down to resuspend the silica resin.
- 21 Close the tube and centrifuge for  00:00:30 at maximum speed to pellet the resin. 30s
- 22 Return wash buffer to freezer.
- 23 Use a pipette with a fresh tip to remove the supernatant, being careful not to disrupt the pellet. Spin the tube briefly to collect any remaining drops of supernatant, and then remove these with a pipette.

Elute DNA


30s

- 24 Add  100 µL of molecule grade water to the silica resin and mix by pipetting up and down several times.
- 25 Incubate the mixture at  57 °C for 5 minutes.


- 26 Centrifuge for  00:00:30 at maximum speed to pellet the resin. 30s
- 27 Label a clean 1.5 µL tube with sample number. Transfer  90 µL of the supernatant to the clean tube, being careful not to disturb the pellet. Discard old tube containing the resin.
- 28 Store sample in freezer until ready to PCR. If going directly to PCR, put sample in refrigerator.







Amplify COI (PCR)

10m

- 29 For each DNA sample, label a PCR microtube with sample ID.
- 30 Turn on PCR thermal cycler and connect to computer.
- 31 Remove molecular grade water from refrigerator. Remove template DNA, PCR master mix, and primer 10m from freezer. Let thaw for  00:10:00 .

A	B	C	D
Primer Name	Direction	Sequence	Concentration
LC01490	Forward	GGTCAACAAATCATAAAGATATTGG	10 µM
HCO2198	Reverse	TAAACTTCAGGGTGACCAAAAAATCA	10 µM

- 32 Add  31 µL molecular grade water to each microtube.

- 33 Mix forward primer by flicking. Add  2 µL forward primer to each microtube.
- 34 Mix reverse primer by flicking. Add  2 µL reverse primer to each microtube.
- 35 Mix template DNA by flicking. Add  5 µL template DNA to each microtube.
- 36 Mix  EZ PCR Master Mix
5X miniPCR by inverting. Add  10 µL PCR master mix to each microtube.
- 37 Put microtubes in PCR thermal cycler and start PCR.
[To do: Add PCR program steps.]
- 38 Remove microtubes from PCR thermal cycler. Verify PCR using gel electrophoresis.
- 39 Store PCR products in freezer ( -20 °C).