



Sep 09, 2020

# RNase One Ribonuclease digestion

Roey Angel<sup>1</sup>, Eva Petrova<sup>1</sup><sup>1</sup>Soil and Water Research Infrastructure

In Development

This protocol is published without a DOI.

SoWa RI Anaerobic and Molecular Microbiology (public)

Tech. support email: [eva.petrova@bc.cas.cz](mailto:eva.petrova@bc.cas.cz)

Eva Petrova

Soil and Water Research Infrastructure

## ABSTRACT

Modified from: <http://www.promega.com/paguide/chap6.htm#RNaseAmount>

This protocol was used to remove RNA from DNA preparations.

The RNases commonly used in an RPA are RNase I, RNase T1 or a combination of RNase T1 and RNase A. RNase I cleaves after all four ribonucleotides, RNase T1 cleaves after G residues and RNase A cleaves after A and U residues. For most efficient cleavage of the single-stranded regions immediately adjacent to the double-stranded RNA, we recommend RNase I (RNase ONE™ Ribonuclease, Cat.# M4261).

Use 1–10u RNase ONE™ Ribonuclease per 10µg of total RNA or 0.1–1 u of RNase ONE™ Ribonuclease per 1µg of poly(A)+ RNA.

## PROTOCOL CITATION

Roey Angel, Eva Petrova 2020. RNase One Ribonuclease digestion. **protocols.io**  
<https://protocols.io/view/rnase-one-ribonuclease-digestion-qdzds76>

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

May 24, 2018

## LAST MODIFIED

Sep 09, 2020

## PROTOCOL INTEGER ID

12441

## GUIDELINES

### Annealing and Digestion Temperatures

RNase ONE™ Ribonuclease digestion works most efficiently at 20–37°C (Brewer *et al.* 1992).

## MATERIALS

NAME	CATALOG #	VENDOR
<a href="#">Roti®-Phenol/Chloroform/Isoamyl alcohol ()</a>	A156.3	<a href="#">Carl Roth</a>
<a href="#">RNase ONE™ Ribonuclease</a>	M4261	<a href="#">Promega</a>
<a href="#">Proteinase K Solution (20 mg/mL)</a>	AM2548	<a href="#">Thermo Fisher Scientific</a>
<a href="#">Gel Loading Dye Purple (6X) - 4.0 ml</a>	B7024S	<a href="#">New England Biolabs</a>
<a href="#">Yeast tRNA (10 mg/mL)</a>	AM7119	<a href="#">Thermo Fisher Scientific</a>

## STEPS MATERIALS

NAME	CATALOG #	VENDOR
RNase ONE™ Ribonuclease	M4261	Promega
Proteinase K Solution (20 mg/mL)	AM2548	Thermo Fisher Scientific
Roti®-Phenol/Chloroform/Isoamyl alcohol ()	A156.3	Carl Roth
Yeast tRNA (10 mg/mL)	AM7119	Thermo Fisher Scientific

BEFORE STARTING

### Materials Required:

- purified total RNA or poly(A)+ RNA from the tissue or cells of interest
- 3.0M ammonium acetate (pH 5.2)
- ice-cold 100% ethanol
- ice-cold 70% ethanol
- 20% w/v SDS
- RPA loading dye

- 1 Incubate 5–10µg of total NA 85°C for 5 minutes to denature the NAs.


 **5 µg total NA**

 **85 °C**

 **00:05:00**

- 2 Add 300µl of RNase digestion buffer and the appropriate amount of RNase ONE™ Ribonuclease (1-10u per 10µg of total RNA). Incubate the samples for 30–60 minutes at 20–37°C.

 **300 µl RNase digestion buffer**



**RNase ONE™ Ribonuclease**  
by Promega  
Catalog #: [M4261](#)

 **00:30:00**

 **20 °C**


- 3 Stop the reaction as follows: Add 10µl of 20% w/v SDS and 2.5µl of 20mg/ml proteinase K. Incubate for 15 minutes at 37°C.

 **10 µl 20% w/v SDS**

 **37 °C**

 **00:15:00**

 **2.5 µl 20mg/ml proteinase K**



**Proteinase K Solution (20 mg/mL)**  
by Thermo Fisher Scientific  
Catalog #: [AM2548](#)

- 4 Extract once with phenol:chloroform:isoamyl alcohol, and remove the aqueous phase to a clean microcentrifuge tube containing 1µl of 10mg/ml carrier tRNA .



Roti®-Phenol/Chloroform/Isoamyl  
alcohol ()  
by Carl Roth  
Catalog #: A156.3

 **1 µl 10mg/ml carrier tRNA**



A yeast tRNA is an effective coprecipitant to aid in recovery of small amounts of nucleic acids.




Yeast tRNA (10 mg/mL)  
by Thermo Fisher Scientific  
Catalog #: AM7119

- 5 Add 825µl of ice-cold 100% ethanol, and chill at –20°C for 30 minutes, then centrifuge at maximum speed in a microcentrifuge for 15 minutes at 4°C to pellet the RNA.

 **825 µl ice-cold 100% ethanol**


 **-20 °C**

 **00:30:00**



#### **CENTRIFUGE**

max. speed

 **00:15:00**

 **4 °C**

- 6 Carefully remove the supernatant. Wash the pellet using 300µl of 70% ethanol and dry.

 **300 µl 70% ethanol**

- 7 Resuspend the pellet.