



MAR 14, 2024

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.5qpvo3jmvz4o/v1

Protocol Citation: Xiaobo Mao, Ramhari Kumbhar, Hanseok Ko, Valina L. Dawson, Ted Dawson 2024. Cylinder Test. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5qpvo3jmvz4o/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Mar 14, 2024

Cylinder Test

Xiaobo Mao^{1,2,3,4}, Ramhari Kumbhar^{1,2,3,4}, Hanseok Ko^{1,2,4}, Valina L. Dawson^{1,2,4,3,5,6}, Ted Dawson^{1,2,4,3,5,6,7}

¹Neuroregeneration and Stem Cell Programs, Institute for Cell Engineering, Johns Hopkins University School of Medicine, Baltimore, MD 21205, USA;

²Department of Neurology, Johns Hopkins University School of Medicine, Baltimore, MD 21205, USA;

³Adrienne Helis Malvin Medical Research Foundation, New Orleans, LA 70130-2685, USA;

⁴Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA;

⁵Department of Physiology, Johns Hopkins University School of Medicine, Baltimore, MD 21205 USA;

⁶Solomon H. Snyder Department of Neuroscience, Johns Hopkins University School of Medicine, Baltimore, MD 21205, USA;

⁷Institute for NanoBioTechnology, Johns Hopkins University, Baltimore, MD, USA

ASAP Collaborative Research Network

Kaplitt Protocols



Eileen Ruth Torres
Weill Cornell Medicine

ABSTRACT

The Cylinder test is used to evaluate locomotor asymmetry and spontaneous movement in rodent models.

MATERIALS

- Small transparent cylinder (height, 15.5 cm; diameter, 12.7 cm)
- Video recorder

Last Modified: Mar 14, 2024

PROTOCOL integer ID: 96716

Keywords: ASAPCRN, behavior,
mouse model

Funders Acknowledgement:

Aligning Science for Parkinson's
Grant ID: 020608

- 1 Set up small, transparent cylinder in front of video recorder.
- 2 Place mouse in cylinder and record spontaneous activity for 10 min.
- 3 Clean cylinder in between mice.
- 4 For all recordings, view in slow motion and count the number of forepaw touches to the cylinder walls, rears, and grooming bouts.