



OCT 19, 2023

Nuclei Isolation for HMBA FACS

Lakme
Caceres¹

¹Princeton Neuroscience Institute



Lakme Caceres
Princeton University

ABSTRACT

This protocol is for purifying nuclei for downstream 10X sequencing.

GUIDELINES

Keep tissue/nuclei on ice as much as possible.

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.kxygx35ywg8j/v1

Protocol Citation: Lakme Caceres 2023. Nuclei Isolation for HMBA FACS. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kxygx35ywg8j/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Oct 17, 2023

Last Modified: Oct 19, 2023

Keywords: nuclei, isolation,
FACS

Prepare Stock Solutions

- 1 Make 50 mL **Nuclear Isolation Media** by filling a 100 mL bottle with 25 mL of MilliQ water and then adding 500 μ L 1M Tris, 1.25 mL 1M KCl, 250 μ L 1M MgCl₂, and 4.28 g Sucrose. Shake until sucrose is dissolved then fill to 50 mL with MilliQ water. Filter it using a stericup into another clean bottle. (4°C - 2 weeks)
- 2 Make 20 mL **10% Triton X-100** by combining 18 mL MilliQ water with 2 mL Triton X-100 in a 50 mL tube. Vortex and then incubate at room temperature for 20 minutes. Filter it through a 0.22 μ m filter with a syringe into a clean 50 mL tube. (4°C - 1 month)
- 3 Make 500 μ L **10 ug/mL DAPI** by adding 5 μ L of 10 mg/mL DAPI stock to 495 μ L of 1X PBS.

Prepare Fresh Solutions

- 4 Make 3 mL **Homogenization Buffer** by adding 2.892 mL Nuclear Isolation Media to a 5 mL eppendorf. Then add 60 μ L protease inhibitor, 30 μ L 10% Triton X-100, 15 μ L RNase inhibitor, and 3 μ L 100 mM DTT.
- 5 Make 15 mL **Blocking Buffer** by adding 13.425 mL 1X PBS, 1.5 mL 10% BSA, and 75 μ L RNase inhibitor to a 15 mL falcon tube.

Homogenization

- 6 Clean dounce, scalpel, and forceps using MilliQ water, ethanol, RNase Zap, then MilliQ again. The red-tape forceps are for unfixed tissue.
- 7 Place your sample on dry ice. Tare a weigh boat or tube on the scale and weigh your sample.

- 8 Add tissue to dounce and push it to the bottom using 1 mL of Homogenization Buffer and the pestle. Homogenize the tissue without creating bubbles. Then add the remaining 2 mL of the Homogenization Buffer and continue to dounce until homogenized.
- 9 Pass all of the nuclei suspension through a 40 μ m filter.
- 10 Centrifuge at 900 g/rcf for 10 minutes at 4°C.

Blocking Buffer & Antibodies

- 11 Add 950 μ L Blocking Buffer to each of your control tubes.
- 12 The amount of Blocking Buffer we will need to resuspend our nuclei pellet with after centrifuging depends on the number of controls we will run. Each control tube will require 50 μ L of nuclei suspension, and our sample needs to have 1,000 μ L of nuclei suspension left over. In this case we are running two controls.
- 13 After the cell homogenate is done centrifuging, discard the supernatant and resuspend the pellet in 1100 μ L of Blocking Buffer. Incubate for 15 minutes on ice.
- 14 Add 50 μ L of nuclei suspension to each control tubes.
- 15 Add the corresponding antibodies to the control tubes:

- 0.5 μ L NeuN-PE (1/2000) to the PE-only control
- 0.33 μ L Olig2-Alexa 488 (1/3000) to the A488-only control

16 Add both 0.5 μ L NeuN-PE and 0.33 μ L Olig2-Alexa 488 to the sample tube.

17 Place all tubes in a rotator placed in a dark 4°C fridge and incubate for 30 minutes.

Wash

18 After incubation, spin down all of the tubes at 500 g for 5 minutes at 4°C.

19 *If there is no visible pellet* after centrifugation, remove 980 μ L of supernatant with a small pipette being careful not to disturb the bottom of the tube. Then resuspend the pellet in 980 μ L Blocking Buffer.

If there is a visible pellet, remove all of the supernatant without disturbing the pellet and then resuspend in 980 μ L of Blocking Buffer.

20 Add 20 μ L 10 μ g/mL DAPI to each sample and control tube and invert to mix.

21 Add 300 μ L Blocking Buffer to falcon FACS tubes and label them A488+, PE+, and Double Negative. We will sort into these.

22 Store tubes on ice until FACS. Bring DAPI stock and Blocking Buffer with you to FACS in case we have to dilute the sample to lower the event count.

