



Feb 14, 2020

Detection of viable *Dichelobacter nodosus* by real-time PCR using PMAxx™

PLOS One

Peter Kuhnert¹, Tobias Hidber¹

¹Institute of Veterinary Bacteriology, Vetsuisse Faculty, University of Bern, Switzerland

1 Works for me dx.doi.org/10.17504/protocols.io.bb9ij96



Peter Kuhnert
Institute of Veterinary Bacteriology, Vetsuisse Faculty, Uni...



ABSTRACT

Dichelobacter nodosus is a gram-negative fastidious anaerobic bacterium and the causative agent of ovine footrot. The disease has a global presence and is endemic in many countries. Clinical symptoms range from mild interdigital dermatitis in benign footrot to severe underrunning and separation of the hoof horn from the underlying tissue in virulent footrot. Clinical symptoms start as early as 2 weeks after first contact and the disease leads to pain, lameness, decreased meat and wool production as well as animal welfare issues. Diagnosis has improved significantly by the development of a real-time PCR to detect and discriminate virulent (*aprV2*-positive) and benign (*aprB2*-positive) *D. nodosus* strains. However, the real-time PCR also detects dead cells making its use for testing e.g. disinfecting agents limited. A PMA (propidium monoazid) real-time PCR using the improved dye PMAxx™ was therefore developed for virulent *D. nodosus* that allows discrimination of viable and dead bacteria. The distinction between viable and non-viable cells is possible, based on membrane integrity. For that purpose, the samples containing *D. nodosus* are treated with the improved nucleic acid intercalating PMA dye PMAxx™ that selectively enters cells with compromised cell membranes, whereas the intact cell membrane presents a natural barrier for this molecule. After exposure to strong light, it covalently binds to the DNA, preventing DNA from being amplified by PCR, thereby enabling differentiation of viable from non-viable cells. The PMA-qPCR proved to be a valid method for comparison of antimicrobial efficiency in *ex vivo* experiments.

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0229066>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

In vitro and ex vivo testing of alternative disinfectants to currently used more harmful substances in footbaths against *Dichelobacter nodosus*

MATERIALS

NAME	CATALOG #	VENDOR
TaqMan™ MGB Probe	4316034	Thermo Fisher
7500 Real-Time PCR System, desktop	4351105	Thermo Fisher
TaqMan™ Fast Advanced Master Mix	4444556	Thermo Fisher
KingFisher™ Duo Prime Purification System	5400110	Thermo Fisher
VetMAX™ Xeno™ Internal Positive Control DNA	A29764	Thermo Fisher
VetMAX™ Xeno™ Internal Positive Control - LIZ™ Assay	A29766	Thermo Fisher
PMAxx™ Dye 20 mM in H2O	#40069	Biotium
LED Floodlight 100 W 6000 K 8500 lm	FL5836	Optonica
MagneSil®RED	A1641	Promega

PMAxx treatment

- 1 place or dissolve your sample containing *Dichelobacter nodosus* in a total volume of 50 µl 0.85% NaCl in a transparent 1.5 ml tube
- 2 in a darkened room add 100 µM (0.25 µl) of PMAxx™ to the tube
- 3 vortex and incubate for 3 min at room temperature in a metal box impervious to light
- 4 after incubation, place tube on ice and expose to LED light (white light, 100 W, 6000 K, 8500 lm) at a distance of 20 cm for 5 min
- 5 centrifuge tube at 15'000xg for 5 min in a microcentrifuge
- 6 to remove the supernatant, the tubes can be held in a horizontal position and twisted tissues used to absorb the liquid
- 7 resuspend pellet in 500 µl SV-lysis buffer (4 M guanidine thiocyanate, 0.01M Tris-HCl, 1% β-mercaptoethanol)

DNA extraction

- 8 add 2 µl of VetMAX™ Xeno™ Internal Positive Control DNA to the sample in SV-lysis buffer (this is the extraction control and later also the internal PCR control)
- 9 add 30 µl MagneSil®RED and resuspend without vortexing
- 10 extract DNA using KingFisher™ Duo Prime Purification System (or an alternative magnetic bead based system) and elute the purified DNA in 60 µl H₂O

real-time PCR

- 11 primers and probes as well as the real-time PCR are described in the publications [116. Dnodocus PCR.pdf](#)
[140. Dnodocus infection SAT.pdf](#)
- 12 prepare 25 µl containing:
 - 1x TaqMan™ Fast Advanced Master Mix
 - 300 nM primer DnApr™-L (CAATAGCCAAATTTCTTTAGATGGTGAT)
 - 300 nM DnApr™-R (CAAGAGCTGTCGCTTCTTTCTTT)
 - 100 nM Probe DnApr™-v (FAM-CGGTGGTTATCCTGAT-MGB)
 - 250 nM Probe DnApr™-b (VIC-TGGTCGTCCTGATC-MGB)
 - pyrogen-free water
 - 1 µl Xeno™ LIZ Primer Probe Mixand 2.5 µl of DNA template
- 13 prepare at least duplicates for each sample and run on an appropriate instrument like 7500 Real-Time PCR-System instrument
- 14 run using cycles of 2 min at 50°C and 10 min at 95°C followed by 40 cycles with 15 s at 95°C and 1 min at 60°C
- 15 results are analyzed using the Sequence Detector 7500 software (v 2.3.) with the threshold set at 0.015
- 16 signals will only result from alive *D. nodosus* and no dead cell will be detected. For absolute quantification a standard curve using defined amounts of *D. nodosus* genome equivalents can be generated in parallel



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited