



Nov 13, 2020

BSCI:414--Lab 11:Confirming Cloning of spike(RBD)-SUMO/pET28

Harley King¹

¹NIST Center for Neutron Research National Institute of Standards and Technology Gaithersburg, MD 20899 Department of Materials Science and Engineering University of Maryland College Park, MD 20742-2115

1 Works for me This protocol is published without a DOI.

USG Fall 2020 BSCI:414



Harley King

NIST Center for Neutron Research National Institute of Stand...

ABSTRACT

We've amplified our target genes (SARS-CoV-2 spike and spike (RBD) and SUMO) with primers. Then we purified and performed Gibson Assembly. After transformation we tested a few colonies and found that they were positive. We sent the purified plasmid away for sequencing and now we wish to confirm if we correctly cloned.

PROTOCOL CITATION

Harley King 2020. BSCI:414--Lab 11:Confirming Cloning of spike(RBD)-SUMO/pET28 . **protocols.io**
<https://protocols.io/view/bsci-414-lab-11-confirming-cloning-of-spike-rbd-su-bpm8mk9w>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 12, 2020

LAST MODIFIED

Nov 13, 2020

PROTOCOL INTEGER ID

44448

ABSTRACT

We've amplified our target genes (SARS-CoV-2 spike and spike (RBD) and SUMO) with primers. Then we purified and performed Gibson Assembly. After transformation we tested a few colonies and found that they were positive. We sent the purified plasmid away for sequencing and now we wish to confirm if we correctly cloned.

Watch Cloning Video

- 1 Review the Benchling notebook regarding Gibson Assembly for cloning the spike-SUMO and RBD genes into pET28.
<https://benchling.com/s/etr-kXqsFUhKxGCbyUGyCx9V>

Creating Spike-related Plasmids for Heterologous E. coli Expression · Benchling.pdf

Create Alignment with Chromatograms


- 2 Open the "[F20_spike\(RBD\)-SUMO/pET28](#)" plasmid from our BSCI:414 plasmids folder. Make a copy into a new folder under "BSCI:414 Lab 11".
- 3 Download the 4 sequence chromatograms from our class folder onto your computer.
 - BSCI414>Files>Lab Results>Lab 11 spike cloning>[RBD_pET28 chromatograms](#)

These chromatograms represent sequence using single primers: T7_promoter and T7_terminator.

Primer	Primer Sequence
T7_promoter	TAATACGACTCACTATAGGG
T7_terminator	gctagttattgctcagcgg

Sequencing primers

- 4 Select the "Alignments" icon on the right side.



Benching "Alignments" icon.

- 5 Add the chromatograms and select "Align".

Are there differences?

What kind of differences?

Which Cell Line Should Be Transformed?

- 6 Next steps: transformation but into which cell lines?