



JAN 25, 2023

OPEN  ACCESS**DOI:**[dx.doi.org/10.17504/protocols.io.3byl4j5xolo5/v1](https://dx.doi.org/10.17504/protocols.io.3byl4j5xolo5/v1)

**Protocol Citation:** Iris Dadole, Kevin Huguet, Didier Blaha, Nicolas Personnic 2023. Protocol for high-throughput isolation of bacterial intracellular nonreplicating persisters.

**protocols.io**<https://dx.doi.org/10.17504/protocols.io.3byl4j5xolo5/v1>**MANUSCRIPT CITATION:**<https://www.nature.com/articles/s41467-019-13021-8>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Jan 12, 2023**Last Modified:** Jan 25, 2023

## Protocol for high-throughput isolation of bacterial intracellular nonreplicating persisters


Iris Dadole<sup>1</sup>, Kevin Huguet<sup>1</sup>, Didier Blaha<sup>1</sup>, Nicolas Personnic<sup>1</sup>

<sup>1</sup>CIRI, Centre International de Recherche en Infectiologie, CNRS UMR 5308, INSERM U1111, Ecole Normale Supérieure de Lyon, Université Claude Bernard Lyon 1, France

Iris Dadole: Equal contribution

Kevin Huguet: Equal contribution

Nicolas Personnic: Corresponding author



Nicolas Personnic

**DISCLAIMER**

There are no competing interests relating to employment, consultancy, patents, products in development, or marketed products. We present a full and objective manuscript that has no interference to its objectivity.

**Keywords:** Persisters, Single-cell, Timer, Growth rate, Infection, Pathogens, Host, Antibiotic, FACS

## ABSTRACT

Failure of antibiotics to clear a bacterial infection is arguably one of the most severe threat to global health. Bacteria ability to withstand drug treatment partly relies on the formation of a subpopulation of persisters. Persisters are rare nonreplicating individuals, transiently refractory to antibiotic treatment and responsible for non-resolving bacterial infection and relapse. Studying persisters is challenging as it is a reversible physiological state occurring in a small fraction of the pathogen population hidden in a complex environment. Persisters have been primarily analyzed in broth, hardly reflecting the conditions prevailing within a host. In addition, indirect read-outs used to detect the persisters, *a posteriori*, limited a comprehensive understanding of the persisters peculiar physiology during the infection. Here, we describe a high-throughput method to track, and collect the bacterial persisters generated during the infection of phagocytic cells, with single-cell resolution. Persisters are detected based on the bacterial division rate measured for each individual bacterium using the Timer fluorescent growth rate reporter and subsequently recover using fluorescence-activated cell sorting (FACS). We will illustrate this method with *Legionella pneumophila*, a facultative intracellular pathogen and the agent of the Legionnaires' disease for which persistence to antibiotics and relapsing infections are clinically documented. Isolated persisters can be analyzed biochemically and functionally, yielding fundamental new insights into the molecular and cellular processes underlying persistence, and into the peculiar physiology of the persisters within the host. This method is adaptable to other microbial pathogens and alternative infection models.

## ATTACHMENTS

[ka3mb2zrp.docx](#)

## GUIDELINES

### SUMMARY:

Persisters represent a subpopulation of nonreplicating individuals transiently refractory to antibiotic treatment. Here we describe a high-throughput single-cell method to track and collect the persisters formed within phagocytes. Cells are infected with bacteria expressing the Timer fluorescent growth rate reporter, and nonreplicating individuals are recovered using fluorescence-activated cell sorting (FACS).

### INTRODUCTION:

During the colonization of the host, the pathogen encounters highly heterogeneous host micro-environments [1], elaborate a complex interplay with the component of the immune systems [2] and/or develops microbial communities with high degree of phenotypic specialization [3] which produce, within the micrometer ranges, a

multitude of genetically identical individuals with alternative physiologies [4]. Bacterial variability is further increased as clonal bacteria constantly develop functionally distinct subpopulations [5]. This reversible and regulated phenomenon termed phenotypic heterogeneity is of major clinical importance as it implies the formation of persisters [6]. Persisters are preestablished or stress induced individuals that transiently evade the bactericidal activity of antibiotics by entering a rare and peculiar physiological state at the cost of abrogation of bacterial division, and finally resume growth on the onset of the treatment termination [7]. The emergence of persisters has been documented for major bacterial pathogens including *Staphylococcus aureus*, *Mycobacterium tuberculosis*, *Escherichia coli*, *Salmonella enterica*, *Pseudomonas* spp., *Listeria monocytogenes* and *Legionella pneumophila* [4,6,8–16]. Traditionally, the presence of bacterial persisters has been revealed by the “biphasic killing” phenomenon of bacterial cultures exposed to bactericidal antibiotics [17]. However, such population-averaging methods, performed a posteriori, are poorly suited to comprehensively and dynamically analyze a phenomenon where reversible cell-to-cell variations plays a key role. The development of high-throughput fluorescence-based single-cell technologies has rendered important breakthroughs reachable. Notably, the microfluidic device called mother machine has been incredibly useful to directly monitor the formation of persisters and the underlying regulatory mechanisms using a combination of genetic and environmental manipulations [6]. Persisters have been mostly studied using non-infectious model organisms cultivated in broth, hardly reflecting the conditions prevailing within a host. Deciphering the biology of the persisters during the infection has remained technically challenging as the persisters represent a rare and unstable subpopulation of the pathogen, that is generated within complex environments.

The ubiquitous environmental Gram-negative bacterium *Legionella pneumophila* is a waterborne human pathogen responsible for a life-threatening pneumonia termed Legionnaires’ disease (LD) [18], the inventoried cases of which have doubled since 2013 [19]. LD is fatal for 5-15% of the patients (up to 50% among the elderly), even when promptly diagnosed and treated [20]. A nation-wide retrospective study describing those cases of slowly or non-resolving LD proposed that recalcitrance to antibiotic treatments or relapsing infections were the expression of bacterial persistence [21]. *L. pneumophila* survives its ingestion by a diverse array of protozoan predators, as well as lung alveolar macrophages, and establishes a replicative membrane-bound compartment termed the *Legionella*-containing vacuole (LCV) [22,23]. LCV formation requires the *Icm/Dot* type IV secretion system (T4SS) and involves approximately 300 different so-called « effector » proteins, which are translocated into host cells [24,25] to hijack numerous host cell functions. The amoebozoan *Acanthamoeba castellanii* has proved a reliable, easy-to-use and cost effective model organism to investigate the molecular interactions of *L. pneumophila* with its host cell, and at large, the pathogenic potential and process of many other human opportunistic pathogens [26].

The plasmid based Timer growth rate reporter [4] is a powerful technology to track the nonreplicating persisters produced during an infection. Timer is a stable fluorescent protein that slowly matures from a green to a red fluorescent protein. As a consequence, in growing bacteria, constitutively producing Timer, green fluorescent Timer dominates over red fluorescent Timer, which is diluted by cell division before maturation, and the individual bacteria show a high [500 nm (green)/600 nm (red)] fluorescence (color) ratio (i.e. “green” fluorescence). By contrast, growth-arrested bacteria accumulate both green and red fluorescent Timer, and the individual bacteria show a low green/red color ratio (i.e., “red” fluorescence).

Here, coupling the protozoan model of infection, the innovative Timer fluorescent growth rate reporter and high-throughput single-cell technologies, we are defining a pipeline to track and collect the persisters produced by *L. pneumophila* with an unprecedented resolution. Isolated persisters can be analyzed biochemically and functionally. To date, using this methodology, we demonstrated that *L. pneumophila* persisters are a preestablished subpopulation of nongrowers and the expression of a regulated pathogen strategy which is applied to evolutionarily distant host cells (i.e., protists and macrophages) [15,16]. Near in situ biochemical characterization of the intracellular *L. pneumophila* persisters unveiled a unique physiology and, remarkably, a specific virulence program in order to survive the combined action of the host defenses and the antibiotics, hence, breaking the dogma that persisters are dormant bacteria [15]. Of note, this approach can be applied to the study nonreplicating persisters produced by other pathogens during the infection or within complex microbial communities.

#### **ACKNOWLEDGMENTS:**

Research at the CIRI is supported by the Inserm, CNRS, University Lyon 1, and ENS Lyon. N.P. was the recipient of an SNF Ambizione fellowship (PZ00P3\_161492, PZ00P3\_185529). The research in the laboratory of N.P. is supported by a CIRI seed fund and the ANR. The funders had no role in study design, decision to publish, or preparation of the manuscript.

#### **DISCLOSURES:**

The authors have nothing to disclose.

#### **REFERENCES:**

1. Bumann D, Cunrath O. Heterogeneity of Salmonella-host interactions in infected host tissues. *Current Opinion in Microbiology*. 2017;39: 57–63.  
doi:10.1016/j.mib.2017.09.008
2. Pham THM, Brewer SM, Thurston T, Massis LM, Honeycutt J, Lugo K, et al. Salmonella-Driven Polarization of Granuloma Macrophages Antagonizes TNF-

- Mediated Pathogen Restriction during Persistent Infection. *Cell Host & Microbe*. 2020;27: 54-67.e5. doi:10.1016/j.chom.2019.11.011
3. Davis KM, Mohammadi S, Isberg RR. Community Behavior and Spatial Regulation within a Bacterial Microcolony in Deep Tissue Sites Serves to Protect against Host Attack. *Cell Host & Microbe*. 2015;17: 21–31. doi:10.1016/j.chom.2014.11.008
  4. Claudi B, Spröte P, Chirkova A, Personnic N, Zankl J, Schürmann N, et al. Phenotypic variation of *Salmonella* in host tissues delays eradication by antimicrobial chemotherapy. *Cell*. 2014;158: 722–733. doi:10.1016/j.cell.2014.06.045
  5. Ackermann M. A functional perspective on phenotypic heterogeneity in microorganisms. *Nat Rev Microbiol*. 2015;13: 497–508. doi:10.1038/nrmicro3491
  6. Balaban NQ, Merrin J, Chait R, Kowalik L, Leibler S. Bacterial persistence as a phenotypic switch. *Science*. 2004;305: 1622–1625. doi:10.1126/science.1099390
  7. Harms A, Maisonneuve E, Gerdes K. Mechanisms of bacterial persistence during stress and antibiotic exposure. *Science*. 2016;354: aaf4268. doi:10.1126/science.aaf4268
  8. Bigger Joseph W. TREATMENT OF STAPHYLOCOCCAL INFECTIONS WITH PENICILLIN BY INTERMITTENT STERILISATION. *The Lancet*. 1944;244: 497–500. doi:10.1016/S0140-6736(00)74210-3
  9. Conlon BP, Rowe SE, Gandt AB, Nuxoll AS, Donegan NP, Zalis EA, et al. Persister formation in *Staphylococcus aureus* is associated with ATP depletion. *Nat Microbiol*. 2016;1: 1–7. doi:10.1038/nmicrobiol.2016.51
  10. Helaine S, Cheverton AM, Watson KG, Faure LM, Matthews SA, Holden DW. Internalization of *Salmonella* by Macrophages Induces Formation of Nonreplicating Persisters. *Science*. 2014;343: 204–208. doi:10.1126/science.1244705
  11. Kortebe M, Milohanic E, Mitchell G, Péchoux C, Prevost M-C, Cossart P, et al. *Listeria monocytogenes* switches from dissemination to persistence by adopting a vacuolar lifestyle in epithelial cells. *PLOS Pathogens*. 2017;13: e1006734. doi:10.1371/journal.ppat.1006734
  12. Manina G, Dhar N, McKinney JD. Stress and Host Immunity Amplify *Mycobacterium tuberculosis* Phenotypic Heterogeneity and Induce Nongrowing Metabolically Active Forms. *Cell Host & Microbe*. 2015;17: 32–46. doi:10.1016/j.chom.2014.11.016
  13. Moyed HS, Bertrand KP. *hipA*, a newly recognized gene of *Escherichia coli* K-12 that affects frequency of persistence after inhibition of murein synthesis. *Journal of Bacteriology*. 1983;155: 768–775. doi:10.1128/jb.155.2.768-775.1983
  14. Mulcahy LR, Burns JL, Lory S, Lewis K. Emergence of *Pseudomonas aeruginosa* Strains Producing High Levels of Persister Cells in Patients with Cystic Fibrosis. *Journal of Bacteriology*. 2010;192: 6191–6199. doi:10.1128/JB.01651-09
  15. Personnic N, Striednig B, Lezan E, Manske C, Welin A, Schmidt A, et al. Quorum sensing modulates the formation of virulent *Legionella* persisters within infected cells. *Nat Commun*. 2019;10: 5216. doi:10.1038/s41467-019-13021-8
  16. Personnic N, Striednig B, Hilbi H. Quorum sensing controls persistence, resuscitation, and virulence of *Legionella* subpopulations in biofilms. *ISME J*. 2021;15: 196–210. doi:10.1038/s41396-020-00774-0

17. Balaban NQ, Helaine S, Lewis K, Ackermann M, Aldridge B, Andersson DI, et al. Definitions and guidelines for research on antibiotic persistence. *Nat Rev Microbiol.* 2019;17: 441–448. doi:10.1038/s41579-019-0196-3
18. Newton HJ, Ang DKY, van Driel IR, Hartland EL. Molecular pathogenesis of infections caused by *Legionella pneumophila*. *Clin Microbiol Rev.* 2010;23: 274–298. doi:10.1128/CMR.00052-09
19. European Centre for Disease Prevention and Control. Legionnaires' disease. In: ECDC. Annual Epidemiological report for 2018. Stockholm: ECDC; 2020.
20. Brady MF, Sundareshan V. Legionnaires' Disease. StatPearls. Treasure Island (FL): StatPearls Publishing; 2022. Available: <http://www.ncbi.nlm.nih.gov/books/NBK430807/>
21. Pouderoux C, Ginevra C, Descours G, Ranc A-G, Beraud L, Boisset S, et al. Slowly or Nonresolving Legionnaires' Disease: Case Series and Literature Review. *Clin Infect Dis.* 2020;70: 1933–1940. doi:10.1093/cid/ciz538
22. Boamah DK, Zhou G, Ensminger AW, O'Connor TJ. From Many Hosts, One Accidental Pathogen: The Diverse Protozoan Hosts of *Legionella*. *Front Cell Infect Microbiol.* 2017;7: 477. doi:10.3389/fcimb.2017.00477
23. Isberg RR, O'Connor TJ, Heidtman M. The *Legionella pneumophila* replication vacuole: making a cosy niche inside host cells. *Nat Rev Microbiol.* 2009;7: 13–24. doi:10.1038/nrmicro1967
24. Burstein D, Amaro F, Zusman T, Lifshitz Z, Cohen O, Gilbert JA, et al. Genomic analysis of 38 *Legionella* species identifies large and diverse effector repertoires. *Nat Genet.* 2016;48: 167–175. doi:10.1038/ng.3481
25. Personnic N, Bärlocher K, Finsel I, Hilbi H. Subversion of Retrograde Trafficking by Translocated Pathogen Effectors. *Trends in Microbiology.* 2016;24: 450–462. doi:10.1016/j.tim.2016.02.003
26. Moliner C, Fournier P-E, Raoult D. Genome analysis of microorganisms living in amoebae reveals a melting pot of evolution. *FEMS Microbiology Reviews.* 2010;34: 281–294. doi:10.1111/j.1574-6976.2009.00209.x
27. Siddiqui R, Khan NA. Biology and pathogenesis of *Acanthamoeba*. *Parasit Vectors.* 2012;5: 6. doi:10.1186/1756-3305-5-6
28. Rao C, Benhabib H, Ensminger AW. Phylogenetic Reconstruction of the *Legionella pneumophila* Philadelphia-1 Laboratory Strains through Comparative Genomics. *PLoS One.* 2013;8: e64129. doi:10.1371/journal.pone.0064129
29. Samrakandi MM, Cirillo SLG, Ridenour DA, Bermudez LE, Cirillo JD. Genetic and Phenotypic Differences between *Legionella pneumophila* Strains. *J Clin Microbiol.* 2002;40: 1352–1362. doi:10.1128/JCM.40.4.1352-1362.2002

## MATERIALS

### Materials

- Isobiotec ball homogenizer (<http://www.isobiotec.com>)
- FACS Aria III sorter
- BBL yeast extract
- Bacto Proteose Peptone
- D(+)glucose monohydrate
- $\text{MgSO}_4$
- $\text{CaCl}_2$
- sodium citrate
- $\text{Fe}(\text{NH}_4)_2(\text{SO}_4)_2$
- $\text{Na}_2\text{HPO}_4$
- $\text{KH}_2\text{PO}_4$
- Bacto™ yeast extract
- agar
- L-cysteine
- $\text{Fe}(\text{NO}_3)_3$
- activated charcoal powder
- KOH
- N-(2-acetamido)-2-aminoethane-sulfonic acid
- glycerol solution
- chloramphenicol

## Cultivation of *Acanthamoeba castellanii*: Axenic cultivation...

1



### Note

Note 1: *A. castellanii* strain Douglas (ATCC 30234) is used. *A. castellanii* is an ubiquitous amoeba with a 2-stage lifecycle: an active trophozoite stage that exhibits vegetative growth and a stress induced dormant cyst stage with minimal metabolic activity [27]. Only *A. castellanii* trophozoites are suited for infection.

Note 2: **Warning.** *A. castellanii* is an opportunistic pathogen producing rare and severe human infections. Therefore, it must be manipulated in a biosafety containment level 2 in a microbial safety cabinet in compliance with local rules.

**Prepare Peptone Yeast Glucose (PYG) broth (see Note 3):**

A	B
---	---

A	B
BBL yeast extract	1 g/L
Bacto Proteose Peptone	20 g/L
D(+)-glucose monohydrate (50 mL of a 2 M	1.8% (w/v)
MgSO <sub>4</sub> (10 mL of a 0.4 M solution)	4 mM
CaCl <sub>2</sub> (8 mL of a 0.05 M solution)	0.4 mM
sodium citrate (3.4 mL of a 1 M solution)	3.4 mM
Fe(NH <sub>4</sub> ) <sub>2</sub> (SO <sub>4</sub> ) <sub>2</sub> (20 mg)	0.05 mM
Na <sub>2</sub> HPO <sub>4</sub> (10 mL of a 0.25 M solution)	2.5 mM
KH <sub>2</sub> PO <sub>4</sub> (10 mL of a 0.25 M solution)	2.5 mM

Add the components (except the glucose) to 950 mL H<sub>2</sub>O, adjust the pH with 1 Molarity (M) HCl to 6.5±0.1 and autoclave. Solubilize 11 g of D(+)-glucose in 50 mL of warmed H<sub>2</sub>O (37 °C) and add to the autoclaved medium (see **Note 4**). Sterilize the medium using a 0.2 µm filter cartouche and store at 4 °C.

#### Note

***Note 3:** For a high reproducibility of growth characteristics, the components should be tested beforehand, and the same suppliers and batches should be used for all experiments.*





***Note 4:** Glucose caramelizes upon autoclaving in combination with the medium. Suspend the D(+)-glucose in warmed 50 mL of H<sub>2</sub>O, filter sterilize (0.2 µm) and add to the autoclaved medium.*

- Thaw an ampule of frozen amoeba at Room temperature and aseptically transfer the contents into a 13 mL sterile Falcon round-bottom polypropylene tubes containing 10 mL of PYG. Centrifuge 250 x g, 00:10:00 and discard the supernatant. Resuspend the pellet in 10 mL of PYG medium to inoculate a tissue culture treated flask (75 cm<sup>2</sup>: 10-15mL medium) for initial amplification.
- Cultivate the amoebae at 23 °C in PYG medium optionally supplemented with antibiotics (e.g. Pen/Strep or Fungizone) to avoid bacterial and fungal contaminations. Use either 100 mm tissue culture treated Petri dishes (minimum of 10 mL medium) or tissue culture treated flasks (25 cm<sup>2</sup>: 5-7 mL medium, 75 cm<sup>2</sup>: 10-15mL medium). Change the medium every 2-3 days. Split the culture when cell reach 80% confluence by repeated pipetting of the media over the



plate.

## Cultivation of *Acanthamoeba castellanii*: Cryopreservation ...

- 4 Remove the PYG medium and harvest the amoeba, when in logarithmic growth phase, by repeated pipetting of fresh refrigerated PYG medium over the plate. Determine the cell concentration by counting the amoeba with a haemocytometer. Adjust the concentration to about  $10^7$  amoeba/mL.
- 5 Prepare freshly a 15% solution of sterile dimethyl sulfoxide (DMSO) in refrigerated PYG medium and mix to the cell suspension in equal portions to reach a final concentration between  $10^6$  and  $10^7$  cells/mL and 7.5 % DMSO.
- 6 Dispense in  0.5 mL aliquots into 2.0 mL sterile plastic screw-capped CryoTube vials. Place the vials in a Nalgene Cryo  1 °C freezing container and store at  -80 °C ,  Overnight .

### Note

For long-term storage, the cell preparation shall be deposited in either the liquid or the vapor phase of a nitrogen freezer.

## Generation of a Timer producing *Legionella pneumophila*

7



### Note

**Note 5: *Warning.*** *L. pneumophila JR32 is used. This lab strains originally derived from a Philadelphia-1 clinical strain isolated from an outbreak in Philadelphia in 1976 [28,29]. L. pneumophila is the causing agent of the Legionnaires' disease and must therefore be manipulated in a biosafety containment level 2 in a microbial safety cabinet in compliance with local rules.*

**Note 6:** *L. pneumophila D1cmT lacks a functional T4SS and therefore cannot grow within infected cells. It is a convenient control to define the Timer spectral properties corresponding to the lack of bacterial division during the infection.*

**Prepare CYE (charcoal yeast extract) agar plates:**

A	B
ACES	10 g/L
Bacto™ yeast extract (see Note 7)	10 g/L
activated charcoal powder (puriss. p.a.)	2 g/L
agar	15 g/L
L-cysteine	3.3 mM
Fe(NO <sub>3</sub> ) <sub>3</sub>	0.6 mM

Dissolve 10 g of ACES and 10 g of yeast extract in 950 mL of H<sub>2</sub>O and adjust the pH to 6.9 with 10 Molarity (M) KOH. Transfer the solution to a 1 L Schott bottle containing 2 g of activated charcoal powder, 15 g of agar and a stir bar. Autoclave and let the agar solution cool down to 50 °C. Add filter sterilized 0.4 g / 10 mL L-cysteine and 0.25 g / 10 mL Fe(NO<sub>3</sub>)<sub>3</sub> solutions (see Note 8). If required add antibiotics. Mix the solution on a magnetic stirrer and pour plates (approximately 40 plates per L of medium). Let the plates dry for 1 day at Room temperature and store at 4 °C.

#### Note












*Note 7: The source and quality of yeast extract and peptone affect the physiology of L. pneumophila. For a high reproducibility of virulence traits and growth characteristics, the components should be tested beforehand, and the same suppliers and batches should be used for all experiments.*

*Note 8: Dissolve L-cysteine and Fe(NO<sub>3</sub>)<sub>3</sub> each separately in 10 mL of H<sub>2</sub>O in a 15 mL tube. Stir the medium and slowly add the L-cysteine solution first, followed by the iron solution to prevent precipitation.*

## 8 Prepare ACES yeast extract (AYE) broth:



A	B
N-(2-acetamido)-2-aminoethane-sulfonic acid (ACES)	10 g/L
Bacto yeast extract	10 g/L
L-cysteine	3.3 mM
Fe(NO <sub>3</sub> ) <sub>3</sub>	0.6 mM

Add  10 g of ACES and  10 g of yeast extract in  950 mL of H<sub>2</sub>O. Add filter sterilized  0.4 g /  10 mL L-cysteine and  0.25 g /  10 mL Fe(NO<sub>3</sub>)<sub>3</sub> solutions. Adjust the pH to  6.9 with  10 Molarity (M) KOH. If required, add antibiotics. Pass the medium several times through a glass-fibers filters, and finally sterilize the medium using a  0.2 µm filter cartouche. Store the medium at  4 °C in the dark (see **Note 9**).

#### Note

*Note 9: Pre-filter the medium 6 to 8 times through a glass microfibers filter paper to remove precipitates. L-cysteine is light sensitive.*


## Generation of a Timer producing *Legionella pneumophila*: P..

9



#### Note


*Note 10: Perform all steps at room temperature unless stated otherwise.*

Streak out *L. pneumophila* wild-type and the isogenic mutant *DicmT* from frozen glycerol stocks onto CYE plates. Colonies or a bacterial lawn will form after 3 days of incubation at  30 °C (see **Note 11**).

#### Note


*Note 11: Electrocompetence of *L. pneumophila* is increased when using freshly prepare bacteria grown at 30°C.*


10

Generously harvest the bacterial lawn with a disposable inoculation loop and inoculate  1 mL of PBS in a 1.5 mL sterile Eppendorf tube.

11




Centrifuge at  5000 rpm, 00:10:00, discard the supernatant and resuspend the bacterial pellet with a 10 % glycerol solution. Repeat the operation 3-4 times. 10m

12 Finally, resuspend the bacterial pellet with  100 mL of a 10% glycerol solution.



13 Directly proceed with the electroporation.

## Generation of a Timer producing *Legionella pneumophila*: E..



14 Add  100 ng of the plasmid pNP107 to the bacterial suspension and transfer into a 2 mm electroporation cuvette.



15 Electroporate using the following settings: U=2.5 kV; C= 25  $\mu$ F; 200  $\Omega$ = $\infty$ ; t.c.  $\approx$ 2.5.



16 Immediately add  500  $\mu$ L of pre-warmed AYE, transfer into a 13 mL Falcon round-bottom polypropylene tubes and incubate on a shaker set at  200 rpm, 37°C, 05:00:00 .



17 Plate on CYE agar supplemented with  5  $\mu$ g/mL chloramphenicol (*see* **Note 12**). Orange colonies (*i.e.*, producing the Timer fluorescent protein) will form after 3 days of incubation at  37 °C .





### Note

*Note 12: Chloramphenicol stock:*  30 mg/mL ethanol. Store the stock solutions at  -20 °C .

## Preparation of infectious *Legionella pneumophila*






18



Streak out Timer producing *L. pneumophila* from frozen glycerol stocks onto CYE plates supplemented with  5 µg/mL chloramphenicol. Colonies or a bacterial lawn will form after 2-3 days of incubation at  37 °C.

19



Prepare planktonic cultures for infection (see **Note 13**). Suspend a loop of the bacterial lawn grown on CYE/Cam agar plates for 2-3 days in  600 µL of AYE medium supplemented with  5 µg/mL chloramphenicol (AYE/Cam) in a 2 mL Eppendorf tube. Mix thoroughly until you reach a homogeneous suspension (see **Note 14**). Measure the OD<sub>600 nm</sub> and adjust the density of the bacterial suspension to an OD<sub>600 nm</sub> of 0.1. Inoculate  9 mL of AYE/Cam with the suspension to obtain an OD<sub>600 nm</sub> of 0.1 and split into 3 x  3 mL in 13 mL Falcon™ round-bottom polypropylene tubes and incubate on a shaker set at  200 rpm, 37°C, 21:00:00.


#### Note

***Note 13:** *L. pneumophila* is most virulent (so-called transmissive) at the transition from the exponential to the stationary growth phase. The appearance of brown pigmentation is an indicator for stationary growth phase.*

***Note 14:** Day 3 bacterial lawn is difficult to dissolve in broth. To prevent the formation of aggregates, first collect the bacterial lawn in a 2 mL Eppendorf tube. Second, add the culture medium and wait for a few minutes prior to mixing.*

20



At 21 h, *L. pneumophila* must be transmissive. Check the bacteria by examining  10 µL of culture loaded on a glass slide covered with a 12 mm round glass coverslip with an inverted light microscope (40 × objective). *L. pneumophila* look like a dense population of small and motile coccobacilli (see **Note 15**). Measure the OD<sub>600 nm</sub>; an OD<sub>600 nm</sub> of 5 corresponds to approximately  $2 \times 10^9$  bacteria/mL (see **Note 16**).

#### Note

***Note 15:** Presence of elongated bacilli or filaments signals sub-optimal culture quality.*

***Note 16:** Under our conditions an OD<sub>600 nm</sub> of 5.0 corresponds to approximately  $2 \times 10^9$  bacteria/mL. At 21 h, the exact correlation of the OD<sub>600 nm</sub> with the bacterial concentration depends on the spectrophotometer used and must be determined experimentally.*

21

Proceed with infection of *A. castellanii*.

## Infection of *Acanthamoeba castellanii*

- 22 Split *A. castellanii* cultures one day before the experiment to obtain cells in logarithmic growth phase, which have been fed the day before (see **Note 17**).


### Note

*Note 17: For best and most reproducible results use A. castellanii growing exponentially. Non-dividing, dense or starved cells will be infected less efficiently and less reproducibly.*

- 23 Wash exponentially growing *A. castellanii* with PYG and detach the cells by repeated pipetting of the media over the plate.





- 24 Determine cell concentration by counting the cells with a haemocytometer.

- 25 Seed cells at a density of  $9 \times 10^5$  amoeba per well of a tissue culture treated 6-well plate. Let the phagocytes settle and adhere to the plastic surface for  24:00:00 under normal culture conditions (see **Note 18**).

1d

### Note

*Note 18: Use  2 mL of PYG medium per well. The next day, cell number will have doubled. Alternatively, cell can be seeded the day of the infection as cell adhesion to the plastic surface takes  00:30:00.*

- 26 Prior to infection, remove the PYG medium and wash out antibiotics if necessary (see **Note 19**). Determine the cell numbers per well using a haemocytometer. Calculate the appropriate number of bacteria to infect host cells with a multiplicity of infection (MOI) of 2.

## Note

*Note 19: Cell culture media used in assays must not contain antibiotics.*

27



Synchronize the infection by centrifugation of the 6-well plate for **250 x g, 00:10:00** and **1d 0h 10m** incubate for 45 to 60 min at **25 °C**. Remove not internalized bacteria by washing.

- Wash the infected cells with PYG. (1/3)
- Wash the infected cells with PYG. (1/3)
- Wash the infected cells with PYG. (1/3)
- Incubate at least **24:00:00** at **25 °C** (see **Note 20 & 21**).

## Note

*Note 20: Subpopulations of proliferative *L. pneumophila* and nonreplicating persisters of are distinctly observed and were characterized at 24 h post-infection (p.i.), in infected *A. castellanii*.*

*Note 21: Infected *A. castellanii* loses motility and develops tomette-like shape.*

28

Proceed with the sample preparation for FACS.

## Sample preparation for FACS

29

## Note


*Note 22: After host cell lysis, maintain the samples at 4°C to slow down the Timer maturation. Thereby the Timer fluorescent ratio will accurately correlate to the bacterial growth rate for extended period of times (several hours).*


*Note 23: Optimal lysis of infected *A. castellanii* requires a 2-step process combining detergent and mechanical forces.*



**Prepare homogenization (HS) buffer:**

A	B
N-2-hydroxyethylpiperazine-N-2-ethanesulfonic acid (HEPES)	20 mM




A	B
Sucrose	250 mM
Ethyleneglycoltetraacetic acid (EGTA)	0.5 mM
pH adjusted to 7.2 with KOH	1 M

Filter sterilize and store at  4 °C .




**30**  On ice , detach the infected amoeba by repeated pipetting of the media over the wells and transfer into 2 mL Eppendorf tubes.

**31** Pellet by centrifugation  250 x g, 4°C, 00:10:00 , remove the supernatant and wash once with  1 mL ice-cold PBS. 10m




**32**  On ice , lyse the cell pellet with  1 mL 0.1 % Triton TX-100 (Sigma) in HS buffer for  00:10:00 . 10m



**33** Centrifuge the cell lysate  250 x g, 4°C, 00:15:00 , remove the supernatant and resuspend the pellet with  1 mL PBS. Keep  On ice . 15m



**34** Assemble and place the Isobiotec ball homogenizer (<http://www.isobiotec.com>) (6 or 8 mm clearance ball)  On ice . Mount a pair of Luer-lock syringes.

**35** Rinse well the ball homogenizer.




- 35.1** Remove the plunger of one syringe and fill the volume with ddH<sub>2</sub>O. Install the plunger and flush. (1/3)
- 35.2** Remove the plunger of one syringe and fill the volume with ddH<sub>2</sub>O. Install the plunger and flush. (2/3)
- 35.3** Remove the plunger of one syringe and fill the volume with ddH<sub>2</sub>O. Install the plunger and flush. (3/3)
- 35.4** Finally flush the homogenizer once with ice-cold HS buffer. Make sure there is no leakage or clogging.

- 36** Install a new pair of Luer-Lock syringes. Remove the plunger of one syringe to load the cell lysate. Install the plunger and press the sample through the homogenizer into the second syringe. Press back and forth nine times (see **Note 24**).



#### Note

***Note 24: Warning.** Homogenization generates aerosol that poses a threat to the operator when working with respiratory pathogens. Wear mask and goggles. Operate in biosafety containment level 2 in a microbial safety cabinet. Notify your colleagues.*

- 37** Collect the homogenized sample and transfer through the cell-strainer cap into a 5 mL polystyrene round-bottom tube (see **Note 25**). Keep  On ice.

#### Note

***Note 25:** The cell-strainer cap retains oversize particles and therefore protect the cytometer microfluidics from clogging.*

- 38** Proceed with the isolation of the nonreplicating persisters.

## High-throughput isolation of nonreplicating persisters by FA.

39

### Note

***Note 26:** The procedure detailed below was established using a high speed multicolor FACSAria III 4L sorter equipped with 4 spatially separated laserlines, for excitation, and the associated detector arrays: 405 nm (violet), 488 nm (blue), 561 nm (yellow/green) and 633 nm (red). Data acquisition using the FACSDiva software.*

***Note 27: Warning.** FACS generates aerosol that poses a threat to the operator when working with respiratory pathogens. Local rules may impose a **biosafety containment level 3** when working on aerosolized *L. pneumophila*. Ensure the sorter is equipped with an aerosol management system to evacuate the sort collection chamber and to trap aerosolized particles during sorting. Wear mask and goggles. Notify your colleagues.*

Start up the sorter.

- 39.1** Decontaminate the sort block, the sort chamber, the sample injection port and the sort collection device (tube holder and lid) with 70% ethanol. Clean the deflecting plates with ddH<sub>2</sub>O. Wipe to dry carefully all surfaces that will be charged.

- 39.2** Assemble and filled the sheath tank with sterile PBS (see **Note 28 & 29**).

### Note

***Note 28:** Use autoclaved tank.*

***Note 29:** The sheath level sensor must be decontaminated with 70 % ethanol.*

- 39.3** Turn on the aerosol management system.

- 39.4** Turn on the sorter (see **Note 30**).

### Note




***Note 30:** Ensure the excitation laser 488 nm and 561 nm are set to "ON".*

- 39.5** In the FACSDiva menu, turn on the fluidic system by selecting “Fluidic startup”(see **Note 31**) and follow the instructions.

**Note**

*Note 31: Ensure that the lid of the sheath tank is properly closed and the tank pressurized.*

- 39.6** In the sort block, remove the close-loop nozzle and insert the nozzle with an orifice of 70 mm at the dedicated position at the lower end of the cuvette flow cell. Set the correct configuration for the 70 mm nozzle in the FACSDiva menu “view configurations”.

- 39.7** Turn on the external temperature control system for the collection tube holder and set the temperature to  4 °C . In the FACSDiva menu “Cytometer”, set the agitation of the sample in the injection chamber to  300 rpm and temperature to  4 °C .

- 40** Establish a stable stream and set up the drop formation.

**Note**

*Note 32: During sorting, drop drive energy is applied to the stream to break it into highly uniform droplets containing particle of interest. Droplets detach from the stream a few millimeters downstream from the nozzle*

*Note 33: In the FACSDiva software, use the breakoff window to control the stream and the drop formation.*

- 40.1** Start the stream that should appear in the breakoff window.


- 40.2** Within the sort block, ensure the stream hit the center of the waste aspirator (see **Note 34**).

**Note**

*Note 34: The entire sort block assembly can be rotated on a fixed pivot point to adjust the position of the stream in the waste aspirator.*

- 40.3** Stabilize the position of the drop breakoff point by adjusting the amplitude (see **Note 35**). Drops before the breakoff point will appear pear-shaped. Ensure satellite drop merge with leading drop.

**Note**

***Note 35:** Each nozzle size configuration is optimized to a preset sheath pressure and also change the related values for drop drive frequency, drop charge levels, laser delay and area scaling. Keep the preset values unless you cannot establish a stable stream (i.e. pressure:  70 undetermined ; gap: 6 and frequency: 87 KHz).*

- 40.4** Copy the generated Drop 1 value (i.e., the distance between the top of the image and the center of the first broken-off drop) in the target field.

- 40.5** Turn on the Sweet Spot to maintain the drop breakoff point during the sorting.

- 40.6** Wait 15-20 min and control the stability of the stream and drop formation. For the Drop 1, a difference of 10 units is acceptable. For the GAP a difference of 1 unit is acceptable.

- 41** Adjust the drop delay.

**Note**

***Note 36:** The drop delay calibrates the time between the particle analysis and the droplet formation in order to charge the drop the sort.*

***Note 37:** Side streams, deflection, and drop delay are defined using the side stream window in the FACSDiva software.*

***Note 38:** Sorting is controlled using the sort layout window in the FACSDiva software.*

- 41.1** Ensure the maximal illumination of the central stream using the micrometer dial.

**41.2** Turn on the deflection plates (voltage) in the side stream window (*see* **Note 39**).

**Note**

*Note 39: Warning. High voltage potential. Do not touch the deflection plates when the voltage is ON (red warning light). Keep the sort block door closed.*

**41.3** In the FACSDiva menu, create an AccuDrop experiment.

**41.4** Place a tube filled with a dilute suspension of AccuDrop beads on the sample injection port and load into the sample injection chamber (*see* **Note 40**). Adjust the flow rate to reach a threshold rate of 2'000-3'000 events per second (*see* **Note 41**).

**Note**

*Note 40: Tube loading, data acquisition and recording as well as flow rate are controlled using the FACSDiva acquisition dash board.*

*Note 41: Avoid applying high flow rate.*

**41.5** Activate the test sort in the side stream window.

**41.6** Activate the optical Filter in the side stream window.

**41.7** Adjust the voltage slider to place the left side stream in the center of the box on the left. Ensure that the central stream is visible in the right box (*see* **Note 42**).

**Note**

*Note 42: Adjust the 2nd, 3rd, and 4th Drop settings to tighten the center stream and fine-tune the side streams, if needed.*

- 41.8** Stop the test sort in the side stream window.
- 41.9** In the FACSDiva sort layout window select “2 tube” as device and “Initial” as a precision mode.
- 41.10** Start sorting the beads by selecting “sort” in the sort layout window (see **Note 43**).
- Note**

*Note 43: Keep the aspirator drawer closed.*
- 41.11** In the side stream window, adjust the drop delay value in 1-drop increments to achieve close to 100% intensity in the left side stream (see **Note 44**).
- Note**

*Note 44: Alternatively, the drop delay can be set automatically using the Auto Delay button in the side stream window.*
- 41.12** In the sort layout window, change the precision mode to “Fine Tune” and continue optimizing the drop delay by small increments until the left side stream intensity is > 95 % (see **Note 45**).
- Note**

*Note 45: Wait a few seconds after each small increment for a complete response to the delay change.*
- 41.13** In the side stream window, inactivate the Optical Filter and set the voltage slider to “0”. In the sort layout, stop the sort. Close the AccuDrop experiment.
- 41.14** Unload the tube filled with AccuDrop beads using the FACSDiva acquisition dash board. Discard the tube or keep it in a fridge for future sort.

## 42 Set up the side streams for the 2-way sort.

**42.1** Insert one 15 mL tube in the 2-way sort tube holder. Close with the universal top. Insert the sort collection device into the sort collection chamber and connect to the temperature control system in order to refrigerate the collection tube.

**42.2** In the side stream window, turn on the deflection plates (voltage) and activate “test sort”. Optimize the position of the first either right or left side stream using the corresponding voltage sliders (see **Note 46**).

### Note

*Note 46: Keep the aspirator drawer closed.*

**42.3** In the side stream window, open the aspirator drawer to confirm the side stream are properly aligned with the collection 15 mL tube.

**42.4** Close the aspirator drawer, inactivate “test sort”, turn off the voltage, remove the sort collection device and trash the 15 mL collection tube.

## 43 Create the gates to sort nonreplicating persisters.

### Note

*Note 47: To ensure that your cytometer is performing consistently over time, you should regularly run performance check.*

*Note 48: Acquisition setup and gating strategy are carried out on the FACSDiva software workplace.*

**43.1** In FACSDiva browser, open a new experiment and click on “Cytometer settings”.

- 43.2** In the Cytometer window, in the parameter tab, delete all unneeded optical configuration. Keep forward scatter (FSC); side scatter(SSC); [Laser 488 nm - 502LP - BP 530/30] (FITC) and [Laser 561 nm - 600LP - BP 610/20] (mCherry).
- 43.3** For each parameter, only record the height (H) of the pulse signal detected.
- 43.4** In the cytometer window, in the threshold tab, set both the FSC and SSC values to the minimum (i.e., 200).
- 43.5** On the global worksheet, create (i) a SSC-H vs FSC-H plot; (ii) a FITC-H vs FSC-H plot and (iii) a mCherry-H vs FITC-H plot.
- 43.6** In the inspector window, select the checkbox for X Axis and Y Axis under biexponential display.
- 43.7** Place the sample tube filled with homogenized lysates on the sample injection port and load into the sample injection chamber using the FACSDiva acquisition dash board and start acquisition (see **Note 49**). Detected signals appear in the plots created on the global worksheet.

**Note**

*Note 49: In the FACSDiva browser / experiment / specimen, make sure that the tube is activated.*

- 43.8** In the cytometer window, parameter tab, set the photomultiplier (PMT) voltage for FSC and SSC in order to center the detected signal on the SSC-H vs FSC-H plot.



**43.9** In the cytometer window, parameter tab, set the PMT voltage for FITC-H in order to separate the green Timer signal from the host debris autofluorescence and the electronic noise on the FITC-H vs FSC-H plot (see **Note 50**). Using the polygonal gating tools, draw the **Gate 1** around the FITC-H positive population.

**Note**

*Note 50: The use of homogenized lysate from cell infected with non-fluorescent L. pneumophila confirm the absence of any background particles with similar fluorescence.*

**43.10** In the mCherry-H vs FITC-H plot, apply Gate 1. In the cytometer window, parameter tab, set the PMT voltage for mCherry-H in order to center the red Timer signal (see **Note 51**).

**Note**

*Note 51: The use of homogenized lysate from cell infected with mcherry or GFP producing L. pneumophila shall confirm the absence of fluorescence overlap in themCherry-H vs FITC-H plot. When using a sorter without spatially separated laserlines and detector arrays possible spectral overlap must be corrected via a process of fluorescence compensation.*

**43.11** In the mCherry-H vs FITC-H plot, two well-separated subpopulations with similar FITC-H signal intensity appear. The population with the higher mCherry intensity corresponds to the nonreplicating persisters. Using the gating tools, draw the **Gate 2** around this subpopulation (see **Note 52 and 53**).

**Note**

*Note 52: L. pneumophila DcmT cannot replicate intracellularly. To control the gating strategy, use homogenized lysate from cell infected with Timer producing L. pneumophila DcmT mutant. All bacteria shall be contained in **Gate 2**.*

*Note 53: On the global worksheet verify the gate are properly hierarchized by showing the "population hierarchy view".*

**43.12** In the FACSDiva acquisition dash board, select Gate 2 as stopping gate and 10'000 as events to record. Start recording by clicking "record" (see **Note 54**).

#### Note


*Note 54: Adjust the flow rate to reach a threshold rate < 12'000 events per second.*

- 43.13** In the mCherry-H vs FITC-H plot, using the population hierarchy view, determine the population size in Gate 2 (% Parent).

- 44** Sort the nonreplicating persisters.

#### Note

*Note 55: Make sure the aerosol management system is ON and both the sort block and sort chamber tightly closed in order to avoid aerosols diffusion.*

- 44.1** Insert one 15 mL tube prefilled with  3 mL ice-cold PBS in the 2-way sort tube holder. In the sort layout window, define the sort location field to align the collection tube with the side stream as setup in 5.3. In the selected sort location field, add the subpopulation to be sorted (i.e., **Gate 2**). Close with the universal top. Insert the sort collection device into the sort collection chamber (as in 42.1.).

- 44.2** In the FACSDiva sort layout window select “2 tube” as device and “purity” as a precision mode (see **Note 56**).

#### Note

*Note 56: Precision modes define the sorting stringency. Sorting in Purity mode results in a sorted sample that is highly pure, at the expense of recovery and yield.*

- 44.3** In the side stream window, turn on the deflection plates (voltage), open the aspirator drawer and initiate the sample flow (“acquire”).

- 44.4** In the FACSDiva sort layout window, lunch the sort by clicking “sort”.

**44.5** In the sort layout window, monitor the sorting by observing the sorting rate, the conflicts rate and the sorting efficiency (ideally above 95 %). Inspect regularly the central and side stream in the side stream window.

**44.6** Sort until the required number of nonreplicating persisters has been sorted (see **Note 57-58**).


**Note**

***Note 57:** For continuous sorting, select “Continuous” from the Target Events menu in the sort layout window. If sample reloading is requested, pause the sorting, by clicking the “Pause” button in the Sort Layout. The counting of sorted cells is retained when you restart sorting by clicking the “Resume” button.*

***Note 58:** Open the sort chamber to visually inspect the 15 mL collection tube filling to avoid spill over. When necessary pause the sort and replace the 15 mL collection tube. Write on the tube the number of nonreplicating persisters sorted and store on ice.*

**44.7** Stop the sort by inactivating “sort” in the sort layout window and save the “sort report”. In the side stream window, control the deflection plates (voltage) are turned off and the aspirator drawer closed.

**44.8** Unload the sample tube using the FACSDiva acquisition dash board and eliminate into the dedicated biological wastes disposal in compliance with local rules.


**44.9** Open the sort chamber, remove and disassemble the sort collection device to collect and the store  On ice the 15 mL collection tube.

**45** Re-analysis of the sorted nonreplicating persisters.

**45.1** In the FACSDiva browser, under experiment and specimen create a new tube (see **Note 59**).

#### Note

*Note 59: A new tube can be created directly from the FACSDiva acquisition dash board.*

- 45.2** Harvest only  100 mL of sorted nonreplicating persisters from the 15 mL collection tube and transfer into a sample tube.
- 45.3** Place the sample tube on the sample injection port and load into the sample injection chamber using the FACSDiva acquisition dash board and activate acquisition. Detected signals appear in the plots created on the global worksheet.
- 45.4** In the FACSDiva acquisition dash board, adjust the flow rate to reach a threshold rate of 2'000-3'000 events per second.
- 45.5** In the FACSDiva acquisition dash board, select Gate 2 as stopping. Start recording by clicking "record". Stop recording and acquisition after 500-1000 events are recorded (see **Note 60**).

#### Note

*Note 60: In the mCherry-H vs FITC-H plot, sorted Timer producing nonreplicating persisters shall be contained in the **Gate 2**.*

- 45.6** In the mCherry-H vs FITC-H plot, using the population hierarchy view, determine the population size in Gate 2 (% Parent) to infer the sort purity.
- 45.7** Unload the sample tube using the FACSDiva acquisition dash board and eliminate into the dedicated biological wastes disposal in compliance with local rules.
- 45.8** Save your data and proceed with the shut-down.

## 46 Shut-down.

**46.1** Clean by sequentially loading and running sample tubes filled with FACSClean FACSRinse, 70% ethanol and finally sterile ddH<sub>2</sub>O.

**46.2** In the side stream window, turn off the stream, depressurized and disassemble the PBS tank (see **Note 61**).

### Note

*Note 61: Autoclaved the tank, decontaminate the sheath level sensor with 70 % ethanol.*

**46.3** Decontaminate the sort block, the sort chamber, the sample injection port, the sort collection device (tube holder and lid) with 70% ethanol. Clean the deflecting plates with ethanol and ddH<sub>2</sub>O. Wipe to dry carefully all surfaces.

**46.4** Assemble and filled the ethanol tank with 70% ethanol and connect to the sorter fluidics system (see **Note 62**).

### Note

*Note 62: Ensure that the lid of the ethanol tank is properly closed and the tank pressurized.*

**46.5** In the FACSDiva menu, turn off the fluidic system by selecting “Fluidic shutdown” and follow the instructions.

**46.6** In the sort block, remove the 70 mm nozzle. Decontaminate the nozzle and its slot at the lower end of the cuvette flow cell with 70 % ethanol. Insert the close-loop nozzle.

**46.7** Turn off the external temperature control system and the aerosol management system.

**46.8** Turn off the sorter.